MEDISANA®

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READ ALL INSTRUCTIONS BEFORE USE.
SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE.

IMPORTANT SAFEGUARDS:

WHEN USING ELECTRICAL PRODUCTS, ESPECIALLY WHEN CHILDREN ARE PRESENT, BASIC SAFETY PRECAUTIONS SHOULD ALWAYS BE FOLLOWED, INCLUDING THE FOLLOWING:

- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
 - NEVER insert pins or other metallic fasteners into the appliance.
 - Use this appliance for the intended use as described in this booklet. DO NOT use attachments not recommended by MEDISANA.
 - Excessive use could lead to the product's excessive heating and shorter life. Should this occur, discontinue use and allow the unit to cool before operating.
 - NEVER drop or insert any object into any opening.
 - DO NOT operate where aerosol (spray) products are being used or where oxygen is being administered.
 - DO NOT operate under a blanket or pillow. Excessive heating can occur and cause fire, electrocution or injury to persons.
 - DO NOT attempt to repair the appliance. There are no user serviceable parts. For service, send to a MEDISANA Service Centre. All servicing of this appliance must be performed by authorised MEDISANA service personnel only.
 - FOR HOUSEHOLD USE ONLY.
 - Do not place or store appliance where it can fall or be pulled into a bath or sink.
- Do not place in or drop into water or other liquid.**

Warranty and repair terms

Please contact your supplier or the service centre in case of a claim under the warranty. If you have to send in the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for MEDISANA products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be corrected free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
 - a. All damage which has arisen due to improper treatment, e.g. nonobservance of the user instructions.
 - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

This guarantee does not affect your statutory rights.

MEDISANA AG
Jagenbergstraße 19
41468 NEUSS
GERMANY

SAFETY CAUTIONS:

PLEASE READ THIS SECTION CAREFULLY BEFORE USING THE APPLIANCE.

- If you are susceptible to allergic reactions or if you have circulatory problems, consult a physician before using this product.
- Do not use on irritated, sunburned or chapped skin, or on any area with dermatological problems.
- Do not apply over skin eruptions, warts, moles or varicose veins.
- Use should be easy and painless. If you experience nail or skin irritation, or hand discomfort when using this product, discontinue use. If a more extreme reaction occurs, consult your doctor
- Individuals with pacemakers and pregnant women should consult a doctor before using this appliance. Not recommended for use by Diabetics.
- DO NOT use on an infant, invalid or on a sleeping or unconscious person. DO NOT use on insensitive skin or on a person with poor blood circulation.
- If you feel any discomfort whilst using this appliance, discontinue use and consult your doctor.
- Ensure that your hands are dry when operating the controls.
- This product should NEVER be used by any individual suffering from any ailments that would limit the user's capacity to operate the controls.
- Please ensure that all hair, clothing and jewellery are kept clear of moving parts of the product at all times.

Safety notes for batteries

- Do not disassemble batteries! • Never leave any low battery in the battery compartment since they may leak and cause damage to the unit. • Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of these parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately! • If a battery has been swallowed seek medical attention immediately! • Replace all of the batteries simultaneously! • Only replace with batteries of the same type, never use different types of batteries together or used batteries with new ones! • Insert the batteries correctly, observing the polarity! • Remove the batteries from the device if it is not going to be used for an extended period! • Keep batteries out of children's reach! • Do not attempt to recharge these batteries! There is a danger of explosion! • Do not short circuit! There is a danger of explosion! • Do not throw into a fire! There is a danger of explosion! • Rechargeable batteries are to be removed from the appliance before being charged! • Do not throw used batteries into the household refuse; put them in a hazardous waste container or take them to a battery collection point, at the shop where they were purchased!

Battery directive

 This symbol indicates that batteries must not be disposed of in the domestic waste as they contain substances which can be damaging to the environment and health. Please dispose of batteries in designated collection points.

WEEE explanation

 This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

PRODUCT FEATURES

INSTRUCTIONS FOR USE:

Inserting the batteries

1. Open the battery compartment on the underside of the device.
2. Insert batteries observing the polarity signs (+ and -) inside the battery compartment.
NOTE: The batteries are stacked with one of the batteries behind the other two.
3. Close the battery compartment.

How to use

1. Ensure that your feet are clean and dry.
2. Use the switch to select a power setting;
For sensitive skin, or softer skin use the LO setting.
For harder skin use the HI setting.
3. Using light pressure place the surface of roller onto the heel or sole of your foot and gradually move the device around to target areas of hard skin. Do not use on a single area for more than 4 seconds at a time.
4. Check the results throughout. Once you are happy with the results turn the switch to OFF.
5. Rinse your feet with clean water & dry before applying your favourite foot cream.

Changing the rollers

To remove the roller.

1. Ensure the switch is in the OFF position.
2. With one hand press the release button on the side of the device.
3. With the other hand lift out the roller. TIP: It's easier to lift the roller out from the end closest to the release button.

To insert a new roller.

1. Loosely position one of the roller pins into the recess on the end furthest from the release button.
2. Push the other end of the roller into place using the groves to guide it into position.
3. The roller will secure into place with a click.

MAINTENANCE AND STORAGE

- Always turn the device off before cleaning it.
- Only clean the device with the supplied brush or a dry cloth. Do not use chemical agents including benzene, alcohol etc.
- Remove batteries if equipment is not going to be used for extended periods of time, and store batteries in a cool dry place away from direct sunlight and metal objects.
- If motor power seems reduced, replace the batteries.
- If roller surface becomes worn or less abrasive, fit a new roller.

Cautions

Before use, make sure that the roller is in good condition. Any malformation, distortion or damage may cause skin injury.
DO NOT APPLY EXCESS PRESSURE OR CAUSE THE MOTOR TO STOP ROTATING. FOR BEST RESULTS, RUN THE UNIT GENTLY OVER THE SKIN. THIS WILL ENSURE SOFT, SAFE AND SMOOTH REMOVAL OF DEAD SKIN.

SERVICE

The service centre address is shown on the attached leaflet.

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