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Massage Gun Pro User manual

Specifications

- A powerful brushless motor
- 3200-3800RPM
- 30 Speed settings
- Weight: 1100 grams
- Noise: 40-60 dB

Vibration frequency
30Hz - 53Hz

Weight
2.5 lbs / 1100 gram

Charge
110 - 240V - 50/60Hz - 5A - 26VDC - 1A

Interchangeable attachments(x7)

Operating instructions

CHARGING

- Fully charge the battery for up to six hours before using or disconnecting the device for the first time.
- To charge, connect the DC end of the supplied 24 V adapter to the battery's charging port and plug in the AC cable in the socket.
- LED band on battery lights pulses to indicate battery level and indicate active charging.
- When the battery is connected, the LED lights on the device will flash to indicate the battery level and indicate active charging.
- Full charge is indicated when all LEDs stop flashing and remain fully lit.
- The battery can be charged at any time and at any battery level.
- It is not recommended to completely discharge the battery to the level of the red LED.
- The average run time is 2+ hours depending on the speed level and the pressure applied during use.

OPERATING

- Align and secure the battery into the base of the device - make sure it is securely locked in place.
- Turn on the device by turning the power switch on the bottom of the battery to the ON position. Turn off the device by turning the power switch to the OFF position.
- With the power switch in the ON position, press the power button once to turn on the device in level one, twice to set the device to level two and three times for level three. If you press the power button again, the device will turn off switched off. The blue LED lights indicate the corresponding power level.
- With the device switched off, fully insert the applicator by pressing it. Remove it by pulling it out.
- Massage the desired body part at the desired speed and apply pressure if necessary.
- To remove the battery, lightly press the battery release button with your thumb. Pull firmly with the other hand at the bottom of the battery/handle.

MAINTENANCE, CLEANING & STORAGE

- Use a slightly damp towel to clean the device and dry it with a soft cloth.
- Remove the battery for storage or travel and keep it in the cardboard box provided.

DISCLAIMER

WARNINGS & CAUTION

- Do not immerse in water. Keep away from liquids or heat sources. Keep ventilation ports free from dust and debris.
- Do not remove screws or attempt to disassemble.
- Do not operate continuously for more than one hour. Allow device to rest for 30 minutes before re-using.
- Unplug the charger from device battery after charging or prior to use.
- For adult use only. Do not use if injured. Consult your doctor before using this product.

SAFETY INSTRUCTIONS

DANGER: TO REDUCE THE RISK OF ELECTRIC SHOCK, FIRE AND PERSONAL INJURY, MUST THIS PRODUCT IS TO BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

- For adult use only.
- Use only on a dry clean surface of the body over clothing by lightly pressing and moving across the body for approximately 60 seconds by regions.
- Use the massage gun only on the body's soft tissue as desired without producing pain or discomfort. Do not use on head or any hard or bony area of the body.
- Use only the applicators that provide the best desired effect.
- Bruising can occur regardless of control setting or pressure applied. Check treated areas frequently and immediately stop at first sign of pain or discomfort.
- Keep fingers, hair or any other part of the body away from the shaft and the back of the applicator, as this could pinch prevent.
- Do not place any objects into the ventilation holes of the Massage Gun
- Do not immerse in water or allow water to enter the ventilation holes of the Massage Gun
- Do not drop or misuse the Massage Gun.
- Only charge with the supplied 26VDC Massage Gun charger.
- Carefully examine the device and battery before each use.
- Do not tamper with or modify the Massage Gun in any way.
- Never let the Massage Gun operate or charge unattended.

DO NOT USE THE Massage Gun OR ANY VIBRATION DEVICE WITHOUT FIRST APPROVAL FROM YOUR DOCTOR. OBTAINED IF ANY OF THE FOLLOWING POINTS APPLY:

Pregnancy, diabetes with complications such as neuropathy or retinal damage, pacemaker wear, recent surgery, epilepsy or migraines, hernias, spondylolisthesis, spondylolysis or spondylosis, recent joint replacements or IUDs, metal pins or plates or other concerns about your physical health. Vulnerable persons and children must be accompanied by an adult when using a vibrator. These contraindications do not mean that you cannot use a vibration or musculoskeletal device, but we recommend that you consult a doctor first.

Research is ongoing into the effects of vibrational exercise for specific medical conditions. This probably will result in a shortening of the list of contraindications as shown above. Practice shows that integration of vibration exercises in a treatment plan is advisable in a number of cases. This should be done on the advice and under the supervision of a doctor, specialist or physiotherapist.

GUARANTEE

IF THIS PRODUCT BECOME DEFECTIVE WITHIN ONE YEAR FROM THE DATE OF PURCHASE BY ANY MATERIAL OR WORKMANSHIP DEFECT, WE WILL REPAIR OR REPLACE THE PRODUCT OR NECESSARY PARTS FREE OF CHARGE.

THIS WARRANTY EXCLUDES:

- damage caused by accident, abuse, mishandling or transportation;
- units subject to unauthorized repair;
- units not used in accordance with instructions;
- damage exceeding the cost of the product.
- deterioration of the delivered product as a result of abnormal storage and/or storage conditions at the customer's premises, and
- failure to present dated proof of purchase.

Some states and countries do not allow limitation of damages, so the foregoing limitation may not apply to you. This warranty guarantees specific legal rights, but other rights may vary from country to country and state to state. This warranty service is available through the operational policy and international distributors.