auronic

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Introduction

- Thank you for choosing our Auronic product! Please always follow the instructions for safe usage. If you have any questions about the product or experience any problems, please contact us by email: service@auronic.nl
- The Auronic Heating Pad quickly heats to 55 - 65 °C and turns off automatically after 90 minutes of continuous use. Choose between 3 settings with the controller with LED display. The Heating Pad is hand washable.

Warnings and Safety Instructions

-  Please ensure that the packaging is kept out of the reach of children! Risk of suffocation!
-  Make sure to keep this user manual. Read the manual carefully before using your Auronic Heating Pad!

General

- Be aware of the risks and consequences resulting from misuse of the Electric Heating Pad. Misuse can result in personal injury, or destruction of the Heating Pad.
- Use the Electric Heating Pad only for purposes described in this user manual.
- If the Electric Heating Pad malfunctions, or appears to be damaged in any manner, do NOT use the Electric Heating Pad and contact Customer Service.
- The Electric Heating Pad is designed and intended for private, non-commercial use only.
- When you are not using the Heating Pad, turn the device off and disconnect the device from the power socket.
- Only setting 1 can be used for continuous operation. ALWAYS slide the slider to setting 1 before going to sleep or when the Heating Pad is warming up.
- Do NOT disconnect the Heating Pad from its power supply by pulling the cable.
- Do NOT carry or pull the device by its cable.
- Make sure that you do NOT trip over the cable. The cable should not be kinked, trapped or twisted.

Persons and Use

- This device should NOT be used by children under the age of 3 or persons with physical, sensory or mental disabilities, or persons with lack of experience with and knowledge of the device, unless they are supervised or have been given instruction concerning the use by a person responsible for their safety.
- Children from the age of 8 should only use the Heating Pad when they are supervised by a person responsible for operating the controller and who is responsible for their safety.
- Keep the Electric Heating Pad away from children under the age of 3 when using and storing. Make sure children do not play with the Heating Pad.
- This electric Heating Pad must not be used by people who are not sensitive to heat or by people who may not be able to react to overheating (diabetics, people with skin conditions or scars, after drinking alcohol, using drugs in intoxicatives, or taking medicine).
- The Heating Pad should only be used on the back and neck. Do NOT place the Heating Pad near your eyes or other sensitive areas. Never place the Heating Pad on or over your head. Do NOT sit on the Heating Pad.
- Prolonged use of the Heating Pad at a high setting could lead to skin burns.
- Stop using the Heating Pad and consult your physician when you experience pain or discomfort when using the Heating Pad.

Medical remarks

- This device is NOT a substitute for treatments administered by doctors or physiotherapists.
- Do NOT use the Electric Heating Pad/consult your physician if you:
 - Are pregnant;
 - Are fitted with a pacemaker or other medical devices.
 - Have circulation disorders, varicose veins, open wounds, bruises, broken skin, inflammation of the veins, or other skin conditions;
 - If you are sensitive to heat.
- Before using the Electric Heating Pad, consult your doctor regarding any unexplained pain, if you are undergoing medical treatment or using medical equipment.

Electronics

- Do NOT reconstruct, repair or disassemble the Electric Heating Pad. This can damage the device. For safety reasons, these parts may only be replaced by authorized service personnel. In order to prevent hazards, always send the Heating Pad to the service centre for repairs.
- Make sure your outlet supplies the correct power voltage to avoid damage to the Electric Heating Pad.
- Always check the power cord and plug before using the Electric Heating Pad. Do NOT operate the Heating Pad if the cord or plug are damaged. A damaged cord or plug must be replaced by a qualified person in order to avoid damages and personal injuries.
- Do NOT use the Electric Heating Pad in or near a bath, shower, swimming pool, or over a basin filled with water when it is turned on. Do NOT operate the device with damp or wet hands and do NOT touch the power cord with damp or wet hands. Make sure to prevent any contact with water. Do NOT submerge the Electric Heating Pad in water. Make sure the controller does not come into contact with water. If water does get inside the casing, turn the device off immediately and contact our customer service.
- If the power fails during usage, turn off the Heating Pad immediately and pull the plug from the socket.
- The Electric Heating Pad heats up during use. Do NOT use the device if you are sensitive to heat.
- Keep the Electric Heating Pad away from moisture, heat and in or close to areas where the surrounding temperature is relatively hot (e.g. stoves or other sources of heat).
- Do NOT cover the controller when it is switched on. Never put Heating Pads or cushions on top of or against the controller. This poses a risk of fire, electric shock or injury.

Other

- Do NOT place a cover over the Heating Pad.
- Do NOT place objects or heat sources (such as hot water bottles) on the electric Heating Pad when it is turned on or stored away.
- The Electric Heating Pad can only be operated using the controller.
- Do NOT use the Heating Pad when it is folded or creased.
- Neither pierce the Heating Pad with, nor attach safety pins or other spiky or sharp objects to the pad.
- Let the Heating Pad cool down completely before folding and storing it. Do not fold it too tight.
- Store the Heating Pad in its original packaging and keep it in a clean, dry place.
- The Heating Pad can be washed by hand. Make sure to ALWAYS disconnect the controller before washing the Heating Pad.

Receiving Your Shipment

- When receiving your shipment, make sure that the contents are consistent with the packing list. Notify your distributor of any missing items. If the equipment appears to be damaged, file a claim immediately with the carrier and notify your distributor at once, giving a detailed description of any damage. Save the damaged packing container to substantiate your claim.

Contents of Package

- Auronic Elektrisch Heating Pad - Model AU3036 - EAN 8720195255536

Included Accessories

- 1x User Manual
- 1x Controller

Product Features

- Technical Data

Size	40x30 cm
Weight	670 gr
Colour	Grey
Material	Fleece
Cable Length	240 cm

Input - Frequency	220V - 240V - 50Hz
Power	100W
3 heat modes	Low: 53±5°C Medium: 59±5°C High: 65±5°C

Other Features

- Automatic shut-off after 90 minutes of continuous use;
- Overheating protection;
- Hand washable;
- Removable controller with LED display.

Operating Instructions

The Heating Pad should only be used on the back and neck. Do NOT place the Heating Pad near our eyes or other sensitive areas. Never place the Heating Pad on or over your head. Do NOT sit on the Heating Pad.

⚠ Prolonged use of the Heating Pad at a high setting could lead to skin burns.

How to Use

1. Connect the plug to a power outlet and slide the slider from 0 to 1. The Heating Pad will start to heat up.
2. Slide the slider to setting 2 or 3 if you prefer a higher temperature. Slide the slider back to a lower setting if the pad becomes too hot.
3. Slide the slider back to position 1 before you go to sleep or when you are using the pad for an extended period of time.
4. Put the slider back on 0 to turn off the pad.

Maintenance and Cleaning

- Remove the plug from its socket and let the Heating Pad cool down for at least 10 minutes. Disconnect the plug from the controller and the heating pad. Hand wash the heating pad. Use lukewarm water and a mild detergent and carefully squeeze the water out. Rinse the heating pad several times to remove all the detergent. Do not tumble dry, or use a (hair) dryer or iron to dry this product. Stretch the heating pad out on a flat surface and dry it naturally. Only use the heating pad once it is completely dry.

Storage

- Let the Heating Pad cool down completely before folding and storing it. Store the Electric heating pad in its original packaging in a clean and dry place. Make sure that you do not damage the power cord. Do not place anything on top of the heating pad. Do not store the device in extreme high or extreme low temperatures.

Disposal - Recycle

- Contribute to a cleaner living environment! This product is subject to the European Directive (EU) 2012/19/EU. Do not dispose of the Heating Pad with the household waste. For proper recycling, please contact your local authority or your household waste disposal service for further details on your nearest designated collection point.

Service and Warranty

- If you need service or information regarding your product, please contact the Auronic Customer Support at service@auronic.nl
- Auronic provides a 2 year warranty on its products. To obtain service during the warranty period, the product needs to be returned with a proof of purchase. Product defects have to be reported within 2 years from the date of purchase. This warranty does not cover: damages caused by misuse or incorrect repair; parts subject to wear; defects that the customer was aware of at the time of purchase; damages or defects caused by customer neglect; damages or defects caused by third parties.

Disclaimer

- Changes reserved; specifications are subject to change without stating grounds for doing so.

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© LifeGoods B.V.