



# Device Description

## 1.1. Instrument description




## 1.2. display interface



## 1.3. Accessories



Plug



Adapter








DC Plug

## While Use




### 1. Selecting the right intensity

- ❖ The device provides 6 different light intensities and advises you on the right setting you should use based on your skin tone. You will always be able to select the light intensity settings as you find comfortable.



- ❖ Press  for 2 seconds to turn the device ON
- ❖ Short press  to adjust the light intensity level, Each time you press  the light intensity level will gradually increase from 1st mode to 6th mode. After 6th mode, press  the light intensity level will gradually decrease until 1st mode.



- ❖ Press and  hold for 2 seconds to switch mode between "Auto" and "Manual"
- ❖ Press  and  hold for 2 seconds to switch mode between "hair removal (HR)" and "skin rejuvenation (SR)"



**Note: There is only one intensity level for skin rejuvenation.**

## While Use

For hair removal, select the desired light intensity settings. Please see the table below to set the light intensity settings.

- ❖ The IPL should never be painful. If you feel uncomfortable, reduce the light intensity settings.
- ❖ After moving to another body area or recent tanning, perform a skin test to determine the right light intensity settings.

Skin Type & Tone	Light intensity Setting
I White: always sunburnt, never tan.	4-6
II Beige: easily sunburnt, tan minimally.	3-4
III Light brown: sunburn easily, tan slowly to light brown.	2-3
IV Middle brown: rarely sunburnt, tan easily.	1-2
V Dark brown: rarely sunburnt, tan very easily.	You cannot use the device
VI Brownish black or darker: rarely or never sunburnt, very dark tanning.	You cannot use the device


### 2. Handling the device



❶. Clean the light exit window before use




❷. Insert the plug



❸. Plug the power supply in the wall socket



❹. Switch on the device and select the right light intensity for your skin tone



❺. Place the device at a 90° angle on the skin so that the skin contact surface of the device.



❻. Leave no gaps between the skin and the device to ensure no light escapes

## Course Of Treatment

### 1st Treatment - (Week - 0)


- ❖ After the first treatment, it can take 1 to 2 weeks for the hair to fall out.
- ❖ During the initial treatments, you might still see some hair growing. These are likely to be the hairs that were in their growing phase during the first treatment.

### 2-3 Treatment - (Week - 2 & Week 4)

- ❖ After 2-3 treatments, you should see a noticeable reduction in hair growth. However, to effectively treat all hair, it is important to keep on treating according to the recommended treatment schedule.

### 4-5 Treatment - (Week - 6 & Week 8)

- ❖ After 4-5 treatments, you should see a significant reduction of hair growth in the areas that you treated with the device.
- ❖ A reduction of hair density should also be visible. Keep on treating with frequent touch-ups (every 4 weeks) to maintain the result.



## Before Use

### 1. Get your skin ready

- ❖ Before you use the IPL device, remove all visible hair in the area you wish to treat.
- ❖ Ensure there are no hair remaining above the surface of the skin. This allows the light to be easily absorbed by the melanin located in the hair roots and ensures effective treatment.
- ❖ You should only shave before using the IPL.
- ❖ Do not wax, epilate or use depilatory creams, as anything else pulls hair from the roots (follicle) and hinders the IPL process, which you don't want.
- ❖ If you choose to wax, please wait 6-8 days before using the device to let your skin rest.
- ❖ We recommend that you take a shower before the treatment to ensure the skin is free from oily substance, dirt and then pat your skin dry.

### 2. Skin test

- ❖ Check the skin tone table to see if your skin type allows the usage of the device.

### 3. Hair test

- ❖ Check the hair color table to see if your hair color allows the usage of the device.

