











Commentary Mode




Features

	IPL rejuvenation		Sonic vibration		Deep Cleaning		Nourish import
	Positive and negative ions		Four-speed adjustment		EMS Microcurrent		Constant temperature hot compress


Packing List



Neck beauty device



Data line



Manual

! Q&A Skincare Small Class

Q: Can it be used every day? How long does it last?
 A: It can be used every day, it is recommended to use it once in the morning and once in the evening for 5-10 minutes each time.

Q: What should be paid attention to before use?
 A: Clean the skin before use. It is better to use with neck cream and massage cream. When using, gently pull and lift each part, do not use too much force, so as not to cause damage to the skin.

Q: Can it be used on acne-prone skin?
 A: It is recommended to use it with caution for acne-prone skin to avoid secondary damage.

Q: How to charge and use?
 A: This product uses a USB charging port to store electricity, which is convenient and fast. Generally, it takes 3 hours to charge, and the machine is fully charged.