



Attention

In order to reduce the risk of failure, electric shock, injury, fire, death and equipment damage, please carefully read all instructions and warning tips in this manual before using this product.

Safety instructions

1. Minors are not recommended to use this product. If you need to use it, please operate it under adult supervision.
2. The massage time of the same part shall not exceed 60 seconds. For sensitive skin, please use it after covering towel or wearing clean and dry clothes
3. It can only be used at high speed and high pressure at the soft tissue part of the body. Please stop using when you feel pain or discomfort.
4. Please do not use it at high speed high pressure on the head and any bone parts of body.
5. Please don't cover the vent of the product.
6. Please don't immerse the product in water, and also do not let water enter the vent of the product.
7. Please don't put the product into the fire or near the fire source.
8. It is forbidden to modify the equipment without permission.
9. Please don't run or charge the product when nobody is around.
10. In case of physical injury, please consult your doctor before use.
11. Please use our original battery and charger. Please charge for 6 hours before using this product for the first time.
12. Keep fingers, hair or other parts of the body away from the massage shaft and cooling holes because of the risk of being caught.

03

Function description

Through strong and powerful high-frequency vibration stimulation, the product can penetrate into the texture, penetrate deep fascia, rapidly disintegrate the accumulation of lactic acid generated after exercise, relax the tight and stiff muscle tissue, comb the muscle nodules, prevent muscle adhesion, release lactic acid completely, and relieve fatigue and discomfort.

This product is only used for physical therapy and health care, not use as medical appliances.

Prohibition

Please don't use this product without the permission of your doctor under any of the following situation:
Pregnancy, diabetic complications (such as neuropathy or retinal damage), pacemaker wearing, recent surgery, epilepsy or migraine, disc herniation, spinal forward, vertebral detachment, or spinal joint stiffness, recent joint replacement surgery or IUD, metal needle or steel plate and any other issues related to your health.

04

Thank you for choosing this product. Before you use it, please read the user manual in detail so as to master the correct use of this product. When referring, please pay special attention to the precautions and keep the operation manual properly for reference.

Note: The company reserves the right to change and explain the design of the product. If there is any change, we will not give any notice. The color of the product is subject to the material object.

Catalog

1. Directions for using 01
2. Buttons and display description..... 02
3. Attention 03
4. Function description..... 04
5. Maintenance and storage..... 05
6. Product parameters..... 05
7. Troubleshooting..... 06
8. Name and content of hazardous substances 07
9. Quality assurance..... 08

Maintenance and storage

1. Wipe the surface of the equipment with a slightly damp towel, then dry it with a dry towel and place it in a dry and ventilated place.
2. Please do not put it in the place where high temperature, humidity or direct sunlight.
3. Please turn off the product for 30 minutes after 1 hour using.
4. When the storage period is more than 6 months, please power on the product for inspection

Product parameters

| | | |
|--------------------|--------------------------|---------------|
| Product name | Massage gun | |
| Product model | KX-2212 | KX-2212D |
| Interface type | DC 16.8V ~ 0.6A | DC 12.6V ~ 1A |
| Rated power | 30W | 25W |
| Battery capacity | 2200mAh | 1800mAh |
| Use time | 6 hours | 4 hours |
| Product weight | 680g | 620g |
| input voltage | 100~240V~50/60Hz | |
| Decibel | ≤55 dB | |
| Product size | 24X18.5 X 6.5cm | |
| Protection time | 15 minutes | |
| Use environment | 0~40°C , Non- waterproof | |
| Executive standard | GB31241-2014 | |

05

Troubleshooting

If there is a fault during the operation of this product, please check it. If the fault cannot be eliminated or the product still cannot be operated, please turn off the power and contact the service center of this product for maintenance.

| Possible causes | Possible reasons | Troubleshooting |
|--|---|--|
| Product cannot work | Check the remaining power. Is the power weak? | Charging |
| Product cannot be charged | Is the power charger connected properly? | Pull out and insert again |
| | Is the charger and socket in good condition? | Check the condition of the socket. If there is any damage, please contact the product service center |
| | Does the adapter meet the requirements? | Use the charger specified by the manufacturer |
| The product can be turned on, but the function cannot work | Check the remaining power. Is the power weak? | Charging |

Please do not disassemble or try to repair by yourself. Once disassembled, we will not be responsible for it.

06

Directions for using

Using direction

1. Power button (speed control button) is located in the middle depression of the control panel. Press for 5 seconds to turn on or off
2. In the start state, it is the first gear. Press the power button continuously, and the speed will cycle from the first gear to the fifth gear and then to the first gear
3. When there is only one power indicator and it is flashing, you need to charge the machine. When it is fully charged, four power indicators will be on. Charging takes 2-3 hours and using time is 4-6 hours
4. In order to avoid muscle and bone damage caused by too long using, the machine will start the intelligent protection function and automatically turn off when it works continuously for 15 minutes, please restart the power button if you need to use it again.
5. The charging port is at the bottom of the machine. If use charging cradle (Selectable), just insert the product into the charging cradle.

Big ball head: Apply to arm, waist, buttock, thigh, calf and other large muscle groups
Air plug head: Suitable for the whole body bones and muscles in sensitive parts
U head: Suitable for neck, spine and Achilles tendon
Bullet head: Applicable to impact deep tissues, such as meridians, palms, soles, etc
Flat head: Suitable for relaxing and shaping all parts of the body

01

Buttons and display description

02

Name and content of hazardous substances

Statement of toxic and harmful substances in products

| Part name | Toxic and harmful substances or elements | | | | | |
|-----------|--|--------------|--------------|------------------------------|-------------------------------|--------------------------------------|
| | Lead (Pb) | Mercury (Hg) | Cadmium (Cd) | Hexavalent chromium (Cr(VI)) | Polybrominated biphenyl (PBB) | Polybrominated diphenyl ether (PBDE) |
| Massager | ○ | ○ | ○ | ○ | ○ | ○ |
| Data line | ○ | ○ | ○ | ○ | ○ | ○ |

This form is prepared in accordance with SJ / T11364-2014
○: It means that the content of the hazardous substance in all homogeneous materials of the part is below the limit requirements specified in GB / T 26572-2011
●: It means that the content of the hazardous substance in at least one homogeneous material of the component exceeds the limit requirements specified in GB / T 26572-2011

This product does not contain harmful substances. It is a green and environmental protection product. It can be recycled after being discarded and should not be discarded at will

07

Warranty card

| | |
|--------------|----------------|
| Product name | Model |
| User name | Contact number |
| Address | |
| Post code | Purchase date |

Distributor

| | |
|-----------|----------------|
| Address | |
| Post code | Contact number |

Customer service instructions

1. Users can enjoy one year free warranty from the date of purchase (appearance is not covered by the warranty), and the purchase time shall be subject to the valid purchase documents issued.
2. Repair with this card and valid purchase documents (please send the machine, this card and valid purchase documents to the repair point together when the foreign customers need to repair).
3. The following situations are not within the scope of maintenance, only the maintenance service is provided, but the labor and cost are charged:
(1) Damage caused by all human factors, including damage caused by use under abnormal working environment and improper use and storage not in accordance with the instructions.
(2) The customer disassembles the machine, repairs or refits it by himself or repairs it by an unauthorized unit of the company.
(3) There is no warranty card and valid documents.
(4) Damage caused by self poor transportation or other irresistible factors (such as flood, fire, lightning, earthquake, etc.) after purchase
4. The company reserves the right of modification and final interpretation of the above terms.

Quality assurance

08