

Regrow^{MD}

MEDICAL DEVICE FOR HAIR GROWTH

How the RegrowMD Laser 272 Device Works

The RegrowMD Laser 272 Device works through a process of photo-biostimulation, stimulating the hair follicle at a cellular level with the precise wavelength of laser light energy for hair growth.

This safe and effective treatment, known as low level laser therapy (LLLT), revitalizes and enlivens weakened hair follicles, restores your natural hair growth cycle, reverses the thinning process and stimulates hair growth.

CONTACT US AT LEXINGTON INTERNATIONAL, LLC:
Corporate Office – Customer Service – Warranty Repairs
 1040 Holland Drive, Boca Raton, Florida 33487 USA
 info@regrowmd.com

Direct line 1-561-417-0200
 United States and Canada 1-866-527-3726
 United Kingdom 0808-134-9868
 Australia +61 7 3198 2484

**Authorized Representative in Europe
 (Regulatory Affairs Only):**

EC	REP	EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands Tel: +31 (0) 70 345-8570 Fax: +31 70 346-7299
----	-----	---

Australian Sponsor:
 Emergo Australia,
 Level 20, Tower II Darling Park,
 201 Sussex Street Sydney,
 NSW 2000 Australia

Storage and transportation temperature/humidity: -20°C to 60°C / 0% to 90%, relative humidity.
 Made in China
 Patent www.HairMax.com/patents
 © 2020 Lexington Intl., LLC. All Rights Reserved



Scan to learn more



Laser Light - Avoid direct eye exposure.



LASER272

HAIR GROWTH LASER LIGHT DEVICE

BIO LIGHT COMFORT
 DESIGN

- This Package Includes:
- Laser 272 Hair Growth Device
 - Sports Cap
 - Battery Pack
 - Belt Clip for Battery Pack
 - Universal DC Adapter
 - User Manual

Silicone Insert for a Comfortable Fit

Lightweight Sports Cap

BioLight Comfort Design

272 High Quality Lasers



LIT-SLVRGMDL272 02 1220



HAIRMAX
 The science of hair growth.
 HAIR GROWTH LASER LIGHT DEVICE

CAT 856823006619 2023-07-27 Lexington International, LLC
 SN 2023031021 1040 Holland Dr. Boca Raton, FL 33487
 (01) 0 0856823 00661 9 (11) 230727 800 973 4769 (USA) 561 417 0200 (All Others) www.hairmax.com
 (21) 2023031021 QTY: 1 EA

