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### IMPORTANT SAFEGUARDS

The Power toothbrush is safe and effective when used according to directions. Use all electrical appliances, however, basic safety precautions should be followed, especially when children are present.

#### READ ALL INSTRUCTIONS BEFORE USING

#### DANGERS-WATER

To reduce the risk of burns, electrocution, fire, or injury to persons:

- Do not place or drop the charger base where it can fall or be pulled into a bathtub, toilet, or sink.
- Do not place or drop the charger into water or other liquid.
- Do not use the product that has been taken into water or other liquid, including immediately.
- Do not use the battery.
- Unplug charger base before cleaning and assure that it is dry before plugging into an electrical outlet.

#### DANGERS-VOLTAGE

To reduce the risk of burns, electrocution, fire, or injury to persons:

- Do not plug the charger base into an outlet with a voltage other than specified on the label of the charger.
- The product is designed to operate within a range of 100-240 volts.
- Voltage conversion is not guaranteed.
- Do not fire the plug into an outlet.
- The plug does not easily fit into the outlet, disconnect use.

#### DANGERS-WARNINGS

To reduce the risk of burns, electrocution, fire, or injury to persons:

- Do not use the product if it is not working properly, if it has been dropped or damaged, or if it has been exposed to water, immediately disconnect use of the product against damage in any way.
- Use the toothbrush, charge base and power cord away from heated surfaces.
- Use the toothbrush, charge base and power cord away from heated surfaces.
- Use the toothbrush, charge base and power cord away from heated surfaces.

### CONFIGURATION

Professional brush head

Color ring

Power button

Mode button

Charging indicator

Charger base

Power cord

Power plug

### SAVE THESE INSTRUCTIONS

### GETTING STARTED

#### Charging the toothbrush

Be sure to turn off the toothbrush before charging. Plug the power plug into the outlet. Then insert the power plug into the charging. The charging indicator light will be on, indicating the toothbrush is charging. After the battery is fully charged, the charging indicator light will turn off, indicating the toothbrush is fully charged.

#### Notes:

- It takes at least 24 hours to charge the battery.
- The toothbrush must be charged for 24 hours, if this is the first time using the toothbrush.
- The toothbrush has not been used for more than three months.
- To keep the battery fully charged, we recommend that you keep the toothbrush in the charger when not in use.
- During charging, the toothbrush will get warm, but this is normal.
- The operating indicator will become red in the battery gate color. If the toothbrush operates only for a short time after being fully charged, this means that the battery is nearing the end of its life cycle. Contact service center for support.

#### HOW TO USE

- Slide the brush head onto the handle. Push brush head down until it locks into place on handle.
- There is a small gap between color ring and the handle (Fig. 1).
- Use the toothbrush and apply a small amount of toothpaste.
- Apply toothpaste to the brush head (Fig. 2).
- Place the brush head against the handle at an angle.
- Apply toothpaste to the brush head (Fig. 2).
- Place the brush head against the handle at an angle.
- Apply toothpaste to the brush head (Fig. 2).
- Place the brush head against the handle at an angle.
- Apply toothpaste to the brush head (Fig. 2).

### TRUBLESHOOTING

Problem	Possible cause	Solution
The toothbrush does not work.	You have not purchased the correct model for your country.	Change the main unit for the correct model.
The toothbrush does not work.	The toothbrush is pressed together against your teeth.	Place the brush head against your teeth.
The toothbrush does not work.	The battery has reached its end of life.	Contact service center for support.
The toothbrush does not work.	Changing the line voltage.	Change the main unit for the correct model.

### WARRANTY

Your toothbrush is guaranteed for 1 year against manufacturer's defects. The warranty is valid for the original purchaser only. The warranty is valid for the original purchaser only. The warranty is valid for the original purchaser only.

### HOW TO USE

Apply light pressure and gently move the brush head across the teeth in a semi-circle.

Apply light pressure and gently move the brush head across the teeth in a semi-circle.

Apply light pressure and gently move the brush head across the teeth in a semi-circle.

### MAINTENANCE

Remove the brush head.

Wipe all surface areas with a damp cloth.

Wipe all surface areas with a damp cloth.

Wipe all surface areas with a damp cloth.

### RECHARGEABLE

Using a soft cloth, wipe off any dirt on the charge and the settings.

Using a soft cloth, wipe off any dirt on the charge and the settings.

Using a soft cloth, wipe off any dirt on the charge and the settings.

### RECYCLING

The entire surface of this product is designed to be recycled. All parts of the product are made from recycled materials.

The entire surface of this product is designed to be recycled. All parts of the product are made from recycled materials.

The entire surface of this product is designed to be recycled. All parts of the product are made from recycled materials.

### BINTOI

Remove the brush head from the handle. Push brush head down until it locks into place on handle.

Remove the brush head from the handle. Push brush head down until it locks into place on handle.

Remove the brush head from the handle. Push brush head down until it locks into place on handle.

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© 2023 Philips. Alle rechten voorbehouden. Philips is een handelsmerk van Philips. Philips is een handelsmerk van Philips. Philips is een handelsmerk van Philips.