

# U-shaped electric toothbrush

Rated voltage: 5V

Body: ABS

Charging time: <120Minute

Rated power: <5W

Brush head: Food grade silicone




Executive standard: GB4706.1-2005/GB4706.59-2008



We recommend to replace the brush head within 3 months

## USE METHOD:



STEP1



STEP2



STEP3

Apply toothpaste on the two  
sides brush a head to  
squeeze toothpaste

Put it in your  
mouth and press  
the switch

Shake the toothbrush  
from side to side to  
clean your teeth

FDA CE FC RoHS

MADE IN CHINA