

Infrared Light Therapy Lamp for Body



✓ Anti-aging



✓ Reduce wrinkles




✓ Improve skin elasticity



✓ Reduce Inflammation



✓ Muscle Recovery



✓ Pain Relief



✓ Wound Healing



✓ Better Sleep

User Instructions :

1. Keep yourself 6-12 inches from the device while exposing your skin for the best results
2. Apply the red light panel for 10-30 minutes of continuous dosing.
3. Use your red light panel device 3-7 times a week depending on your tolerance.

Contact Us : Rxc-service@outlook.com

Made in China
 FC CE FDA RoHS