FOR A GOOD REASON

GRUNDIG

**INSTRUCTION MANUAL
BEDIENUNGSANLEITUNG
MANUEL D'UTILISATION
GEBRUIKSAANWIJZING
MANUALE DI ISTRUZIONI
MANUAL DE INSTRUCCIONES**

UK: Digital body scale

DE: Digitale Körperwaage

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NL: Digitale personenweegschaal

IT: Bilancia digitale pesapersona

ES: Báscula digital

ITEM: 16426

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[UK] INSTRUCTION MANUAL

Dear User,

Thank you for choosing our body scales.

For the purpose of maintaining product longevity, weighing accuracy and correct use of the body scales, please read the following instructions carefully before use.

SPECIFICATIONS:

Auto off: power off automatically after no operation for approximate eight seconds.

Maximum weight capacity: 180 kg

Minimum measurable weight: 2 kg

Graduation: 50 g

User-selectable units: kg / lb

Low-battery indication: LO

Overload indicator: oL

Error display: Err

Operating temperature: 5 °C – 35 °C

Power supply: 2xAAA battery

INSTRUCTIONS FOR USE

Depending on the particular type of your scales:

1. Place 2 x AAA battery in the battery compartment in the back of the scales.
2. Choose your preferred weight units (kg / lb) by pressing the unit conversion button.

At the bottom of the scale you can find a unit conversion button as shown in this image

3. When weighing, place the scales on a hard, flat surface. Do not use the scales on a carpet or any other soft surface. When measuring, remove your shoes, socks etc. to ensure that the measurement is as accurate as possible. Remove any items from your person that could affect the weight. Make sure the bottoms of your feet are clean

before you step onto the scales barefoot.

4. Turn the scales on by standing on them. The display will show 0.0 kg. Now you can start weighing yourself.
5. Always stand in the centre of the platform, distributing your weight evenly.
6. Within seconds, your weight will be established and displayed.
7. A few seconds after you have stepped off the scales, they will switch themselves off automatically .

TROUBLESHOOTING:

1. Display shows 'oL'.

If the weight exceeds the maximum capacity of 180 kg, 'oL' appears on the display. Step off the scales immediately, otherwise they may be damaged.

2. Display shows 'Lo'.

Low battery power. Replace the batteries.

PRECAUTIONS:

1. Avoid overloading the scales, otherwise permanent damage will occur.
2. Do not disassemble the product. Only replacing the battery is allowed.
3. Clean the scales after use with a slightly damp cloth. Do not use solvents or immerse the product in water.
4. Do not inflict any excessive impact or vibration on the scales, for example by dropping them on the floor.

NOTICE:

1. The platform will be slippery if it becomes wet. Keep it dry!
2. Do not place any objects on the glass platform while you are not using the scales.

IMPORTANT: Remove the plastic packaging under the battery compartment before use

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Ⓢ The symbol above and on the product means that the product is classed as Electrical or Electronic equipment and should not be disposed with other household or commercial waste at the end of its useful life.

The Waste of Electrical and Electronic Equipment (WEEE) Directive has been put in place to recycle products using best available recovery and recycling techniques to minimise the impact on the environment, treat any hazardous substances and avoid the increasing landfill. Contact local authorities for information on the correct disposal of Electrical or Electronic equipment.

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