


Is the HairMax right for your Hair Loss?

HairMax can be used by men and women to help grow hair and treat hair loss. Doctors use the following scales to describe different hair loss and skin types. The shaded areas show the type of hair loss and skin types that can be treated with Laser Light Therapy. The HairMax has not been clinically tested for hair loss or skin types outside the shaded areas.

Females



Males



Fitzpatrick Skin Type Classification Guide:

A way to describe your skin type based on your complexion and what happens to your skin in sunlight.

Clinically Tested for these skin types:

- Type 1 Very white or freckled skin, always burns with sun exposure (very fair; often in people with red or blond hair and blue eyes)
- Type 2 White skin, usually burns with sun exposure (fair; often in people with red or blond hair and blue, green or hazel eyes)
- Type 3 White or olive skin tone, sometimes burns with sun exposure (fair; seen in people with any hair or eye color)
- Type 4 Brown skin, rarely burns with sun exposure (common in people of Mediterranean descent)

Not Clinically Tested for these skin types:

- Type 5 Dark brown skin, very rarely burns with sun exposure (common in people of Middle-Eastern descent)
- Type 6 Black skin, rarely burns with sun exposure

Chart adapted from Fitzpatrick's Dermatology in General Medicine. McGraw-Hill Professional, 5th edition

*HairMax laser devices are indicated to treat Androgenetic Alopecia, and promote hair growth in females who have Ludwig (Savin) I-4, II-1, II-2, or frontal patterns of hair loss and in males who have Norwood Hamilton Classifications of IIIa to V and who both have Fitzpatrick Skin Types I to IV.

©Copyright 2019 Lexington Intl., LLC. All Rights Reserved

Storage and transportation temperature/humidity: -20°C to 60°C / 0% to 90%, relative humidity

Made in Taiwan

Patent www.HairMax.com/patents

