# ACHATÉ ANTI AGING

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## **INTRODUCTION**

---

Thank you for purchasing one of our Achaté products. We recommend that you read the user manual carefully before using the product. Always follow the safety instructions for safe use.

Do you have questions about the product? Please feel free to contact us by email.

The Achaté Light Therapy Facial Mask is suitable for every skin type thanks to the 7 different LED functions. This device treats acne, scars, wrinkles, pigmentation spots and other imperfections.

## **HOW TO USE**

---

Step 1. Connect the face mask to a power network.

Step 2. Press the on/off button.

Step 3. Press the fan button (heat dissipation)

Step 4. Adjust the time to your own needs.

Step 5. Choose the type of treatment you want. Keep pressing the yellow button until you get the position you want. You can also switch on double light.

Step 6. Turn off the device after use and unplug the face mask from the power.

## **WARNINGS AND SAFETY INSTRUCTIONS**

- Read the instructions below thoroughly before using the Achaté LED Light Therapy Facial Mask.
- Follow the steps for using the device at all times.
- Following the steps is important because of the risks and consequences associated with misuse of the device.
- Misuse of the device may result in personal injury or damage.
- This appliance is not suitable for people under the age of 16, people with physical or mental disabilities. Is the user under 16 years old? Then make sure that there is a supervisor to guarantee the safety of the user.
- Only use the face mask for what the product is intended for.
- Make sure that the appliance is always connected to the correct voltage.
- If the device is connected to a power network, make sure that you do not leave the device unattended.
- Regularly check that the plug and cord are not worn. Is there damage? In that case, stop using the appliance and have the cord and plug replaced.

- In all cases of damage, you can no longer use the device.

- You should never repair the device yourself, this is due to the risk of personal injury or damage.

- This product is not water resistant, so keep it away from damp places or water.

Should the device have come into contact with water;

1. Remove the plug from the socket immediately.
2. Do not touch the device again.
3. Have the product checked before using it.

- If you are going to clean the product, please disconnect it from the power supply first.

- Only use the device indoors.

- Only use the device for personal use.

- Keep the face mask out of the reach of children.

## **LIGHT THERAPY MODES**

**Red light:** Stimulates blood circulation, accelerates production of collagen and makes skin more elastic, smooth and even

**Green light:** Face for combination skin, heals wounds, reduces skin wrinkles and improves pigmentation (mainly anti-aging)

**Blue light:** Soothes sensitive skin, tightens skin and breaks down bacteria

**Yellow light:** Promotes lymphatic drainage, smooths skin, reduces redness and improves oxygen in cells

**Purple light:** Improves cell renewal, collagen production and cleanses the skin

**Blue & Green light:** Oil secretion is reduced, less shiny skin and accelerates metabolism

**White light:** Penetrates deepest into skin; accelerates skin metabolism, tightens skin, reduces fine lines and wrinkles, reduces hyper-pigmentation and reduces inflammation & scarring

## **CONTENT PACKAGING**

- 1x Achaté LED Light Therapy Face Mask
- 1x Manual English

### **Product Features**

## **MAINTENANCE AND CLEANING**

The LED Light Therapy Face Mask from Achaté is almost completely maintenance-free. If you want to clean the mask, first remove the plug from the socket. Clean the device with a slightly damp cloth.

## **STORAGE**

When you are not using the device, we recommend that you store the face mask and unplug it from the mains. Keep it in a dry and clean place. Make sure that the device is not stored in extremely high or low temperatures.

## **ENVIRONMENT**

When the device has reached the end of its life, dispose of it at the designated collection points.

## **SERVICE AND WARRANTY**

If you have any further questions or would like information about the product, please do not hesitate and contact us directly by email; [info@achate.nl](mailto:info@achate.nl)  
All Achaté products come with a 2-year warranty.

Warranty does not apply when:

- Damage that has been caused by yourself
- Improper use or repair
- Wear on sensitive parts
- Due to damage from third parties
- Defects due to customer's own fault
- In case of defects that the customer is aware of

# **ACHATÉ ANTI AGING**