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
## NON-CONTACT INFRARED THERMOMETER

*Dear customer,*

*You have just bought this product and we thank you for your purchase. We hope you enjoy using the product. In order to be completely satisfied, we recommend that you read these directions carefully.*

### MAINTENANCE

- Clean the thermometer with a soft, dry cloth.
- Clean both lenses with a cotton swab moistened with alcohol and immediately wipe dry with a clean cotton swab. After cleaning, allow at least 20 minutes drying time before taking temperatures.
- Do not touch either lens with your fingers.
- Do not submerge the thermometer in water or any other liquid.
- Do not rinse the thermometer in water.
- Do not use water directly on the LCD screen.
- Store the thermometer in the dry, clean place.
- If the thermometer is not used for a long time, remove the batteries.

### SAFETY

- Keep out of reach of children and to avoid inhalation or swallowing of small parts. Do not allow children to take their temperatures unattended.
- Do not use this appliance for any other purposes than those for which it was designed (human body temperature).
- Use of this thermometer is for monitoring purposes only and is not intended as a diagnostic tool. Make sure that any unusual readings are communicated to your physician. Under no circumstances should you change your prescribed medicine doses.
- Use of this thermometer is not intended as a substitution for consultation with your physician. Measurement results are for reference only. Contact your physician if you have or suspect any health.

### RECOMMENDATION

- Store and use the thermometer in the same room.
- Do not use in areas where it may be subject to room temperature variations.
- Determine your normal temperature range by taking your temperatures when healthy. These readings can be used as a form of reference to compare with future readings when ill.
- Before taking your temperature, make sure your forehead is free of sweat and hair, and make sure the infrared sensor is clean.

## NON-CONTACT INFRARED THERMOMETER

- Do not take temperature within 30 minutes of eating, exercising, or taking a bath.
- Wait one minute between each temperature reading. Two consecutive readings may give slightly different measurements. Should this occur, take the average of the two.
- Do not attempt to disassemble the thermometer.
- Keep thermometer away from direct sunlight and store well away from any source of heat in a clean, dry location.
- Do not drop or break the thermometer. Do not use if broken. Do not disassemble.
- Don't use a mobile phone nearby when the thermometer is working or do not use in an area with possible electromagnetic interferences.
- If these recommendations are not respected, the readings may be inaccurate.
- Temperature will vary depending on the part of the body. This is why you should not compare rectal and forehead temperature readings, for example.
- Temperature will also vary depending on the person, age, and time of day.
- Don't store the thermometer to extremes temperature or extremes humidity. Failing to do so may cause inaccuracy. Don't operate this thermometer under extremes temperature or extremes humidity. Failing to do so may cause inaccurate measurement.

### BODY TEMPERATURE

Normal body temperature is a range. A person's normal temperature range tends to decrease with age. The following table shows normal temperature ranges by age:



<b>0 – 2 years</b>	36.4 – 38 °C	97.5 – 100.4 °F
<b>3 – 10 years</b>	36.1 – 37.8 °C	97.0 – 100.0 °F
<b>11 – 65 years</b>	35.9 – 37.6 °C	96.6 – 99.7 °F
<b>&gt; 65 years</b>	35.8 – 37.5 °C	96.4 – 99.5 °F

The range of normal temperature varies from person to person and can be influenced by many factors such as time of day, measure from different sites, level of activity, medications, emotion and so on. So we recommend that you practice with the Thermometer on yourself and family members when you are healthy. This way you know how the thermometer works and can feel more confident of the measurements you take when a family member is ill.


## NON-CONTACT INFRARED THERMOMETER

## PRODUCT DESCRIPTION

## INSERTING &amp; CHANGING BATTERIES

- To insert or change the 2 LR3 batteries (supplied with the thermometer), open the battery compartment by sliding off its cover, insert the batteries ensuring that the negative and positive ends are in the right place, and replace the cover before use.
- When the low battery indicator  appears on the LCD, or the thermometer does not function at all, you should replace new battery as soon as possible.

## CHOOSING THE TEMPERATURE MEASURING UNIT (C° OR F°)

- Press one of both buttons to turn on thermometer. Then, to switch from one measuring unit to the other, press around 5 seconds [M].

## NON-CONTACT INFRARED THERMOMETER


## USING THE PRODUCT

## FOREHEAD TEMPERATURE

Place the thermometer in your hand as shown with your thumb on the [measure] button.

- Press the [measure] button  to turn on the thermometer.
- Place the thermometer from 0 to 5 cm of the forehead and press the [measure]  button again. A double beep sound indicates that the distance is correct.
- If the position is OK, release the button, but keeping the thermometer flat until you hear another double beep sound indicating that the temperature has been taken. The screen lights up displaying temperature. The measurement is complete and the thermometer may be removed from the head.
- You can read the temperature on the display and the color indicator (red or green) around the [measure] button helps to analyze the temperature, if fever or not.

Green: 34.0°C (93.2°F) ~ 37.9°C (100.3°F)

Red: 38.0°C (100.4°F) ~ 42.9°C (109.3°F)

- Temperature is automatically memorised.

**CAUTION:**

After pressing the temperature measurement button, the distance sensor will estimate the distance between the forehead and the thermometer.

If the distance between the forehead and the thermometer is adequate (between 0 and 5 cm) you will hear a first double beep while your temperature is displayed followed immediately by a 2nd double beep.

If the thermometer is too far from the forehead (more than 5 cm) short beeps will follow one another indicating that the thermometer must be brought closer to the face, until you hear the distance validation double beep.

However, if the distance is still not adequate for body temperature measurement, the thermometer will display the room temperature.

## NON-CONTACT INFRARED THERMOMETER

**CAUTION:**

External factors may influence forehead temperature particularly when an individual has:

- been exposed to very hot or very cold temperatures.
- been recently swimming or bathing.
- had their forehead covered.

In these cases, let the individual stay away from the above scenario and station in a room for 30 minutes prior to taking a temperature.

**OBJECTS OR LIQUID TEMPERATURE**

- To measure the temperature of the surface of an object or liquid, put the thermometer in object mode.
- Press the [Measure] and [Mem] buttons together in order to display like right screen.

**CAUTION:**

Make sure don't use this mode to measure body temperature.

- Hold the thermometer and ensure the probe is around 0-5 cm from the liquid or surface whose temperature you want to measure. Press the [measure] button  to start the measurement.
- A beep sounds signals the end of measurement. Remove the thermometer from the object. The LCD displays the measured temperature.

**CAUTION:**

With object measurement, the current surface temperature of the object is displayed. It can be different from its internal temperature, especially if the surface is exposed to direct sunlight or a draught.

- The thermometer will automatically return back to standby mode and to body temperature mode after 1 minute.

## NON-CONTACT INFRARED THERMOMETER

**CHECK THE MEMORY**

- There are a total of 9 memory data for recording forehead temperature measurements.
- The current measurement is always stored in the last storage space. When all storage spaces have been occupied, the oldest measurement is always deleted from the memory.
- Press [M] button to turn on thermometer and [M] button again to enter in memory data .

The LCD displays memory number followed by measured temperature like figures. Press the [MEM] button repeatedly to go to the memory number you need (9 is the youngest, 1 the oldest).

**PROTECT THE ENVIRONMENT**



Dispose of used batteries in a designed container so that they can be collected and recycled. Do not mix different types of batteries. Do not mix new and used batteries. After life, dispose of this product at a designated waste recycling point.

**GUARANTEE**

This product is guaranteed against defects in materials or manufacturing. During the guarantee period, any such defects will be replaced free of charge (proof of purchase must be shown in the case of a complaint under guarantee). The guarantee does not cover damage resulting from improper use or negligence. In the event of a claim, first contact the shop where you purchased the product.

## NON-CONTACT INFRARED THERMOMETER

## TROUBLESHOOTINGS

Troubles	Checklists	Counter measures
	<p>The measured temperature is lower than 32.0°C/89.6°F or higher than 42.9°C/109.3°F Please check the operation method</p>	<p>Follow User's Manual for proper measurement.</p>
	<p>Hardware problem</p>	<p>Contact your distributor</p>
	<p>Operating temperature is out of the range</p>	<p>Use the thermometer in the range of operating conditions.</p>
<p>The measurement is not accurate or if there is any doubt on the measured result</p>	<p>Please check if the infrared sensor is clean or not</p>	<p>Clean the infrared sensor with cotton swab according to user manual</p>
	<p>Please check the measuring way is correct or not</p>	<p>Ensure you have read the manual and know how to use the thermometer properly. Please keep the thermometer and patient in the measuring room at least 30 minutes before using Please take the measurement indoor Do not held the thermometer in your hand too long. It could affect the accuracy</p>
	<p>Please check if there is low battery symbol on LCD</p>	<p>Change new battery</p>

## NON-CONTACT INFRARED THERMOMETER

## ELECTROMAGNETIC COMPATIBILITY(EMC)

- This device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used.
- Use of accessories other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- **When the operating environment is relatively dry, strong electromagnetic interference usually occurs. At this time, the device may be affected as follows:**
  - the device stops output;
  - the device turns off;
  - the device restarts;

**The above phenomenon does not affect the basic safety and essential performance of the device, and the user can use it according to the instruction. If you want to avoid the above phenomenon, please use it according to the environment specified in the manual.**

Table 1

Declaration - electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The thermometer use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The thermometer are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Not applicable	

## NON-CONTACT INFRARED THERMOMETER

Table 2

Declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315°  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°  0 % UT; 250/300 cycles	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

## NON-CONTACT INFRARED THERMOMETER

Table 3

Declaration - electromagnetic immunity			
The thermometer are intended for use in the electromagnetic environment specified below. The customers or the users of these thermometer should assure that it is used in such environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3V 0.15 MHz to 80MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of device, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ 150 KHz to 80 MHz $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 80 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10 V/m	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range.  <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 

## NON-CONTACT INFRARED THERMOMETER

Note1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note2: These guidelines may not apply in all situations.

Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

<sup>a</sup> Field strengths from fixed RF transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which device is used exceeds the applicable RF compliance level above, device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating device.

<sup>b</sup> Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4

<b>Recommended separation distances between portable and mobile RF communications equipment and device</b>			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and device, as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer			
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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## SPECIFICATIONS

Power supply:	DC 3V (2 x AAA battery)
Mode of operation:	Adjusted mode
Reference body site:	Core
Body measuring range:	32.0°C ~ 42.9°C (89.6°F ~ 109.3°F)
Laboratory accuracy:	±0.2°C(±0.4°F) for 35.0°C (95.0°F) ~ 42.0°C (107.6°F) ±0.3°C (±0.5°F) for other range
Object measuring range:	-22°C ~ 80.0°C (-7.6°F ~ 176.0°F)
Object measuring accuracy:	±2.0°C(±4.0°F)
Displayed room temperature range:	5.0°C ~ 59.9°C (41.0°F~139.8°F)
Room temperature accuracy:	±1°C (±2°F)
Measure time:	1 second
Measure interval:	1 second
Dimension:	130.1(L)x45.0(W)x50.3(D)mm
Weight:	About 62.0g (without battery)
Forehead and object operating condition:	15.0°C~40.0°C (59.0°F~104.0°F) with a relative humidity of 15%~85% Atmosphere pressure: 700hPa-1060hPa
Storage and transport condition:	-20°C~55°C (-4°F~131°F) with a relative humidity of 15%~90% Atmosphere pressure: 700hPa-1060hPa

Service life: 2 years

Service life of the battery: With a new battery (Carbon) approx. 6 months based on the use frequency at 5 times/day

### Caution:

This thermometer has been thoroughly tested and inspected to assure proper performance and operation!

## NORMALIZED SYMBOLS

	Applied part of type BF (The applied part is the probe tip)
	Disposal in accordance with Directive 2012/19/EU (WEEE)
<b>CE0197</b>	Complies with the European Medical Device Directive (93/42/EEC) and amended by directive 2007/47/EC requirements. Notified body TÜV Rheinland (CE0197)
	The name and the address of the manufacturer
	The name and the address of the Authorized EC-representative in Europe
	Caution
	Refer to Instruction Manual
<b>IP22</b>	The first number 2: Protected against solid foreign objects of 12,5 mm Ø and greater. The second number: Protected against vertically falling water drops when enclosure titled up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is titled at any angle up to 15°, on either side of the vertical.
	Date of Manufacture

Shenzhen Dongdixin Technology Co., Ltd.  
Floor 1-2, No.3 Building, Fanshen Xusheng Industrial Estate  
Xilixiaobaimang 518108 Nanshan District, Shenzhen P. R. China

MedEnvoy Switzerland  
Gotthardstrasse 28, 6302 Zug, Switzerland

**EC** **REP**

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestraße 80, 20537 Hamburg Germany

**CH** **REP**

MedNet SWISS GmbH  
Bäderstrasse 18, 5400 Baden, Switzerland

Model: TH1009N / THERMO DISTANCE

# Terrillon

**Terrailon SAS France & Headquarters**

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**Terrailon Asia Pacific Ltd**

4/F, Eastern Centre  
1065 King's Road  
Quarry Bay – Hong Kong  
Tel: +852 (0)2960 7200  
customerservice@terraillon-asia.com

**Terrailon Corp USA**

contact@terraillon.fr

**Terrailon**  
[www.terraillon.com](http://www.terraillon.com)

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**Terrailon UK Limited**