



VI SAFETY INSTRUCTIONS

Read all safety instructions before using and use correctly.

Follow the operating instructions to use this device.

1. Read all safety instructions before using and use correctly.

2. Follow the operating instructions to use this device.

3. Do not pull down or remake the device, please consult the distributor for repairing.

4. Only use the device for its intended purpose as described in this manual. Do not use the device in the high temperature and humidity place

5. Do not store or use the device in the high temperature and humidity place such as under sunshine or near fire or in the bathroom, etc. to avoid fire or electric shock. Do not use or place it in or near the high temperature appliance.

6. Do not use and place the device in the humid condition, in case of any short circuit, electrical shock or fire.

7. Keep away from children and physical disabled people allow to use, close supervision is necessary when this device is operated by user under 18 years of age.

8. Keep out of reach of the children.

9. Pregnant, and people with hypertension, heart disease are not allow to use.

10. The person who is sensitive to the skin diseases excessively is not allowed to use.

WARNING

1. To avoid malfunction or fire, please stop use in case of unusual hot.

2. Please do not drop or clash the device.

3. Please apply serum constantly during use.

4. Please do not use it if the power cord is damaged.

5. Below conditions are not allow to use:

Heart disease, fever, communicable disease, cancer, open cuts, acute disease, Hypotension, hypertension, install pacemaker, artificial heart or lung, during medical treatment, skin disease, sunburned, over sensitive people, anapetia.