C CARE

SMC 120

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Service adres:

Cresta International B.V
Rolbrugweg 4 | NL-1332 AS ALMERE
info@cresta-international.com
www.cresta-international.com

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INSTRUCTIONS

We congratulate you on purchasing the C Care Wireless -, Rechargeable Neck - and Shoulder Massager SMC120 and thank you for the confidence placed in our brand and products.

We are convinced that with this purchase you have acquired a high quality product that you will enjoy using for many years to come.

Please read these operating instructions carefully before using the product, even if you are familiar with the type of product. It gives you useful tips to get the most out of your product.

If, after reading this manual, you still have questions about this product and how to use it, please send an e-mail to the address below; we will answer your question immediately

CRESTA INTERNATIONAL BV

Rolbrugweg 4

NL-1332 AS ALMERE

The Netherlands

info@cresta-international.com

www.cresta-care.com

Tel: +31 36 3030063

PRODUCT DESCRIPTION

This electronic wireless - , rechargeable neck and shoulder massager will provide you a penetrating and relaxing massage with the intention to relax and loosen your muscles. The use of the product will relieve muscle aches and pain, will reduce stress, relieve tension and will improve your blood circulation.

SAFETY INSTRUCTIONS

To ensure a safe and proper use of this product read this instruction manual carefully first

USAGE ENVIRONMENT

- Do not use outdoors
- Do not use the appliance in a wet or dusty environment as this may result in malfunction or an electric shock.
- Use this appliance only for its intended use as described in this manual
- Do not use the product if you notice any damage on the cord or plug or if it is not working properly.
- Do not use the product if it has dropped in the water
- Keep the cord away from heated surfaces
- Do not operate under a blanket or pillow as excessive heat may occur and may cause fire, electrocution or any other physical injury.
- If you are pregnant consult your doctor before using this product
- If you have a pacemaker or other aid consult your doctor before using this product.
- If you have any concerns about your health consult your doctor before using this product.
- This product is not a substitute for any medical treatment
- Do not use this product directly on swollen or inflamed areas or skin eruptions.
- Do not use on insensitive skin or with persons with a poor blood circulation.
- We recommend you not to use the massager for more than 15 minutes at a time.
- This product is a non-professional appliance designed for home use and to provide a soothing massage to worn muscles.
- Do not use the product if it is wet
- Do not insert any objects like pins - into the appliance
- This product should not be used by individuals suffering from any physical ailment that could limit the user's capacity to operate the control or who has sensory deficiencies in the lumbar
- Children, invalids or disabled should only use this product under supervision.
- Do not leave the appliance unattended when it is plugged in. Unplug it from the socket when it is not in use.
- Always unplug the appliance from the socket after using and before cleaning
- Do not use while bathing or in the shower
- Do not place or store the appliance where it can fall or be pulled into a tub or sink. Store it in a safe location.

C CARE

- Do not expose the appliance to a heater or direct sunlight.
- Do not drop or insert any object into any opening
- Handle the heated surfaces with great care . The unattended use of the heated surfaces by children or disabled may be dangerous
- Do not lift or carry this appliance by the power cord
- Do not cover the appliance during the operation
- The appliance is not intended for persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge. They can use the product under supervision of a person responsible for their safety.
- Extensive use can lead to excessive heating and physical damage. Should this occur discontinue the use and allow the unit to cool down before operating it again.
- For people under medication or medical treatment please consult your doctor before using the product. Do not use the product if you do not feel well.
- Do not open the product, do not repair by yourself.

MAIN UNIT

OPERATION AREA

ON/OFF BUTTON
REVERSE BUTTON
SPEED BUTTON
HEATING BUTTON

CHARGING

This neck and shoulder massager is a rechargeable -, wireless product with integrated rechargeable batteries included. Charging the batteries to 100% will take about 2 hours.

With fully charged batteries the product can operate 75 to 90 minutes.

To charge connect the adapter to the appliance and the socket.

Attention:

The product can not be used during charging

We advise you to fully charge the batteries before the first use.

OPERATING INSTRUCTIONS

Switch the power on by pushing the power ON/OFF button for 3 seconds. The massage heads will start rotating and the heating will be switched on

If you wish to increase the massage speed push the SPEED button. There are 3 speedlevels selectable - low - middle and high.

If you wish to stop the heating during the massage please push the HEATING button. This button has the ON or OFF function.

By pushing the REVERSE button you can select the massageheads direction - forwards or reverse.

Switch off the product by pushing the ON/OFF button for 3 seconds.

This product has an integrated timer of 15 minutes. After 15 minutes of usage the product will be switched off automatically.

PRODUCT MAINTENANCE

STORING THE UNIT

- Place the massager preferably in the original giftbox and store it in a safe, dry and cool place. Avoid contact with sharp edges or pointed objects that might cut or puncture the surface.
- Do not wrap the power cord.

CLEANING

Unplug the unit and allow it to cool down before cleaning

Clean only with a soft, damp sponge

Do not allow water or any other liquid to go inside the unit

Avoid immersing the unit into any liquid during the cleaning

Do not use abrasive cleaners, brushes, gasoline, kerosene, glass/furniture polish or paint thinner to clean the product

PRODUCT SPECIFICATIONS

Product name	:	3D Wireless Neck and Shoulder Massager
Reference	:	SMC120
Power	:	25,2 W
Power Frequency	:	50-Hz
Adapter output	:	DC 12,6V, 2A

ENVIRONMENT FRIENDLY DISPOSAL

The European Directive 2002/96/EC on Waste Electrical and Electronic Equipment - WEEE - requires that old household electrical appliances should not be disposed of in the normal unsorted municipal waste stream. Old appliances should be collected separately in order to optimize the recovery and recycling of the materials they contain and reduce the impact on human health and the environment.

The crossed out "wheeled bin" symbol on the product reminds you of your obligation that when you dispose of the appliance, it should be separately collected. Contact your local authority for information concerning the correct disposal of your old appliance.

WARRANTY

CRESTA INTERNATIONAL B.V. warrants this product to be free from manufacturing 'defects and material errors for a period of 36 months from the date of purchase. In the event of a defect within this period, CRESTA will at its discretion repair or replace the device free of charge for the same or an equivalent device. This warranty does not apply in the event of normal wear and tear or if, in the opinion of CRESTA INTERNATIONAL B.V. the defect was caused by improper use, falling or impact, water damage, leaking batteries, if applicable, and the like. The warranty is only applicable if the warranty claim can be demonstrated with the purchase receipt.

Do not try to open the product and repair it yourself. Your warranty claims will lapse

Requirements:

If a repair is necessary, you can return the product to the store where you bought it. You can also send the product directly to us. Send a message in advance with your complaint description and do not forget to state your full address details - not a PO box - when sending so that we can return the product to you after repair or replacement. If you send the product to us, you must frank it sufficiently. After repair or exchange you will receive your product from us delivered free.

Service address:

Cresta International B.V

Rolbrugweg 4 | NL-1332 AS ALMERE

info@cresta-international.com

www.cresta-international.com

PRINTING ERRORS AND MODIFICATIONS RESERVED

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