FITAGE

USER MANUAL

English	6
Nederlands	19
Français	33
Deutsch	47
Español	61
Italiano	74

SCAN THE QR CODE

to view the digital manual
in all languages

EN

WE WANT TO LISTEN TO OUR CUSTOMERS

Your opinion is everything for FITAGE. Your experience drives us to help contributing to improve people's health and fuels our community. If you like your FITAGE product, please share your experience with us!

1. Go to Amazon / Bol.com / fitage.com
2. Click "Your Orders" and search for the FITAGE product.
3. Click on "Write A Review" and tell us about your experience!

ANY PROBLEMS?

If you encounter any problems, please contact us via Email.

We are available on Monday - Friday 09.00-17.00 (+1 GTM) in English, Dutch or German. Support Email:

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EN: Scan the QR code and discover more about our FITAGE community!

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WELCOME TO **FITAGE** GET IN SHAPE

Introducing the FITAGE mini massage gun, a groundbreaking device crafted to enhance your health and well-being. This massage gun is designed to deliver deep tissue massage and promote relaxation, aiding in quicker recovery and improved performance. Whether you're an athlete striving to enhance your abilities or someone seeking relief from everyday stress and tension, the FITAGE massage gun is the ideal solution to support your goals. With its advanced features and easy-to-use design, this massage gun is essential for anyone aiming to reclaim their health. Don't hesitate – invest in your wellness today with the FITAGE mini massage gun!

SPECIFICATIONS AND CERTIFICATIONS

Product Name	FITAGE Mini Massage Gun
Product Model	FITAGE-23
Massage heads	8 Pieces
Motor Type	Brush Motor
Battery Type	18650 Lithium Battery
Power Supply	USB-C
Exterior Shell	ABS
Input Voltage	7.4v
Net Weight	640 gram
Rotation Rate	1800 – 3200 rpm
Noise Level	45-50dB
Certifications	CE, RoHS

The product and its packaging are made from materials that can be recycled, which helps reduce waste and protect the environment. Please dispose the packaging properly by utilizing local facilities designated for collecting paper, cardboard, and lightweight packaging.

Proper disposal of your massage gun is essential to protect the environment and human health. Reach out to local facilities to dispose this electronic device. The Waste Electrical and Electronic Equipment (WEEE) Directive provides guidelines for the responsible handling and disposal of electronic devices.

WHAT'S INCLUDED

- 1x FITAGE Massage Gun
- 8x Massage Heads
- 1x Instructions Manual (EN/NL/FR/DE)
- 1x USB-C Cable with EU Power Plug
- 1x Suitcase

INSTRUCTIONS FOR USAGE

The FITAGE mini Massage Gun alleviates discomfort in sore and stiff muscles by focusing on specific muscle groups using 8 customized massage heads, providing deep relief and reducing recovery time.

PRODUCT BREAKDOWN

1. Massage head
2. Air Outlet
3. USB-C Charging Port
4. Touch Display
5. Handle

TOUCH DISPLAY

1. Electric Power Display & Speed Level
2. Decrease Speed Level
3. USB-C Charging Port
4. Increase Speed Level
5. On/Off Switch

OPERATIONS

1. To start the massage gun, press and hold the On/Off button for approximately two seconds. Similarly, to turn off the machine, press and hold the On/Off button for about two seconds.

**When the touch display screen lights up, it indicates that the device is on standby.*

2. Increase or decrease the speed setting from 1 to 30 by tapping the "+" or "-" icons.
3. Begin massaging desired body parts at a low speed, gradually increase the speed if you want the massage to be stronger. You don't have to apply a lot of pressure; the massage gun does most of the work.
4. Massage back and forth across the muscle and avoid keeping the massage gun in the same spot for longer than five seconds. If you are massaging a more tense area, you can hold a little longer. Be sure not to get too close to harder bone areas.

REMINDERS

- Do not use this massage gun while it's charging. Please disconnect it from the charger before using.

- Note that applying too much pressure can lead to muscle damage.
- After every 20 minutes of use, the massage gun will pause automatically. Press the "+" icon to resume.
- If the device is not charging, check the connection between the USB cable and the device. Ensure that you're using an original 5V-2A 7.4v charger provided by the supplier.

INSTALL AND SWITCH MESSAGE HEADS

1. Press the massage head gently into the massage gun
2. Gently pull out and remove the massage head

DIFFERENT MASSAGE HEADS

ROUND HEAD

Palm,
Back,
Buttocks
and Legs

LONG TAPER

Eases deep tis-
sue pain in areas
like the knees,
joints, palms,
and soles

FLAT HEAD

Back,
Buttocks,
Legs,
Calves

Y-SHAPE

Alleviates
stiffness in
the spine and
Achilles tendon

GYRO HEAD

Joints and
Sensitive
Muscles

SHORT TAPER

Eases deep
tissue pain
in areas like
the knees, joints,
palms, and soles

SPADE HEAD

Back,
Buttocks,
Legs and
Calves.

THUMB HEAD

Targets the Palm
and Plantar
for hand
and feet.

MAINTENANCE

- To clean the massage gun and its massage heads, use a slightly damp cloth followed by gentle drying with a soft cloth.
- Ensure you regularly clean the air outlet to keep it free from debris and dust.
- When storing or traveling with the massage gun, ensure all items are placed inside the travel case and securely close it.

FREQUENTLY ASKED QUESTIONS

What are the benefits of a massage gun?

Massage guns help enhance mobility and support overall well-being. They've been demonstrated to improve range of motion, boost athletic performance, sharpen mental focus, enhance blood circulation, alleviate back and muscle discomfort, reduce stress, and aid muscle recovery post-workout, among other benefits.

The massage gun automatically turns off after 20 minutes of use. Is there any way to prevent this?

Yes, the massage gun is designed to automatically turn off after running for more than 20 minutes or if placed on a hard surface with too much pressure. To prevent this, simply press the "+" button on if you need to continue using it.

How loud is the massage gun when set to its highest level?

The FITAGE Massage Gun is designed to provide high power while keeping noise levels low, ranging from 35dB to 50dB. It is more quiet than standard massage guns on the market.

How frequently should you use a massage gun?

It varies for each person. However, using it too often, directly on bones, or for too long in one area could exacerbate the issue. If you're new to massage guns, gauge your tolerance by starting with no more than one session per day.

Is it safe to use a massage gun for neck pain?

With caution, yes, you can use a massage gun on your neck. However, it's essential to assess your tolerance level first. Applying too much pressure could potentially exacerbate neck pain.

Is it safe to use a massage gun during pregnancy?

If you're pregnant, always consult your doctor before using a massage gun or any other therapeutic device.

Can a massage gun alleviate lower back pain?

Yes, a massage gun can be effective for relieving lower back pain. However, it's important to consult your doctor before using one.

What should I do if my massage gun won't turn on?

1. First, make sure you have properly charged the device

using the included power charger. The charger light will turn full on once the battery is fully charged.

2. If that doesn't work, check that you press longer than two seconds on the On/Off button.
3. If the massage gun still won't turn on, please contact our customer support team at support@fitage.nl

SAFETY WARNING

Always seek advice from a doctor prior to using the massage gun or any vibration device.

May not be appropriate if you are pregnant, diabetic, have neuropathy or retinal damage, wear a pacemaker, have recently undergone surgery, experience migraines, or have epilepsy, among other conditions.

1. To prevent damage to your property or injury to yourself and to ensure that the product continues to function properly over time, please read and follow these instructions carefully while using the product.
2. This unit is not suitable for individuals who are 12 years old or younger.
3. In order to minimize the risk of injury, do not directly massage the skin with the device. Wear clean and dry clothing over the massage area and use gentle pressure for no more than 60 seconds to prevent abrasions.

4. Use high speed and high pressure settings only on soft tissue areas of the body to prevent pain and discomfort. Do not use high speed and pressure on the head or any bones to avoid injury.
5. Keep your fingers and hair away from the telescopic rod and motor heat dissipation hole to avoid pinching injury.
6. Use the device in a dry environment with normal room temperature to maintain its long-term capability.
7. Handle the device gently and avoid dropping or hitting it, as this can damage the device.
8. This device is for household use only and should not be used commercially.
9. Do not run or charge the device in an unmanned state.
10. Do not use the device to diagnose or treat any medical issues, and consult a doctor before using it if you have any health problems.
11. Do not use the device if you are pregnant, have diabetes, neuropathy, retinal damage, wear a pacemaker, or have undergone surgery.
12. If the device malfunctions, check the battery level first and charge it if necessary. Do not attempt to open or modify the device. Use the original charger and remove it after the battery is fully charged.

13. To avoid damaging the device, do not put it in water or let water enter the openings. If the device becomes wet, use a dry cloth to clean it and do not use soap or any other cleaning agents. Also, remember to clean the motor heat sink frequently to keep it free of dust and debris.
14. Do not use the unit continuously for more than one hour.
15. To prevent any damage to the device, do not use or charge it in high-temperature environments. Also, make sure to remove the charger from the unit before using it and after it has fully charged.

WARRANTY

We provide a two-year repair warranty on our product, with exceptions including:

- Damage resulting from incorrect operation or transportation by the user.
- Unauthorized tampering or repair of the equipment.
- Failure to consult the instruction manual for assistance.
- Overuse of the product beyond its recommended capacity.
- Damage due to inadequate storage or maintenance by the customer.
- Lack of proof of purchase provided by the customer.

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