Handleiding Deluxe exact-75

Technische Specificaties

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User Manual Deluxe exact-75

Technical Specifications

- Working temperature: 10 ~ 40 ° C
- Working humidity: <75%
- Max. capacity: 20kg
- Min. capacity: 50g
- Graduation: 5g
- Power supply: 2 x AAA (LR03) batteries
- Low power indication: Lo
- Unit change: kg, lb or oz
- Turn off automatically after use
- 150 cm soft measuring tape

Operation

- Make sure the scale is placed on a flat, hard and stable surface.
- After inserting the batteries, press the "ON/OFF/TARE" button on the scale and wait until "0.00" appears on the display.
- Press the "UNIT" button to select the desired unit of weight (kg/lb/oz).
- Carefully place the baby on the scale, making sure the baby lies still in the middle of the scale. The weight of the baby will appear.
- Press and hold the ON/OFF/TARE button for 3 seconds to turn the scale off.

Keys' functions

- ON/TARE: Press to turn on the scale and initiate the startup process. Press again to set the scale to zero. Press for 3 seconds to turn the scale off.
- UNIT: Use this button to adjust the weight unit. You can toggle between kg, lb and oz. The weight unit is displayed on the right side of the screen. Long press this button to lock or release the measurement.

Device is made in class III of insulation.
Device is compliant with EU directives:
- electromagnetic compatibility (EMC)
Device marked CE mark on rating label

 To protect your environment: please separate carton boxes and plastic bags and dispose them in corresponding waste bins. Used appliance should be delivered to the dedicated collecting points due to hazardous components, which may affect the environment. Do not dispose this appliance in the common waste bin.

Cautions

- Do not expose the scale to moisture.
- Do not drop the scale.
- Do not immerse the scale in water.
- Do not expose the scale to extreme heat.
- Do not overload the scale, this will permanently damage the scale.
- Do not leave weight on the scale when not in use.
- Never leave the baby unattended on the scales.
- Before use, place the scale on a hard and flat surface for a good and safe measurement.
- Replace the batteries immediately when the display shows "Lo".
- Insert the two AAA batteries correctly into the battery compartment, make sure the plus and minus poles are correctly aligned, and firmly close the battery cover again.
- Keep the scale clean to ensure normal operation.
- Remove the batteries from the scale when you do not use it for a longer period of time.
- Do not place more than 20 kg on the scale.

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