VRL *SPORT*

MASSAGE GUN HANDLEIDING

(NL) (EN) (DE)

VRL Shop
Postbus 4189,3130 KD, Vlaardingen, Nederland
info@vrl-shop.nl | www.vrl-shop.nl

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Gewicht
1.1 kg / 2.5lbs

Opladen
110-240V 50/60Hz-5A-26VDC-1A

(EN)
- Powerful brushless high-torque motor with Quiet Gbde technology.
- Removable battery.
- Rechargeable lithium-ion battery. (8 hours of use per charge)
- Multiple speeds for optimal pressure.
- Cordless and lightweight (2.5 lbs) including AC wall charger.
- On/off switch.
- Interchangeable massage heads (x6).

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(EN) 1: Keep away from liquids or heat sources. 2: Do not remove screws and do not attempt to disassemble. 3: Do not use for more than one hour at a time. Let the device rest for 30 minutes before using again. 4: Disconnect the charger from the battery after charging or before use. 5: For adult use only. Do not use if you are injured. Consult your doctor before using this product.

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(EN) CAUTION: TO REDUCE THE RISK OF ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

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(EN) 1: For adult use only. 2: Only use on a dry and clean surface of the body. Use over clothing by lightly pressing and moving over the desired area for about 60 seconds. 3: Use the VRL Massage Gun only on soft body tissue to avoid causing pain or discomfort. Do not use on the head or any hard part of the body. 4: Only use the attachments that give the best desired effect. 5: Bruising can occur regardless of the setting or pressure applied. Check the treated areas regularly and stop immediately at the first sign of pain or discomfort. 6: To avoid pinching, keep fingers, hair, and other body parts away from the shaft and back of the VRL Massage Gun. 7: Do not place objects in the ventilation openings. 8: Do not immerse in water, and do not let water or dirt enter the ventilation openings. 9: Do not drop or misuse the VRL Massage Gun. 10: Charge only with the provided charger. 11: Carefully inspect the device and battery before each use. 12: Do not modify the VRL Massage Gun in any way. 13: Never operate or charge the VRL Massage Gun unattended.

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(EN) 1: Fully charge the battery up to six hours before the first use, both connected and disconnected from the device. 2: To charge, plug the end of the charger into the charging port of the battery and plug the side adapter into the socket. 3: Battery lights flash to indicate battery level and active charging. 4: When the battery is fully charged, all lights stop flashing and remain on. 5: The battery can be charged at any time and at any battery level. 6: It is not recommended to fully discharge the battery until the red indicator light is on.

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(EN) 1: Attach the battery to the base of the device. Make sure it is securely locked in place. 2: Turn the device on by setting the on/off button at the bottom of the battery to 'on'. Turn off the device by setting the on/off button to 'off'. 3: With the on/off button in the 'on' position, press the on/off button once to set the device to level 1. Press the 'on' button twice to set the device to level 2. Press three times for level 3. If you press the on/off button again, the device will be turned off. The blue LED lights indicate the power level. 4: Massage the desired body part with the preferred speed, applying pressure if necessary. 5: To remove the battery, lightly press the battery release button with your thumb. Pull the battery down firmly.

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(EN) Use a slightly damp cloth to clean the device and dry it with a soft cloth. Remove the battery during storage or when traveling and store it in the supplied cardboard box.

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(EN) DO NOT USE THE VRL MASSAGE GUN OR ANY OTHER VIBRATION DEVICE WITHOUT APPROVAL FROM YOUR DOCTOR, IF ANY OF THE FOLLOWING APPLY: 1: Pregnancy, diabetes with complications such as neuropathy or retinal damage, wear and tear of pacemakers, recent surgery, epilepsy or migraine, hernia, spondylolisthesis, spondylolol-ysis of spondylosis, recent weight replacements or intrauterine devices, metal pins or plates, or any concerns about your physical health. 2: Vulnerable individuals and children must be supervised by an adult when using a vibration device. These contraindications do not mean that you are unable to use a vibration or training device, but it is recommended to consult a doctor first. 3: Ongoing research is being conducted into the effects of vibration training on specific medical Conditions that is likely to result in shortening the contraindication list as described above. Practice shows that integration of vibration exercise in a treatment plan is recommended. This should be done under the advice of a doctor, specialist, or physiotherapist.

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(EN) PLEASE CONTACT VRL IF YOU ARE NOT SATISFIED WITH THE PRODUCT OR IN CASE THERE IS ANY DEFECT IN THE PRODUCT.

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