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560*100mm

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Product structure

Skin head
Ionic conductor
Status display area
Power switch
Mode select key
USB charging interface

◀Other accessories▶ ●The instructions ●USB charging line

Product characteristics

Ultrasonic wave is a kind of high frequency mechanical vibration wave, which is widely used in medical and beauty. The mechanical wave has the mechanical action, the temperature work and catalysis, at low energy high frequency vibration can be on the human skin produce micro resonance effect, enhance the vitality of cells and cortical groups. The permeability between the weaves, dredge sweat glands, shrink pores.

Ultrasonic cleansers and beauty meters use ultrasound and nanotechnology to improve the look. That's 24,000 vibrations per second, which makes the molecules violently rub against each other. Vibrations can go deep into the skin, causing intense friction between molecules and raising the temperature in the deeper layers. High cell metabolism is exuberant, microcirculation is accelerated, promote metabolic to return flow. Through the conversion of ion flow, the permeability of cell membrane is enhanced, so that the semi-permeable membrane produces intracellular or protoplasmic microflow, which is the dissimulation of fatty acids into CO₂ and H₂O fat can reduce the delicate volume, and at the same time will in certain pertinence. Products, beauty effect is more significant.

Through high-frequency vibration for deep massage, in the same fat movement. When improving microcirculation, make fat cell combustion. I ysis into metabolism. The product is discharged from the body to achieve the cleaning effect. Special shovel type a probe, which can effectively enter 4-7 cm under the skin and pass through the acoustic wave. The vibration on fibrils and atomizes the beauty lotion, removing the dirt and excessive burn to be emulsified, so that the plug in the pores of cosmetics residue, sweat. The scale emulsifies and, with the four spatulas removed, this cleanser is more effective. Faster and more efficient. This is ultrasound washing like a precision washing machine. This happens to the skin as well as a s the purifier. To remove. Conspicuous layer, adipose grain, which, let your m ainain the nutrient of article. From now on no longer waste, so that the sk in better absorption.

Product features

- 1 Deeply clean and remove aged cutin.
- 2 Pigment decomposition, fade spot.
- 3 Deeply remove wrinkles and enhance skin elasticity.
- 4 Promote blood circulation and lymphatic metabolism.
- 5 Enhance skin cell regeneration and permeability.
- 6 Massage skin, nutrition import, mottled export.
- 7 No skin damage or side effects.
- 8 Beautiful appearance, easy to use, simple operation, easy to carry.

Product parameters

Product name: ultrasonic beauty peeling machine
Model: C-1
Rated voltage: ac100v-240v 50HZ/6HZ
Output voltage: 5V/600mA
Charging mode: rechargeable (lithium battery)
Net content: body is about 99g
Power: 2 w
Size: 55 * 165 * 15.5 MM

The charging method

- 1 USB charging cable is connected to USB interface
The battery is under full state before shipment. Please charge it before use.
- 2 Plug in the computer, adapter (sold separately) or the phone's lithium battery is charged via a USB port
The LED light flashes continuously during charging, and lights up after charging. Please confirm the status of lithium battery before charging. Before use for the first time, please tear off the sealing film of skin brightening head before use.

Charging precautions

The captain presses the power button and does not show the starting state, please charge. When the product is running, a "beep-beep" note will appear, and then the machine will be shut down. Status indicator that the battery is not enough. Please charge the product immediately. The product will automatically be in standby state after working for 10 minutes, if not used for one minute, it will automatically shut down. It takes about 150 minutes for the product to be fully charged. If it needs to be charged for more than 5 hours in the first use, it can be used continuously for about 80 minutes after being fully charged. It is recommended to charge the product regularly.

Method of use

Cleaning mode

- 1 Long press the power button for 2 seconds to switch on the power
- 2 Select the shaving mode on a one-time basis
It is recommended to use facial cleanser first. Use after cleaning your face
- 3 Apply the product to your face
- 4 As shown in figure 2, there is a brightening head for cleaning

Protect wet model

- 1 Long press the power button for 2 seconds to switch on the power
- 2 Click Moisturizing mode twice
For cleaning export mode, it is recommended to use facial mask first. Use after abs orbing nutrition
- 3 Apply the lotion to your face
- 4 As shown in figure 2, import with the skin brightening head

Method of use

Ascension patterns

- 1 Long press the power button for 2 seconds to switch on the power
- 2 Upon triple selection upon Lifting
Tight lifting mode, recommended with the concentrate. Use together, the effect is better
- 3 Apply the serum to your face
- 4 As shown in FIG. 2, pull in with a brightening head

Maintenance and storage methods

- 1 Make sure the power is off
- 2 Please use wet wipes or alcohol for hair and host stains. Fine wipe. Do not use diluents, solvents and other corrosive chemicals to wipe the product.
- 3 After cleaning, please keep it in a dry place. Do not put in the bathroom and other wet places for safekeeping. Keep out of reach of children.

Ultrasonic beauty peeling machine

Operating instruction

Product structure

Skin head
Ionic conductor
Status display area
Power switch
Mode select key
USB charging interface

◀Other accessories▶ ●The instructions ●USB charging line

Product characteristics

Ultrasonic wave is a kind of high frequency mechanical vibration wave, which is widely used in medical and beauty. The mechanical wave has the mechanical action, the temperature work and catalysis, at low energy high frequency vibration can be on the human skin produce micro resonance effect, enhance the vitality of cells and cortical groups. The permeability between the weaves, dredge sweat glands, shrink pores.

Ultrasonic cleansers and beauty meters use ultrasound and nanotechnology to improve the look. That's 24,000 vibrations per second, which makes the molecules violently rub against each other. Vibrations can go deep into the skin, causing intense friction between molecules and raising the temperature in the deeper layers. High cell metabolism is exuberant, microcirculation is accelerated, promote metabolic to return flow. Through the conversion of ion flow, the permeability of cell membrane is enhanced, so that the semi-permeable membrane produces intracellular or protoplasmic microflow, which is the dissimulation of fatty acids into CO₂ and H₂O fat can reduce the delicate volume, and at the same time will in certain pertinence. Products, beauty effect is more significant.

Through high-frequency vibration for deep massage, in the same fat movement. When improving microcirculation, make fat cell combustion. I ysis into metabolism. The product is discharged from the body to achieve the cleaning effect. Special shovel type a probe, which can effectively enter 4-7 cm under the skin and pass through the acoustic wave. The vibration on fibrils and atomizes the beauty lotion, removing the dirt and excessive burn to be emulsified, so that the plug in the pores of cosmetics residue, sweat. The scale emulsifies and, with the four spatulas removed, this cleanser is more effective. Faster and more efficient. This is ultrasound washing like a precision washing machine. This happens to the skin as well as a s the purifier. To remove. Conspicuous layer, adipose grain, which, let your m ainain the nutrient of article. From now on no longer waste, so that the sk in better absorption.

Product features

- 1 Deeply clean and remove aged cutin.
- 2 Pigment decomposition, fade spot.
- 3 Deeply remove wrinkles and enhance skin elasticity.
- 4 Promote blood circulation and lymphatic metabolism.
- 5 Enhance skin cell regeneration and permeability.
- 6 Massage skin, nutrition import, mottled export.
- 7 No skin damage or side effects.
- 8 Beautiful appearance, easy to use, simple operation, easy to carry.

Product parameters

Product name: ultrasonic beauty peeling machine
Model: C-1
Rated voltage: ac100v-240v 50HZ/6HZ
Output voltage: 5V/600mA
Charging mode: rechargeable (lithium battery)
Net content: body is about 99g
Power: 2 w
Size: 55 * 165 * 15.5 MM

The charging method

- 1 USB charging cable is connected to USB interface
The battery is under full state before shipment. Please charge it before use.
- 2 Plug in the computer, adapter (sold separately) or the phone's lithium battery is charged via a USB port
The LED light flashes continuously during charging, and lights up after charging. Please confirm the status of lithium battery before charging. Before use for the first time, please tear off the sealing film of skin brightening head before use.

Charging precautions

The captain presses the power button and does not show the starting state, please charge. When the product is running, a "beep-beep" note will appear, and then the machine will be shut down. Status indicator that the battery is not enough. Please charge the product immediately. The product will automatically be in standby state after working for 10 minutes, if not used for one minute, it will automatically shut down. It takes about 150 minutes for the product to be fully charged. If it needs to be charged for more than 5 hours in the first use, it can be used continuously for about 80 minutes after being fully charged. It is recommended to charge the product regularly.

Method of use

Cleaning mode

- 1 Long press the power button for 2 seconds to switch on the power
- 2 Select the shaving mode on a one-time basis
It is recommended to use facial cleanser first. Use after cleaning your face
- 3 Apply the lotion to your face
- 4 As shown in figure 2, import with the skin brightening head

Protect wet model

- 1 Long press the power button for 2 seconds to switch on the power
- 2 Click Moisturizing mode twice
For cleaning export mode, it is recommended to use facial mask first. Use after absorbing nutrition
- 3 Apply the lotion to your face
- 4 As shown in figure 2, import with the skin brightening head

Method of use

Ascension patterns

- 1 Long press the power button for 2 seconds to switch on the power
- 2 Upon triple selection upon Lifting
Tight lifting mode, recommended with the concentrate. Use together, the effect is better
- 3 Apply the serum to your face
- 4 As shown in FIG. 2, pull in with a brightening head

Maintenance and storage methods

- 1 Make sure the power is off
- 2 Please use wet wipes or alcohol for hair and host stains. Fine wipe. Do not use diluents, solvents and other corrosive chemicals to wipe the product.
- 3 After cleaning, please keep it in a dry place. Do not put in the bathroom and other wet places for safekeeping. Keep out of reach of children.

Ultrasonic beauty peeling machine

Operating instruction