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EN DERMAWAV has been developed for cosmetic facial application. With the help of its innovative combination of different functions, it supports the visual appearance of the skin.

FUNCTIONS

- Microcurrent
- Vibration
- Heat
- Red light

Smart-touch technologie

Easy to use. Intelligent touch activation with an electronic intelligent sensor chip: the Derma Wav is activated when its head comes into contact with the skin. Independently adjustable head: **Closed:** for the entire face - larger areas such as forehead, as well as the lips and neck. **V position:** for eyes, cheeks, chin, neck

CHARGING BEFORE USING THE DEVICE

To charge the device, plug one end of the charging cable into the charging port on the underside of the product and the other end into the USB port on your computer. The indicator light will flash, meaning that the device is charging. Leave the device attached until the indicator light stops flashing and stays continuously lit. The charging time for a full charge is 90 minutes. **Note:** The device will not work when it is connected for charging. If using the device expends a complete battery charge, a full charge cycle of 90 minutes is recommended. When using for the first time, peel the transparent protective film from the product head. Pay attention to the direction of the rotating head. Do not turn the head in the opposite direction. This will damage the head.

CLEANING

Cleaning the device: To clean your device, use a soft, lint-free cloth slightly dampened with water or rubbing alcohol. Never immerse the device or hold it under running water.

USE:

For best results, it is recommended that you use Derma Wav after first cleaning the skin. Cleanse your face and neck with an oil-free cleanser. Note: Products containing oil reduce the penetration of the stimulation current. For best results, use Derma Wav on well-moisturised skin. The best way to do this

is to use a moisturising serum or lotion, apply a few minutes before treatment, and allow it to absorb. Derma Wav uses intelligent smart-touch activation. To turn Derma Wav on and off, simply place the treatment head on the skin of the face to turn it on, and remove it to turn it off. Gently stroke the treatment head over the face or the areas you want to treat with medium pressure. When used correctly, the treatment feels very gentle, beneficial and mild.

RECOMMENDATIONS

Use the following methods: 

FREQUENCY OF TREATMENT • Use the Derma Wav for at least 10 minutes each day. The daily treatment can last for as long as up to 30 minutes. • For the first 60 days, you should use Derma Wav at least 5 times per week. • After that, you should use the product 2-3 times a week to maintain the achieved results.

WARNINGS

- Please follow the instructions / operating manual.
- Never forcibly turn the head.
- Keep the device away from water (this includes bathtubs, showers and sinks).
- If the device starts to overheat, use immediately.
- Stop using the device immediately if you suspect that it is not working properly or has become defective. Do not disassemble or modify the device as this could result in damage, malfunction or personal injury. The device contains no user-serviceable parts.
- Do not use on swollen, inflamed areas of skin, open wounds or rashes, or on broken capillaries, varicose veins, etc.
- Cease treatment with the Derma Wav immediately if you feel unwell.
- Keep the device out of the reach of children.
- A slight tingling sensation may be felt when used around the eyes; this is a normal reaction.
- Do not use if you have recently had facial or other related surgery.
- If you have any concerns, speak to your doctor before using the Derma Wav.

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