

U-shaped electric toothbrush

Rated voltage: 5V

Body: ABS

Charging time: <120Minute

Rated power: <5W

Brush head: Food grade silicone




Executive standard: GB4706.1-2005/GB4706.59-2008



We recommend to replace the brush head within 3 months

USE METHOD:



STEP1



STEP2



STEP3

Apply toothpaste on the two
sides brush a head to
squeeze toothpaste

Put it in your
mouth and press
the switch

Shake the toothbrush
from side to side to
clean your teeth

FDA CE FC RoHS

MADE IN CHINA