
# CONTACTLESS INFRARED THERMOMETER

— User instructions v1.1 —

(Model: DT-8806C)

- ENGLISH(P01-12)
- DUTCH(P13-25)
- Please read this operation instruction carefully before use.
- Please keep it for future reference.

# CONTENTS

---

1、 Product Introduction-----	02
2、 Working Principle -----	02
3、 Packing List -----	02
4、 Product Structure-----	03
5、 Display Icon Definition-----	04
6、 Technical Specifications-----	05
7、 Battery Installation-----	06
8、 Power On-----	07
9、 Shutdown-----	07
10、 How To Use -----	07~09
11、 Parameter Setting -----	09~10
12、 Maintenance & Storage-----	11
13、 Precautions-----	11
14、 After Sales & Warranty-----	12

## 1、 Product Introduction

Thank you for purchasing DT-8806C contactless infrared thermometer! This product is composed of high-precision infrared sensor, microprocessor, memory, power supply, electroacoustic element, LCD and shell. It can measure the temperature of human body or object quickly and accurately


## 2、 Working Principle

The product uses infrared temperature sensing method to measure the temperature to display or output the measured data. To show the temperature of the subject by measuring the forehead or the thermal radiation of the object.

## 3、 Packing List


User operation guide\*1、 AAA Battery\*2、  
Inspection & Guarantee Card\*1、 Cloth bag\*1

## 4、Product Structure

- ① Probe
- ② Measuring button
- ③ Battery case
- ④ LCD display screen
- ⑤ “+” Button
- ⑥ “MODE” Button
- ⑦ “💡” Backlight Button
- ⑧ “-” Button

## 5、Display Icon Definition


- ① Body temperature
- ② Object temperature
- ③ Temperature value
- ④ Memory
- ⑤ Memory number
- ⑥ Lights on/off
- ⑦ Low battery warning
- ⑧ Centigrade
- ⑨ Fahrenheit
- ⑩ Voice on / off

## 6, Technical Specifications

Measuring distance	3cm~5cm	
Measuring range	Body	34 °C~42.9 °C(93.2 °F~109.2 °F)
	Object	1 °C~99 °C(33.8 °F~210.2 °F)
Measurement accuracy	Body	(35-42) °C $\geq \pm 0.3$ °C/0.54 °F (35-42) °C $\leq \pm 0.2$ °C/0.4 °F
	Object	$\pm 1.0$ °C/1.8 °F
Display resolution	01 °C/0.1 °F	
Operating environment	16 °C~35 °C(68.8 °F~95 °F)RH $\leq$ 85%(Non condensing condition)	
Storage environment	-25 °C~55 °C(-13 °F~131 °F)RH $\leq$ 90%(Non condensing condition)	
Power source	DC 1.5V x 2(2*AAA )	
Energy consumption	Standby status < 8uA	
	Power on measurement status < 50mW	
Electricity indication	Low battery prompt(Voltage below 2.5V)	
Memory function	40 groups of measurement values can be memorized	
Unit conversion	°C / °F	
Turn off time	Automatic shutdown without any operation within 20 seconds	
Product size	116.6mmx45.4mmx155mm	
Main material	ABS	
Weight	113g(excluding batteries)	
life span	5 years	

## 7、 Battery Installation

Open the battery case, Install two AAA batteries.

\*Pay attention to the positive and negative poles of the battery

\*Please recycle the used batteries according to the local laws and regulations.

## 8、 Power On

Press the measurement button once, the product will automatically enter the self-test state, and enter the temperature to be measured state after the self-test


## 9、 Shutdown

After 20 seconds of inactivity, the product will automatically shut down and enters standby mode.


## 10、 How To Use

### 1) Human body temperature measurement


Align the probe of the thermometer with the forehead of the human body, with a distance of about 3-5cm. Press the measurement button once, hold for about 1~2 seconds, and the product will emit one "beeps" to complete the measurement. The LCD screen will display the temperature value.


- Normal temperature state:  
The green backlight of the display is on, and the product emits one beeps.

- Higher temperature:  
The display screen light orange and the product beeps five times as prompts .  
If the temperature exceeds the test range, the display screen shows "HI" and beeps three times.

- Body temperature not detected  
The display screen light green, shows "LO" and beeps three times.

### 2) Object temperature measurement

Short press "MODE" button to switch to object temperature MODE and display screen shows "Object TEMP".

Aim the thermometer probe at the object to be measured, the distance is about 3~5cm.

Press the measurement button for 2 seconds, the screen displays the object temperature.

### 3) View temperature records

In the displayed state, press (+ /-) button to preview the temperature records from 01 to 40 times.

4)When the battery is low, it means need to replace with new battery.

5)Short press "MODE" key to switch between Body Temperature Mode and Object Temperature Mode

### 2)"F2" Aleft value setting


In F1 state, press "MODE" button for about 3 seconds, "F-2" will be displayed on the screen, and press "+ / -" button to adjust between (36-40 degrees).

\* Note: after the function setting of F2 is completed, it needs to be press cycled once (from F2 F3 F4 F1 F2) to remember the setting of F2.

### 3)"F3"Temperature unit setting( °C/ °F)


In F2 state, press and hold "MODE" button for about 3 seconds, "F-3" will be displayed on the screen, and press "+ / -" button to switch between Celsius and Fahrenheit

## 11、Parameter Setting


### 1)"F1"Temperature trim setting

In the power on state, press and hold the "MODE" button for about 3 seconds, and the "F-1" word appears on the screen,Press the "+ / -" button to adjust.(Range:  $\pm 2^{\circ}\text{C}$ )

### 4)"F4"Voice on/off setting

In F3 state, press "MODE" button for about 3 seconds, the screen displays "F-4", press "+ / -" button to turn on or off.

## 12、Maintenance & Storage

- 1) Please store the product in a well ventilated and dry indoor environment;
- 2) Avoid dust, pollution and direct sunlight.

## 13、Precautions

- 1) This product can't replace the doctor's diagnosis, only for screening, if there is any abnormality, please see a doctor in time;
- 2) Do not use the product in the environment with humidity greater than 85%;
- 3) Do not repair the fault by yourself, please contact the manufacturer or agent;
- 4) To ensure the best measurement effect, please make sure that the probe is 3-5cm away from the main rest surface to be measured;
- 5) Transfer to the environment with large temperature difference in a short time will lead to measurement error.
- 6) When not in use for a long time, please take out the battery and store it separately.

## 14、After Sales & Warranty

This product has after-sales warranty service.  
If there is any quality problem, please contact the seller for free replacement within 15 days after purchase.  
If any quality problem occurs within 1 year after purchase, the product can be repaired for free.  
If there is more than one year, the cost of repair parts will be charged.

**Note: Return, replacement and warranty do not cover the following situations:**


- 1) Products have been artificially modified, damaged, or disassembled.
- 2) Incorrect operation lead to damage.
- 3) Irresistible factors damage (such as geological disasters, wars, etc.).

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