**Spoused massage gun  
Model Num: SP931**

---

**NL: p.3-6**

**EN: p.7-10**

**DE: p.11-14**

# Inhoud

Spoused massage gun

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## EN

1. Warnings
2. Product details
3. Installing
4. Attachments
5. User manual
6. Specifications

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# Warnings

Safety precautions before use

1. Read the instructions before use.
2. Keep massage gun out of reach of water.
3. Do not dismantle any parts of the massage gun.
4. Do not use the massage gun for more than 60 consecutive minutes. After prolonged use, wait 30 minutes until new use.
5. Do not use the device while charging.
6. Prevent the cord from getting hot.
7. Do not massage bones and joints, only massage muscles.
8. Stop using immediately if it feels painful.
9. Keep out of reach of children.
10. Avoid getting the motor and charging cable wet.
11. The massage gun is made to be used under normal household conditions.
12. Consult a doctor or professional if you are unsure if you should use the device.
13. Consult a doctor before use during pregnancy.

# Product details

Product compartments

**1** Speed level

---

**2** Display

---

**3** Pressure

---

**4** Slower

---

**5** Faster

---

# Installing

1. Easily change between attachments, easily pull out the current massage attachment and insert the desired one.

# Attachments

1. **Round head** - sensitive muscles, back, legs, waist
2. **Long taper** - to deeply penetrate the muscle
3. **Spade head** - hamstrings and glutes
4. **Thumb head** - hands and forearms
5. **Flat head** - back
6. **Y-Shape** - for stiff muscles

# User manual

1. Press the power button for 2 seconds.
2. Press the + button to set the desired speed.
3. Massage the desired body parts and easily determine the desired speed for each body part.
4. After every 20 minutes, the massage gun pauses automatically, press + to continue.
5. Easily charge the massage gun with the included USB cable. It can take up to 4 hours to fully charge the massage gun.

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