





# PerfectSkin



includes:  
beauty instrument  
charging cable  
manual



## Product Features

### Bio-microelectronics

The multi-functional beauty instrument uses the latest research and development technology - biological micro-electricity. The conventional ion function only has positive ions or negative ions, while the biological micro-electricity has both positive and negative ions, similar to RF radio frequency, which is equivalent to twice the effect of conventional micro-current.



### Warming

Warming feature can effectively relieve skin fatigue and promote blood circulation. Persistent use can relieve skin problems such as dark circles, eye bags, and puffiness. It can also be used with serums and creams to promote absorption.



### Resurfacing

The multi-function beauty instrument uses 620-630nm red light. The cells in the dermis can absorb red light and accelerate the synthesis of cells and collagen, so as to achieve the effect of firming the skin, lightening acne marks and melanin.



### Sonic vibration

The use of sonic micro-vibration can convert the serums and creams into tiny molecules, which are then introduced into the skin to promote absorption.



### Smart touch on, smart power saving

When you hold the Wand, the massage head touches the skin, the multi-functional eye beauty instrument starts to work, and stops when you take it off.



### 270° rotating guide

The innovative 270° rotation can make the Wand better adapt to different parts such as the eye area and the lip area, as well as the entire face and neck skin.

## Use

Cleanse the face to prevent facial oil from affecting the use experience.

Long-term use can effectively relieve skin problems such as eye fine lines, dark circles, eye bags, nasolabial lines, neck lines, etc., and firm and smooth skin. For different parts of the face, you can switch the appropriate angle to use according to different parts. The use method is shown in the figure below. It is recommended to use the bottom-up method to achieve the effect of lifting.



## Charging

Insert one end of the USB charging cable into the charging port at the bottom of the eye beauty instrument, and connect the other end to the charger or computer USB port.

- When the battery is low, the indicator light flashes quickly.
- When charging, the indicator light flashes slowly.
- When fully charged, the indicator light is always on. The charging time is 100 minutes to fully charge. The battery will be exhausted after 60 minutes of continuous use.
- Charging time should not exceed 24 hours.

## Precautions

- Do not disassemble and modify this product without permission.
- Do not select the guide head in the opposite direction,
- Due to the deviation of the resistance value of the human body, if there is any discomfort, stop using it immediately.
- After use, wipe with a clean cloth or wet paper towel.
- If you find any abnormality such as vibration during use, please check whether the battery is low or shut down.
- When not in use for a long time, turn off the machine and store it in a dry place to avoid sunlight exposure.
- Place this product out of the reach of children. This product cannot be used as a toy.

## Product Parameters

Item model number: RL05-C Red Light Therapy Wand  
Product material: aluminum alloy/PC  
Power: 1.5W  
Voltage: 3.7V  
Battery capacity: 350mAh  
Charging input: 5V=1A  
Size: 150\*15\*12mm