PRODUCTHANDLEIDING

 REFECTIO

Refectio
MASSAGE GUN

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VEILIGHEIDSINFORMATIE

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GARANTIE

INHOUD

PRODUCTSPECIFICATIES

TECHNISCHE SPECIFICATIES

AAN DE SLAG

OPLADEN

GEBRUIKEN

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VEILIGHEIDSINFORMATIE

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ATTACHMENTS

SAFETY INFORMATION

Read all safety information before using the Refectio. To reduce the risk of electric shock, fire, injury or other damage, the Refectio massage gun must be used according to the following safety instructions: This product is to be used by adults only. If in doubt about the product, consult a physiotherapist or doctor. The massage gun can be used responsibly for a maximum of 60 consecutive seconds in the same muscle group. Only use the massage gun on dry body surfaces, never on wet surfaces. The massage gun can be used over clothing. Only use the Refectio on the soft muscle tissue with the correct attachment. Never use the massage gun on the head, bones or sensitive body parts. Using the massage gun should not hurt or feel uncomfortable. Bruising can be the result of incorrect use of the massage gun. Stop immediately if a bruise develops.

Only use the attachments on the recommended muscle groups. Do not put objects or other things in the ventilation hole of the massage gun. Never allow the ventilation hole to come into contact with water or other liquid agents. Never drop the massage gun.

IMPORTANT NOTES

Only charge the massage gun with the original supplied adapter. View and examine the massage guns for any damage before use. Never change the parts inside the massage gun yourself. Never leave the massage gun unattended behind while charging. Only use this Refectio device with the permission of a doctor or physiotherapist in the following cases:

- Pregnancy
- Diabetes
- Other complications such as:
 - neuropathy, retinal wear, wear and tear pacemakers, recent surgery, epilepsy, migraine, hernia, spondylolisthesis, spondylo analysis, spondylosis.
- Recent joint replacements, IUDs, metal pins or plates.
- Concern about physical health.

These contraindications do not mean that you are unable to use a vibration or musculoskeletal device, but we recommend that you consult a doctor first. Research is ongoing into the effects of vibration exercise on specific medical conditions. This will likely result in a shortening of the list of contraindications as shown above. Practical experience has shown that integration of vibration exercises into one treatment plan is advised. This must be done on the basis of the advice and in the supervision of a doctor, specialist or physiotherapist.

Read all important notes before using the massage gun. The massage gun must not come into contact with liquids or heat sources. Keep the ventilation port free of dust and debris at all times. Do not remove stuck parts such as screws. Do not use the massage gun for longer than 1 hour in a row. Let the massage gun cool down at least 30 minutes before using again.

Important comments

Remove the charger from the battery after the massage gun is fully charged and is ready for use. The charger must also be removed during use from the charging port.

Product is only suitable for adults.

Consult a doctor or physiotherapist on how to use the massage gun to prevent (worse) injuries. Transport the massage gun only in the included storage bag.

GUARANTEE

If this product becomes defective within one year from the date of purchase due to a manufacturing defect in the material or if the device does not work properly, Refectio will repair this defect free of charge, replace it or find a suitable solution.

These are excluded from the guarantee:

- (a) Damage has been caused by accident, abuse, misuse or during transportation.
- (b) Units are subject to unauthorized repair.
- (c) Units not used in accordance with the safety instructions and the instructions for use of the Refectio.
- (d) Damages exceed the product cost
- (e) Deterioration of the delivered product due to abnormal storage / or protective conditions of the customer.
- (f) Failure to provide dated proof of purchase.

In some states and Countries it is not allowed to have a limitation of damage, so the foregoing limitation may not apply to you. This warranty guarantees specific legal rights, but other rights may vary from country to country and state to state.

 REFECTIO

***PRODUCT
INFORMATION***

CONTENT

1X REFECTIO
MASSAGE GUN

1X CHARGER

8X ATTACHMENTS

1X LITHIUM BATTERY

* WITH 1 USER MANUAL
AND STORAGE BAG

TOUCHSCREEN
LCD DISPLAY

8 DIFFERENT
ATTACHMENTS

POWERFUL MOTOR
WITH 30 SPEEDS

SILENCE TECHNOLOGY

LITHIUM BATTERY

TECHNICAL SPECIFICATIONS

PRODUCT FEATURES

| | |
|---------------|------------------------------|
| WEIGHT | 1.9kg |
| BATTERY | Rechargeable lithium battery |
| BATTERY LIFE | 6 hours |
| CHARGING TIME | 2 hours |
| CAPACITY | 2500mAH |
| VOLTAGE | 25.2V |

PRODUCT PERFORMANCE

| | |
|-----------------|--------------------------|
| VIBRATIONSTANDS | 30 |
| SPEED | 3200RPM |
| FREQUENCY | 53 vibrations per second |
| LOUDNESS | 45db |
| ATTACHMENTS | 8 |
| AMPLITUDE | 12mm |

LET'S START

TURN ON: BUTTON UNDER

Switch on the massage gun by pressing the button under the massage gun from "off" to "on".

Step 1

DISPLAY ON

Click on the red power button to vibrate. You can change the speed with + and -.

Step 2

START / STOP | VIBRATE

Higher position | Lower position

*TIP: switch to "off" when you are not using the massage gun. This extends the life of the battery.

PAIR ATTACHMENT

Attach the attachment to the massage gun by means of a rotating movement.

Step 3

CHECK ATTACHMENT

Check that the attachment is properly attached by gently pulling it. If the motor moves with you, the attachment is properly attached.

LET'S START

Step 4

YOU CAN NOW USE YOUR MASSAGE GUN!

CHARGING

Step 1

CHARGING GOES FAST

1. Connect the Refectio to the adapter.
2. The green light turns on.
3. The Refectio is fully charged after 2 hours.

Step 2

CHECK THE BATTERY

Turn on the massage gun. Is the right number set at "99"? Then the Refectio is fully charged.

LET'S START

USE MASSAGE GUN

Now you use the massage gun on your muscles. Slowly move back and forth and massage the length of the muscle.

Step 1

RESISTANCE

To get deeper into the muscle or connective tissue, you can press the massage gun a little harder against the muscle.

Step 2

YOU CAN NOW USE YOUR MASSAGE GUN!

ATTACHMENTS

DISCOVER ALL 8 ATTACHMENTS

Your Refectio comes with 8 unique attachments. Each attachment has a unique function to treat all muscle groups, connective tissue in the body.

LARGE BALL

Suitable for every muscle group. Works optimally with large muscle groups such as the legs, buttocks, shoulders and arms. Avoid the neck muscles. Recommended maximum position: up to 30

FLAT HEAD

Suitable for every muscle group. Works optimally with small muscle groups such as the neck muscles, arms, calves and tendons and tendon plates. Recommended maximum position: up to 30

THESE ARE THE MOST COMMONLY USED ATTACHMENTS.

ATTACHMENTS

AIR CUSHION HEAD

Suitable for flat muscles. Works optimally on the abdominal, back and chest muscles. Recommended maximum position: up to 20

SHOVEL HEAD

Suitable for arm and leg muscles. Works optimally for attachments, such as the calf muscle and the lower leg. Use in the length of the muscle. Recommended maximum position: up to 30

COLUMN HEAD

Suitable for the whole body. Works optimally in areas surrounding the joints and to massage deeper muscle tissue. Recommended maximum position: up to 15

ATTACHMENTS

U-FORK

Suitable for the back and arm muscles. Works optimally in the length over the back along the spine. And along the length of the inside and outside of the arm. Recommended maximum position: up to 20

THUMB HEAD

Suitable for the whole body. Works optimally in treating of trigger points in the deeper muscle tissue. Recommended maximum position: 15
Note: use this attachment with care.

BULLET HEAD

Suitable for the whole body. Works optimally in treating of trigger points in the deeper muscle tissue. Recommended maximum position: 15
Note: use this attachment with care.