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
# Congratulations on your new SANBO Acupressure mat

---

You are one step closer to improving your sleep, recovery, reducing stress and preventing and healing aches and pains.

With the SANBO Acupressure mat you can create a moment of peace where and when you want.

In this manual you can read what the acupressure mat is good for, among other things. And how best to use it. Thank you for your purchase and have fun with the SANBO Acupressure mat!

The logo for SANBO, featuring the brand name in a bold, white, rounded font with a light blue outline, set against a white background.

## QUICK START

# Start your session

1

### **Be prepared, it can hurt!**

If you are not yet used to lying on an acupressure mat, it can be experienced as painful, especially in the beginning.

This is normal and not bad! You will notice that your body will react positively to this after a couple of sessions.

3

### **Listen to your body carefully**

If you are just starting out with the acupressure mat, we recommend building up use slowly.

Always listen carefully to your body. Start with a few minutes first and extend it as soon as it feels right. For optimal results, we recommend a minimum of 10-20 minutes. Your session can be longer than 20 minutes if you feel like it.

2

### **Create a calm environment**

It is important to really see your session as a moment for yourself.

If you want to benefit from all the positive points of acupressure therapy, we recommend lying down in a quiet environment without any disturbing factors from outside.

**4**

## **Sit down and slowly roll on the mat**

Distribute the weight evenly. For the best results, use the mat directly on bare skin.

To get used to the spike mat, you can start by wearing a thin layer of clothing until you get used to it.

**6**

## **Enjoy the relaxation**

Within minutes, your body will warm up, circulation is activated and muscles relax.

Your body's self-healing ability is activated and waste products are removed faster and better.

Enjoy the relaxation as long as it feels good for you.

**5**

## **Focus on your breathing**

Experience what happens in your body. At first, you may experience a twinge of pain. Your body produces endorphins and oxytocin (two happiness hormones).

Hold on, within minutes this "pain stimulus" will subside and you will slowly become completely relaxed.

**7**

## **After care**

Always drink plenty of water after a session on the mat! We recommend drinking at least 1 liter of water. Do this spread out over 45-60 minutes.

# The history of acupressure

The doctrine of acupressure is ancient. There are scholars who say that acupressure was used in Chinese alternative medicine as early as 5,000 years BC. In Chinese history, people were very concerned with the cosmos and the forces and rhythms of the earth. According to Chinese medicine, life energy, called qi, flows through pathways in our body. These pathways are connected to each other and are called meridians. In the skin of your head, back, chest, arms and legs, there are a dozen to fourteen energy pathways that end in your fingertips and toes. Each energy pathway transmits energy to your main organ and other muscles and organs. The energy pathways each have different pressure points.

Physical and mental complaints result from a blockage in the meridians. Acupressure is a way to release blockages in the energy pathways by manipulating the pressure points and bringing life energy back into the pathways. The knowledge from Chinese medicine about relieving pain by applying pressure to certain body points was gathered and spread throughout Asia. In India, Fakirs then made beds out of a wooden board with nails. Fakirs were known as miracle workers who also often performed as magicians and walked on a bed of nails during a show. They also used the bed of nails outside of shows to sleep on.

According to the Fakirs, the pain felt by the brain causes a natural production of a painkiller called Endorphins. This is a happiness hormone that causes you to become relaxed which allowed them to lie on the bed of nails longer. In the 1980s, there was a Russian doctor who had researched the principle of a bed of nails. He came to the conclusion that a bed of nails not only produced endorphins, it also stimulated blood flow, lowered blood pressure and produced Oxytocin.

Science has long denied the existence of meridians. This slowly changed in 1978 after a study using radioactive tracing methods. After injecting radioactive tracing fluid into acupuncture points, a clear line could be observed. After further research, it was determined in 1992 that meridians do exist.

## What happens to your body?

**0 - 30 seconds.** Wow this hurts  
Should it be like this? This can't be good for a person....

**30 - 60 seconds.** Your body stimulates blood flow and you feel it getting warm. The real work can begin.

**1 - 3 minutes.** The heat and pain can feel like you've been burned by the sun. Focus on your breathing and don't let it drive you crazy!

**3 - 6 minutes.** A tingling / itching sensation may occur.

**6 - 10 minutes.** Your natural relaxation mechanism takes over. The relaxation and release of stress can finally begin.

**10 - 20 minutes.** The natural blood flow is stimulated, your muscles are relaxed, the stress in your body slowly flows away. A calming feeling is created. The body needs to be hydrated!

SANBO App

# The massage guide in your pocket

Get guidance for your  
session.

Read about information  
related to acupressure.

SANBO Care

# A team that sits ready especially for you!

Please contact our team if anything is unclear or if there are questions about the product.