# lumie<sup>®</sup> vitamin L<sup>™</sup>

# Safety

If the unit is damaged, please contact Lumie before use.

Keep away from water and damp.

For indoor use only.

The unit gets warm when the light is on so make sure that it is not covered.

Use only the mains power adaptor supplied.

Keep the power supply and cable away from heated surfaces. If either are damaged they should be replaced with a Lumie approved product.

If you need to disconnect the unit from the power supply, remove the power adaptor from the mains socket or remove the power cable from the back of the unit.

Children using the light, or playing near it, should be supervised.

Do not attempt to service or modify any part of the unit.

If the unit has been stored outside the recommended temperature range (5 °C to 35 °C), allow it to stand for at least 1 hour before applying power to avoid condensation or overheating.

Not to be used by or on individuals who are unconscious, anaesthetised or incapacitated e.g. under the influence of alcohol.

## Electromagnetic compatibility

This product conforms to the electromagnetic compatibility (EMC) standard for medical devices and is unlikely to cause interference affecting other electrical and electronic products in its vicinity. However, radio transmitting equipment, mobile phones and other electrical/electronic equipment can interfere with this product; if this happens simply move the product away from the source of interference.

For further detailed information on the electromagnetic environment in which the device can be operated refer to the website: [www.lumie.com/pages/e-m-c](http://www.lumie.com/pages/e-m-c).

## Disposal

At the end of its life, do not throw the unit out with normal domestic waste. Please take it to an official recycling point to be disposed of responsibly.

 SN	Serial number		For indoor use only		Double-insulated (Class II electrical)
	European Declaration of Conformity: Medical Device Class IIa		Do not place in household waste. Dispose of at a recycling centre.		Standby
	Manufacturer: Outside In (Cambridge) Ltd, 3 The Links, Bar Hill, Cambridge, CB23 8UD, UK.		Read instructions for use		Member of the 'Green Dot' recycling scheme
	Keep dry		Caution		Recyclable

# Precautions

Research shows that light therapy is very safe. However, please check with your doctor first if any of the following apply:

- You have a major eye problem or have had eye surgery;
- You have been advised to avoid bright light or are taking medication that may cause photosensitivity;
- You have been diagnosed with depression or other mood/sleep disorders. This is especially important if you are also taking medication for this, since light therapy can also affect your mood.

A few people experience side effects but these are generally mild (including headaches, eye strain, nausea and hyperactivity) and can usually be solved simply by switching off the light. Over the next few days, gradually build up the time you spend in front of it to find a treatment time that works best for you.

Some people complain of difficulty sleeping and often this is down to using the light too close to bedtime.

If you have any problems while using Vitamin L, or if you've been using the light for a week but haven't noticed any benefit, please contact Lumie or speak to your doctor.

**IMPORTANT:** the guidelines in this leaflet are not medical directions for treating any condition. If your doctor has suggested using light therapy, follow their instructions and discuss any effects with them.

# Introduction

Using Vitamin L will provide you with the bright light your body needs during the dark winter months. It can help to boost your energy levels, put you in a better mood and make you feel more awake<sup>1</sup>. Light therapy is safe and is proven to help alleviate the winter blues and SAD (Seasonal Affective Disorder)<sup>2</sup>.


Please read these instructions carefully to get the most out of your Vitamin L. 

1. Acta Psychiatr Scand. 2001 Apr;103(4):267-74.

2. Am J Psychiatry. 2005 Apr;162(4):656-62.

## Using Vitamin L

- Connect the mains power adaptor to the socket at the back of the unit.
- Plug in the unit and switch on at the mains.
- Press  to turn on/off.
- Slot the stand into the back of the unit. For extra stability, reverse the stand when switching between portrait and landscape use:

- Position Vitamin L at arm's length (about 50 cm) so the light is shining towards your face.

Although you don't need to stare at it, the light must reach your eyes to have any effect. This means you can use Vitamin L whilst you're working

at the computer, watching TV, reading, etc. Dark glasses, tinted lenses or closing your eyes will all limit the effects of light therapy.

You should notice positive benefits after using Vitamin L for three or four days in a row.

## When should I use it?

It's best to tackle the winter blues and SAD as soon as possible. Look out for early signs - often in September or October - and if you start to feel lethargic, anxious or irritable, start using Vitamin L.

Most people can skip a day here and there, especially when the weather is bright. As spring comes around you will probably start to cut back your light therapy; your body will tell you if it's too soon and if you notice your symptoms returning simply use your Vitamin L again for a few days.

Light therapy is proven to put you in a better mood, to boost alertness and to help you feel more awake. Use Vitamin L wherever and whenever you need to restore concentration and give yourself a mental boost.

## What time of day?

If you're using Vitamin L to help with the winter blues or SAD, bright light in the morning has been shown to be most effective so try that first. For some people, early evening light exposure is beneficial.

Whether you're using Vitamin L for winter blues, SAD or as a general pick-me-up, once you've found a pattern that's right for you, try to use it around the same time every day.

### **I find it difficult to get going in the morning**

Switch on as soon as you can. If you don't have much time in the morning, then you should aim for a top-up or full session in the afternoon.

### **I feel sleepy or fall asleep too early**

Use Vitamin L in the late afternoon/early evening. It's best to avoid bright light too close to bedtime - within about three hours - as you may find it difficult to get to sleep.

## How long?

Everyone's different, so position Vitamin L at arm's length (about 50 cm) and start with 30 minutes a day. Some people will find this is enough for effective SAD treatment. If you need to, increase your exposure each day until you find a treatment time that works for you.

Distance (cm)	Light received (lux)	Treatment time (minutes)
20	10,000	30
30	5,000	30–60
50	2,500	60–120

You can get your daily dose of Vitamin L all in one go or use it on-and-off throughout the day for a cumulative effect. If you're more comfortable with Vitamin L further away, then less light will reach your eyes and you'll need to use it for longer.

There's no harm in using the light for longer than suggested. If you find your eyes ache slightly, try using it further away or switch it off for a while.

## LED lifetime

The LEDs in this product have a typical life of 30,000 hours. Although it's unlikely, it's possible for individual LEDs to fail. This will make only a negligible difference to light output and you can continue to use the product.

## Cleaning

Clean only with a dry or lightly moistened cloth.

# Warranty

This unit is covered under warranty against any manufacturing defect for 3 years from the date of purchase. This is in addition to your statutory rights. If possible, use the original packaging (including internal pieces) to return your Vitamin L for service. Make sure the plug is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first. Please contact Lumie for further advice.

## Technical description

Subject to change without notice.

### **Lumie Vitamin L**

100-240 Vac, 50/60 Hz, 0.5A max power adaptor (supplied) 

Operating temperature: 5 °C to 35 °C

Storage/transport temperature range: -25 °C to +70 °C

Medical Device Class IIa  0120

## Contact us

Lumie is Europe's leading specialist in light therapy. Our products are based on the latest scientific studies and we focus solely on light therapy and its benefits. If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on +44 (0)1954 780500 or email [info@lumie.com](mailto:info@lumie.com).

Lumie, 3 The Links, Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, UK.

Lumie is a registered trademark of:

Outside In (Cambridge) Limited, incorporated in England and Wales.

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© Lumie 2019
Outside In (Cambridge) Limited, constituida en Inglaterra y Gales.