

## USER MANUAL

### Blackhead Remover Device



## Product Information

The blackhead remover device is a high-tech skin care device that uses a combination of different technologies to effectively remove blackheads. It is designed for use on the face, neck, and chest. The device is easy to use and does not require any special skills or training. It is suitable for all skin types and is a safe and effective way to remove blackheads. The device is made of high-quality materials and is designed to be durable and long-lasting. It is a great addition to your skincare routine and will help you achieve clearer, smoother skin.

### Features

- 1. **LED Light Therapy:** The device uses LED light therapy to stimulate collagen production and reduce inflammation. It is safe and effective for all skin types.
- 2. **Microcurrent Therapy:** The device uses microcurrent therapy to stimulate the skin's natural healing process and reduce the appearance of fine lines and wrinkles.
- 3. **High-Frequency Therapy:** The device uses high-frequency therapy to stimulate the skin's natural healing process and reduce the appearance of acne and blemishes.
- 4. **High-Power Plasma Therapy:** The device uses high-power plasma therapy to stimulate the skin's natural healing process and reduce the appearance of scars and wrinkles.
- 5. **Easy to Use:** The device is easy to use and does not require any special skills or training. It is suitable for all skin types and is a safe and effective way to remove blackheads.

## Suitable For All Skin Types



## Replacement Heads 7x



## Using Steps

1. **Preparation:** Cleanse your face thoroughly with a gentle cleanser. Pat your face dry with a clean towel.
2. **Device Setup:** Insert the replacement head into the device. Turn on the device and select the desired treatment mode.
3. **Treatment:** Gently glide the device over your skin. Do not press or rub the device against your skin. The device will automatically adjust the intensity of the treatment based on your skin type.
4. **Post-Treatment:** After the treatment is complete, apply a soothing toner or serum to your skin. Avoid using harsh skincare products for 24 hours.
5. **Frequency:** Use the device 2-3 times per week for best results. Avoid using the device on the same area of your face more than once per day.
6. **Precautions:** Do not use the device on broken, irritated, or sunburned skin. Avoid using the device on the eye area.

### Notes

The device is not intended to be used on the eye area. It is not intended to be used on the neck or chest area. It is not intended to be used on the face of children. It is not intended to be used on the face of pregnant women. It is not intended to be used on the face of people with a history of seizures. It is not intended to be used on the face of people with a history of heart disease. It is not intended to be used on the face of people with a history of high blood pressure. It is not intended to be used on the face of people with a history of diabetes. It is not intended to be used on the face of people with a history of kidney disease. It is not intended to be used on the face of people with a history of liver disease. It is not intended to be used on the face of people with a history of lung disease. It is not intended to be used on the face of people with a history of cancer. It is not intended to be used on the face of people with a history of HIV/AIDS. It is not intended to be used on the face of people with a history of tuberculosis. It is not intended to be used on the face of people with a history of syphilis. It is not intended to be used on the face of people with a history of hepatitis. It is not intended to be used on the face of people with a history of malaria. It is not intended to be used on the face of people with a history of dengue fever. It is not intended to be used on the face of people with a history of Zika virus. It is not intended to be used on the face of people with a history of Ebola virus. It is not intended to be used on the face of people with a history of COVID-19.

## FAQ

- Question 1: Can I use the device on my neck?**  
Answer: Yes, you can use the device on your neck. The device is designed to be used on the face, neck, and chest.
- Question 2: Why does my skin feel red after using the device?**  
Answer: The device uses LED light therapy, microcurrent therapy, high-frequency therapy, and high-power plasma therapy. These treatments can cause temporary redness and irritation. This is a normal reaction and will subside within 24 hours. You can use a soothing toner or serum to help reduce the redness.
- Question 3: Can I use the device on my eye area?**  
Answer: No, you should not use the device on your eye area. The device is not intended to be used on the eye area.
- Question 4: Can I use the device on my chest?**  
Answer: Yes, you can use the device on your chest. The device is designed to be used on the face, neck, and chest.
- Question 5: Why does the device stop working after a few days?**  
Answer: The device may stop working if the battery is low or if the device is not properly charged. Check the battery level and charge the device for 24 hours before using it again.

## Product Parameter

Model Name	Blackhead Remover Device
Material	ABS Plastic
Power	5W
Color	White
Weight	100g

## Packaging Accessories

