# *Parvum*

ENGLISH

DEUTSCH

NEDERLANDS

FRANÇAIS

POLSKI

## Rollator **SUPERLIGHT**

TABLE OF CONTENTS

Parts overview .....3

Indication.....4

Contraindication .....4

Safety instructions .....4

Height adjustment .....4

Braking and parking .....4

Adjusting the brake .....4

Sitting.....5

Folding .....5

Unfolding .....5

Cleaning and maintenance .....5

Assembly .....5

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Indicatie .....8

Contra-indicatie.....8

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Hoogteverstelling .....8

Remmen en parkeren .....8

Rem afstellen .....8

Zitten.....9

Opvouwen .....9

Uitvouwen .....9

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Montage .....9

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Indikation.....6

Kontraindikation.....6

Sicherheitshinweise .....6

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Sitzen .....7

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Aufklappen .....7

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Montage .....7

TABLE DES MATIÈRES

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Indication.....10

Contre-indication.....10

Consignes de sécurité .....10

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S'asseoir .....11

Pliage .....11

Déplier .....11

Nettoyage et entretien.....11

Montage .....11

SPIS TREŚCI

Części .....3

Wskazania .....12

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Siedzenie.....13

Składanie.....13

Rozkładanie.....13

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Montaż .....13

Parts overview/ Onderdelenoverzicht/ Teileübersicht/ Vue d'ensemble des pièces / Części

- (1) Handles (Grips) / Handvatten / Griffe / Poignées / Uchwyty
- (2) Brake levers / Remhendels / Bremshebel / Leviers de frein / Dźwignie hamulca
- (3) Backrest / Rugsteun / Rückenlehne / Dossier/ Oparcie
- (4) Seat surface (Seating area) / Zitvlak / Sitzfläche / Assise / Siedzisko

- (5) Bag / Tas / Tasche / Sac / Torba
- (6) Brakes / Remmen / Bremsen / Freins / Hamulce
- (7) Rear wheels / Achterwielen / Hinterräder / Roues arrière / Tylne koła
- (8) Front wheels / Voorwielen / Vorderräder / Roues avant / Przednie koła

Indication

This rollator is intended for use both indoors and outdoors to support people with limited mobility or improve their mobility. The purpose of this rollator is to assist the user’s mobility when moving around, for example by making it easier to walk further or to move around more easily.

Warning:

Users with limited hand and/or arm function may require special adjustment or modification of the handles or supports. Consult your healthcare professional. Rollators may not be suitable for users with balance disorders. This should be assessed individually by your healthcare professional.

The rollator is intended for use in normal, everyday activities only.

Contraindication:

The use of the rollator is not recommended for persons who are unable to sit or who have perception and/or balance disorders.

Preface

We are delighted that you have chosen this product from Parvum.

We recommend that you read these instructions carefully before using the product for the first time. These instructions contain essential safety instructions and useful advice on the correct use and maintenance of the product. If you have any questions or require additional information, please contact the store where you purchased this product.

Safety instructions

- Use the rollator exclusively as a walking aid.
- Ensure that the rollator is only used on flat, solid surfaces.
- The rollator is not intended to be used as a seat in (moving) vehicles such as buses or trains.
- Take into account the maximum permissible weight of the user. This information can be found in the technical specifications in this user manual or on the product sticker.
- Only use the rollator’s bag to transport objects and do not overload it. The maximum load can be found in the technical specifications in this user manual or on the product sticker.
- It is strictly forbidden to use the rollator on escalators and similar moving surfaces.
- Do not hang bags on the rollator, as this significantly increases the risk of tipping.
- When using the rollator on sloping roads, always use the walking brake by squeezing the brake lever.
- When sitting on the rollator, it is important to first set the brakes to the parking position so that the rollator does not roll away, causing you to fall.
- Do not lean too heavily on the handles, as this significantly increases the risk of the rollator tipping backwards.
- Be careful when unfolding or folding the rollator to prevent clothing or body parts from getting caught.
- Protect your rollator from prolonged exposure to sunlight, as certain parts may become hot and cause skin injuries when touched.
- Before each use, ensure that your rollator is in good condition and that all adjustable parts are securely fastened. Also check that the brakes are functioning properly.

- Regularly check the screws on the rollator and tighten them if necessary.
- The effectiveness of the brakes may decrease over time due to wear and tear on the wheels. Therefore, regularly check both the functioning of the brakes and the condition of the tires.
- The rollator must not be lifted by the seat, removable parts or brake cables.
- Any serious incident related to the use of the product should be reported to the manufacturer and to the competent authorities of the Member State where the user of this product is located.
- Always sit with your back to the front wheels when using the seat.

Height adjustment

Adjust the handles to wrist height. Have the user stand behind the rollator, upright with their arms relaxed at their sides. Then adjust the height of the handles so that they are at wrist height. Depending on the user, a higher or lower setting may be more comfortable. Consult your healthcare provider for the optimal setting.

Braking and parking

To brake while walking with the rollator, squeeze both brake levers upwards. Make sure you brake slowly and evenly with both brake levers at the same time.

To park or secure the rollator, press the brake levers down until you hear a click, which indicates that the parking brake has been activated. To release the parking brake, squeeze both brake levers until they are in the middle position and thus unlocked.

**Please note:** the effectiveness of the brakes may decrease over time due to wear and tear on the wheels. Therefore, regularly check both the functioning of the brakes and the condition of the tires.

Adjusting the brake

If the brake on the rollator is too weak or too powerful, you can adjust it by loosening or tightening the adjustment screw on the brake pad. If you are unable to do this yourself, please contact the shop where you purchased the rollator.

A local bicycle repair shop can often also assist with adjusting the brake.

The brake is correctly adjusted when the rollator brakes effectively without making any dragging noises when walking and braking with the rollator.

Sitting

Place the rollator on a flat, firm surface. Engage the parking brake by pushing the brake levers down. Important: Never sit on the seat if the brakes are not in the parking position. Make sure the rollator is stable and cannot slip away. When sitting down, sit upright and do not lean back. Always sit with your back facing the front wheels.

Folding

To fold the rollator, pull the cord on the seat. This unlocks the folding mechanism. You can now push the handles toward each other. The rollator is now folded.

Unfolding

Push the handles toward each other and push the seat down until you hear a click.

Cleaning and maintenance

Clean this product by hand every 3 months. Use a damp cloth for this.

Avoid using aggressive cleaning agents such as chlorine and vinegar for cleaning and disinfection.

Do not use abrasive cleaning agents.

Check the screws on the product every month and tighten them if necessary.

The product should not be exposed to temperatures above 60°C for more than one minute.

Check the condition of the tires every month. The tires are puncture-proof and do not require checking of tyre pressure or inflation. However, it is important to check them for cracks or other damage. In case of cracks/damage, contact the shop where you purchased this product.

It is advisable to have your rollator serviced at least once a year. Ask the shop where you purchased this product about the options.

In case of defects, please contact the shop where you purchased this product.

Assembly

The rollator is delivered fully assembled in the box, after which you can make it ready for use in a few simple steps. No tools are required. Instructions are provided below. A video tutorial is also available; please consult the shop where you purchased the rollator.

Remove the rollator from the box.

Slide the front forks into the frame, as shown in image 1. Make sure you attach the correct front fork on each side. You can tell which is which by the stickers on the front forks and the frame.

Image 1.

Attach the rear wheels by pushing the rear wheel tube into the frame at an angle and gradually pushing the tube into the frame. You will feel some resistance from the brake cable.

Image 2.

Install the back strap by attaching the plastic connectors to the handles.

Image 3.

You attach the cane holder by pushing the clamps firmly over the frame.

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#### Declaration of conformity

The manufacturer confirms that the rollator complies with EU Regulation 2017/745 for medical devices. A declaration of conformity can be requested from Parvum.com customer service.

#### Warranty

We supply this product with a warranty period of 24 months from the date of purchase. If you wish to make a warranty claim, you can use the invoice as proof of purchase and warranty. In the event of a material or manufacturing defect, we will replace parts free of charge.

Wear and tear on parts is generally not covered by this warranty, except when repair or replacement of these parts is necessary due to an obvious manufacturing or material defect.

We reserve the right to make changes based on technical improvements or design changes.

The warranty does not apply to the following parts: wheels (including tires, rims, bearings, etc.), brakes (brake levers, brake cables, brake pads, etc.), handles, backrest and seat.

#### Lifespan

The manufacturer estimates the lifespan of this product to be 4 years, provided it is used as intended and the maintenance and service instructions in this manual are followed.

The lifespan can be significantly extended if the product is handled, maintained, cared for and used with care within the technical limits of current science and technology.

Please note that intensive use or improper application or misuse can shorten the service life. This estimate of the products lifespan by the manufacturer does not imply any additional warranty.

#### Reuse

This product is intended for reuse. When passing the product on to a new user, ensure that you include the technical documentation, which is important for safe use, when passing this product on to a new user.

Before reuse, the product must be checked, cleaned and restored to good condition by the shop where you purchased it. Accessories for this product are available from the shop where you purchased this product.

#### Disposal

If you decide that you no longer need the product, please contact the shop where you purchased it. If the product has reached the end of its life, take it to a local waste collection point.

#### Contents

- Rollator
- Bag
- Backrest
- Cane holder
- Instructions for use

The rollator is delivered in a cardboard box. If possible, keep the cardboard box as it can be used for storage later. Upon receival, immediately check that the product is complete and intact.

Please note that the colour of the product shown may differ from the colour of your product. All colour variations of this model have the same technical specifications.

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Parvum B.V.
Ralon Medical And Science CO., LTD
SUNGO Europe B.V.