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
LED-gezichtsmasker - MJ01
LED face mask - MJ01
LED-Gesichtsmaske. - MJ01

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NL

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
Dear Customer,

Thank you for your order with **Maxi Meilleur®**

In this manual, you will find instructions on how to use and care for your product.

Please read the manual carefully and keep this manual for future use. If you have any questions, our team is always here to help, we'll be happy to answer them!

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**“GLOW WITH CONFIDENCE,
LIGHT UP YOUR BEAUTY”**

Maxi Meilleur® technology offers advanced innovation developed by experts. The mask is designed for home use as part of a beauty and self-care routine to support the appearance of the skin. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.

Model	MJ01
LED Lights	480 LEDS
Recommended session Time	10 to 15 minutes (user adjustable)
Charging Voltage	DC 5V  1,0A / 3,7V 3000mAh (Li-ion)
Net Weight	320g
Rated Power	4W
Product Size	Full set: 34 x 23 x 4.5 cm
Material	Certified, food-grade silicone

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
The Netherlands

Dr. Lelykade 44K, 2583 CM Den Haag

Package contents

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

Decolleté Nova


Eye protection


USB - C


Adapter


Attachment Straps


Remote control

User manual

Wavelengths

1. Red light (630nm):

Helps the skin look fuller, firmer and more refreshed. Ideal when the skin appears dull or fine lines seem more visible than desired.

2. Yellow light (590nm):

Yellow light is known for its calming presence in skin care routines. It supports a more even and balanced-looking complexion – a great choice when the skin needs balance.

3. Blue Light (460nm):

Supports skin that tends to look impure and needs clarity. Pleasant during periods when the skin feels unbalanced and in need of brightness.

4. NIR Infrared Light (850nm)


Near infrared is invisible to the eye but widely used in light applications. It supports a relaxed skin feeling and a softer-looking appearance.

First-Time Setup



Carefully remove all items from the box and check if all parts are included. You can refer to the “**Package contents**” section for the complete list.

First-Time Setup



1. Remove the seal on the back

2. Attach the eye protection to the inside of the mask

3. Attach the correct straps to the mask.

Remote Control & Charging

Charging

1. Make sure the remote control is not connected to the LED mask.

Warning: Do not use while charging.


2. Insert the charging cable into the bottom of the remote control and connect it to the supplied adapter.

3. Plug the adapter into the power outlet.

Warning: Use only the supplied adapter.

Remote Control

1. **Power Button:** Press and hold for 2–3 seconds.
2. **Intensity:** Adjust the brightness (intensity) of the lights by pressing 'E'
3. **Modes:** Use the "M" button to select the different light modes.
4. **Timer:** Use the "T+" button to add 5 minutes and use 'T-' for minus 5 minutes.
5. **NIR:** Choose On/Off for Near Infrared. Can be used with any light.

Instructions for use

The 4 Steps for Use

1. Cleanse Your Skin

Rinse your face to remove any makeup, sunscreen, or skincare products, and dry your skin.

2. Set up your treatment

Always inspect the device, remote control, and power cables for damage before each use. Do not use if it is damaged, malfunctioning, or if the unit has been exposed to water. Start with a lower intensity. Gradually increase the setting, especially if you have sensitive skin.

3. Disinfect After Use

Always disconnect the power adapter before cleaning the mask. Disinfect the mask after every session and remove oil and residue to keep it as hygienic as possible and maintain the quality of the LED therapy. Do not use abrasive, harsh, or solvent-based cleaning agents. After cleaning with alcohol, wipe the mask thoroughly with a damp cloth to remove any residue.

For stains: You can gently remove stains with 70% isopropyl alcohol.

4. After care

If desired, apply a hydrating or nourishing cream after the treatment.

Frequently Asked Questions

How often should I use the LED mask?

We recommend use 3 to a maximum of 5 times per week. Read the full instructions carefully before use. Use the build-up schedule for our best practices.

How do I gradually build up usage?

We recommend the following if you are using LED therapy for the first time: During the first week, use it for a maximum of 3 days for 10 minutes on the lowest setting. In the second week, use it every 2 days for 10 minutes on the lowest setting, and in the third week, use it every 2 days for 10 minutes on the medium setting. If the skin shows no irritation, you can gradually increase the frequency. If irritation does occur, taper off and then gradually increase the frequency again so that your skin can get used to the LED therapy.

The device won't turn on

- Press and hold the ON/OFF button for 3 seconds
- The device may need charging
- The cable might be in the wrong port

How long does it take to see results?

With consistent use, you can expect visible improvement within four weeks. Many clients notice the first changes after just two weeks: calmer skin and a healthier glow. It is not about rapid transformations, but about gradual improvement.

What if I skip a week?

Skipping a week does not mean you lose all your results. The positive effects gradually diminish with longer breaks, but your skin does not revert to zero. LED therapy works best as a regular part of your weekly routine.

Still have questions? Don't hesitate to contact us via email, we're happy to help!

General Safety

- Do not use the product if any parts are damaged.
- Do not disassemble or repair the product yourself.
- This product is not fire-resistant. Keep it away from fire and corrosive liquids.
- Do not immerse the product in water or other liquids.
- Do not use this device in a bathtub, shower, or water-filled reservoir. If the device is used in a bathroom, unplug it after use, as the proximity of water presents a hazard even when the device is switched off.
- Store the product in a place out of reach of children.
- Store the product in a dry, clean place away from moisture and direct sunlight.
- Always unplug the product after use and before storing it.
- Always follow the instructions in this manual. Improper use may lead to malfunctions or accidents.
- Do not use the product for purposes other than those described in the manual.
- Protect your eyes during use. Wear the provided eye protection and keep your eyes closed during the session.
- If you experience severe discomfort while using the product, stop using it immediately and consult a doctor.
- Keep the product out of reach of babies and children.
- Always perform the LED light sensitivity test before use to ensure compatibility with your skin.
- This device has been tested in accordance with EN 62471 for photobiological safety and classified as Exempt Group. We recommend keeping your eyes closed during use. For sensitive users and when using the blue (460 nm) program, we strongly advise using the provided eye protection.
- Do not use the product if you are under 18 years of age unless supervised by an adult.
- Avoid using this product on the same day as products containing photosensitizing ingredients.
- Do not use the product if you are taking photosensitizing medication (photosensitizers), if you have a genetic eye condition, and/or if you are sensitive to light. If in doubt, consult your doctor or pharmacist.
- Do not use the product if you have severe skin reactions to light or if you experience discomfort.

- Photosensitive conditions: If you have a condition that causes sensitivity to light, do not use this product.
- Do not use the product if you are pregnant or breastfeeding.
- This product is not suitable for children.
- The product is not intended for use in rainy, wet, or damp environments.
- The product is not intended for professional, industrial, or commercial use.
- If you have a physical, sensory, or mental condition, or if you are undergoing treatment for dermatological conditions (such as lupus erythematosus, photosensitive eczema, or albinism), consult a doctor before use.
- Note: Always consult a specialist before using the product if you have an underlying medical condition.

Warranty

We offer a 12-month warranty from the date of purchase for this product. This warranty covers defects in materials or workmanship under normal use.

To claim your warranty, you must keep your proof of purchase.

If the seal has been removed, the product may no longer be eligible for return for hygiene reasons. If you experience any issues during the return or warranty period, we recommend contacting our customer service team.

Covered by Warranty: Material or manufacturing defects & functional issues caused by manufacturing defects

Not Covered by Warranty: Damage caused by misuse, improper care, or accidents, Normal wear and tear, Removal of the seal during the return period, Unauthorized repairs or modifications carried out by yourself

The warranty does not cover damage caused by overheating or any other use outside the intended purpose as described in the manual.

Hygiene seal

You may open the outer packaging and review the product information without losing your right of withdrawal. For reasons of hygiene and health protection, the Maxi Meilleur LED face mask is packaged with a double seal: an outer seal on the box and a separate hygiene seal on the product itself. Because the device comes into direct contact with the skin of the face during use and the silicone surface can no longer be restored to a saleable condition in a hygienically responsible manner once the seal has been removed, the right of withdrawal expires as soon as the hygiene seal on the product is broken. Breaking this seal therefore results in the loss of the right of withdrawal.

Battery Safety (Lithium-ion)

This product contains a built-in rechargeable lithium-ion battery. Read the following warnings carefully to prevent the risk of fire, explosion, leakage, or personal injury.

Use and handling

- Never attempt to open, pierce, crush, heat, disassemble, or modify the battery. Doing so may cause fire, explosion, or leakage of toxic substances.
- Never expose the battery to open flames, sparks, or high temperatures. Never dispose of the product by burning.
- Do not use the product if it has been dropped, is visibly damaged, or if you notice swelling, deformation, liquid leakage, unusual odour, or excessive heat. Stop using the product immediately and contact customer service.
- Avoid short-circuiting the battery contacts with metal objects.
- If battery fluid comes into contact with skin or eyes, rinse immediately with plenty of water and seek medical attention.

Charging

- Use only the original power adapter and USB-C cable supplied with the product. Using other adapters may cause overheating, battery damage, or fire.
- Only charge the device at an ambient temperature between 0 °C and 40 °C (32 °F and 104 °F). Charging outside this temperature range may permanently damage the battery or create a safety hazard.
- Do not charge the product in direct sunlight, near a heat source, in a closed vehicle, or on a soft surface such as a bed or sofa.
- Do not leave the device unattended for extended periods while charging, and do not charge it while you are sleeping.
- Disconnect the adapter as soon as the battery is fully charged. Do not overcharge the product.
- Do not use the product while charging.

Storage

- Store the product at an ambient temperature between 10 °C and 30 °C (50 °F and 86 °F) in a dry, clean place, away from direct sunlight, moisture, and heat sources.
- For long-term storage (longer than 3 months): charge the battery to approximately 40-60 % and recharge every 3 to 6 months to prevent deep discharge.
- Avoid storage at temperatures below 0 °C or above 45 °C (32 °F or 113 °F). Extreme temperatures significantly reduce battery life.

Battery Disposal

- This product contains a rechargeable lithium-ion battery (Li-ion) (3.7 V, 3000 mAh).
- Do not attempt to remove the battery yourself. The battery is built in and may only be removed by qualified personnel at an authorised collection point.
- At the end of its service life, return the complete product to an authorised collection point for batteries and electrical equipment. Never dispose of the product or battery with regular household waste.
- A short-circuited, damaged, or leaking battery can cause fire. When disposing, cover the contacts with insulating tape.

WEEE

Separate Collection of Electrical Equipment

The crossed-out trash can symbol (WEEE/AEEA mark) means that electrical and electronic equipment must not be disposed of with residual waste. This product falls under the European WEEE Directive (2012/19/EU). Do not dispose of this product with regular household waste. Dispose of the product in accordance with the applicable collection and recycling regulations. Use designated collection points for batteries. Do not incinerate the product. For information about disposal, contact the manufacturer or the EU Responsible Person listed on the product or packaging. This product complies with the applicable EU safety regulations. The EU Declaration of Conformity is available upon request.

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Maxi Meilleur HealthTech B.V.