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Electric Oral Dental Cleaner Set User Manual


Precautions

1. This product is for indoor use only.
2. The product should not be stored in a place where direct sunlight can reach. It should be placed in a dry or good-ventilated space.
3. This product is designed only for adults due to sharp accessories. Please place it out of reach of children and animals to avoid any accidents.
4. This product is equipped with LED light. Please avoid direct exposure to the eyes for a long time.
5. This product is only for cleaning teeth, gums and tongue. Usage for other purposes is prohibited.
6. If you have any dental problems, you should directly consult your dentist before using this product.
7. This product is only a personal cleaning product. It is not recommended to share this with another person.
8. For personal hygiene, please regularly replace the toothbrush head every 3 months.

Specifications

Working modes: 3 (low, medium and high)
 Charging mode: DC 2.0 charging cable USB
 Charging voltage: DC 5V
 Rated voltage: 3.7V
 Rated current: 70 - 130ma
 Motor speed: 11500rpm
 Waterproof grade: IPX6
 Battery type: Polymer 602040
 Battery capacity: 400mAh
 Product size: 26.9 x 26.9 x 224 (mm)

Description

Instructions

1. Charge the device for 3 hours for first-time use.
2. When installing the toothbrush head, align the brush head so that the bristles face the handle and insert the brush head into the transmission shaft.
3. The default is low for the first time, the key is switched in 10 seconds and the button is off after 10 seconds.
4. If you need to turn on the LED light, press the light button. Press it again to turn off the light.

Features

1. Equipped with three dental cleaning heads. Easy to operate and suitable for multiple oral dental cleaning purposes according to personal needs:
 Dental pick: Place pick at the gum line between teeth and gently pull the device away from the gum line. Use it weekly. Bristle toothbrush: Clean the teeth.
 Tongue scraper: Gently glide along your tongue from back to front and rinse. Daily cleaning of the tongue surface will help promote fresher breath. This improves taste sensitivity and a cleaner healthier mouth.
2. LED spotlight: White illumination makes mouth cleaning much more easier.
3. Ergonomic and aesthetic design - integrated induction of silicone case body with leather imitation effect. Non-slip handle with anti-fog protection.
4. The product is lightweight with storage box for 3 head parts. The cleaning head can be changed at will for convenient storage.
5. Battery is durable and has a long time usage.

Charging instructions

1. Take out the charging cable from the package. Insert the round probe of the cable into the DC charging jack. The other end belongs into the USB cable adaptor of 5V and 1A.
2. During charging, the white light indicator is flashing. The white light is keeping bright on when it's fully charged. When the battery is low, the indicator is continuously blinking to remind it needs a recharge.

Replace oral cleaning heads

This product is equipped with 3 cleaning heads as a kit. When you change the head, you should hold the ring-part of the cleaning head and pull out vertically shown as below:

Common troubleshooting problems

Malfunction	Cause	Solution
Weak swing intensity	Press too hard when using	Lightly place the brush head on your teeth
LED light brightness is weak	Low battery	Please charge
Toothbrush doesn't work properly	The toothbrush body is low in battery	Please charge. If it is still not working after fully charging, please consult the retailer