EN 1

DA 13

DE 25

FI 40

FR 53

NL 68


NO 81

PL 93

SE 107

sunbooster

POWERED BY SUNLED TECHNOLOGY

SLS2000

Content

- 1x SunBooster SLS2000
- 1x USB-C cable
- 1x User manual

Manufacturer

SunLED Life Science Research B.V.

Science Park 106

1098 XG Amsterdam

The Netherlands

Email: info@sunled.health

Website: www.sunbooster.health

Intended purpose

This screen-mounted wellness device is intended to energize and improve the user's mood and physical well-being. It is not a toy and must not be used by children. The product is designed for indoor use only and must be used in accordance with the instructions provided in this manual. Any other use is considered improper. The manufacturer is not liable for damages resulting from misuse or use contrary to these instructions.

Intended users



- This product may be used by persons with reduced physical, sensory, or mental capabilities, or by persons lacking experience and knowledge, provided that they have been given supervision or instruction concerning the safe use of the product and understand the hazards involved.
- Children shall not play with the product, and cleaning shall not be carried out by children without supervision.
- This product shall not be used by:
 - individuals with sunburn;
 - individuals under medical supervision for conditions associated with photosensitivity;
 - individuals undergoing treatment with photosensitising medication.

Unpacking

1. Remove and properly dispose of all packaging material.
2. Make sure that all parts are present and undamaged.

Overview

1. Near-infrared LED lights
2. Display
3. TOF sensor
4. Mount
5. Return button
6. Previous and next buttons
7. Confirm button

Safety

Before use

- Thoroughly read this user manual before use.
- Follow all safety instructions.
- Keep this manual for future reference.
- This manual is part of the product and must be given to the new owner upon resale or exchange.

Operation and environment


- For indoor use only. 
- Protect the product from contact with water and moisture. Do not use the product in bathrooms or other high-humidity areas.
- The product may feel warm to the touch.
- Ensure the product is safely mounted to a monitor or laptop screen using the provided mount.
- This product is intended to be powered directly through a USB-C port. If an external power adapter is used, it must comply with the applicable regulatory certifications EN 61558-2-16 or EN 60335-2-29.
- This product is a class I consumer laser product and contains a Time-of-Flight (TOF) distance sensor. The sensor uses a Class I laser, which is safe under all conditions of normal use. The laser component is compliant with EN 50689.
- Only use the product with sufficient ambient light.

- Only use the product at a minimum distance of 30 cm.
- If unexpected side effects occur within 48 hours after first use, stop using the product and seek medical advice before further exposure.

Maintenance and service

- Do not use the product if it is damaged.
- Clean the product only as described in this manual.
- Never attempt to open, disassemble, or service this device. Unauthorized access may result in hazardous situations.
- Maintenance or repair may only be performed by a professional electrician or the manufacturer.
- The battery is not replaceable.
- Disposing of the battery in a fire or heater, subjecting it to mechanical crushing or cutting, leaving the battery in an extremely high temperature environment, or subjecting the battery to extremely low air pressure, may result in an explosion, leakage of flammable liquid or gas or other hazardous conditions.
- The light source of this product is not replaceable; when it reaches its end of life the entire product must be replaced.

Using the SunBooster for the first time

1. Connect the device to a USB-C port, such as a laptop, monitor, etc.
2. Use ⏪ to select a language and press ⏹ to confirm the language setting and go to date settings.
3. Set the correct date and time using the ⏪ buttons and the ⏹ button to confirm.
4. An overview of the settings is shown, press ⏹ to confirm and go to the start screen.

Note: Press ⏮ during set-up to return to the previous settings screen.

Turning the device on/off

- Press and hold ⏹ for 3 seconds to turn the device on.
- Press and hold ⏹ for 3 seconds to turn the device off.

Note: If the device is unplugged without being turned off, it will automatically power on and continue when reconnected to the power source.

Starting a cycle

1. Connect the device to a USB-C port, such as a laptop, monitor, etc.
2. The start screen is displayed. Press ⌂ to start a new cycle. A progress bar displays the cycle's progress.

Note: If a cycle is started between 20:00 and 04:00, the device displays a late-start warning. Press ⌂ to start a cycle.

Pausing a cycle

- Press ⌂ to pause the cycle.
- Press ⌂ again to resume the cycle.
- If the cycle is paused for more than 5 minutes, the device enters sleep mode and the display turns off.
- Press any button to exit sleep mode and resume the cycle.

Progression of a cycle

The screen displays the cycle's progress as a percentage. Each cycle lasts 2 to 4 hours.

Note: One cycle per person per day is sufficient, although additional cycles do not pose any risk of harm.


Resetting a cycle

1. To reset the current cycle, press ⌂.
2. Press ⌂ to confirm and return to the start screen. If no action is taken within 10 seconds, the cycle resumes automatically.

Distance to the product

The recommended operating distance is between 30 and 80 cm. The duration of a complete cycle varies with the distance to the SunBooster: the closer the user is seated, the shorter the cycle time.

1. < 30 cm: the near-infrared lights will be switched off.
2. > 100 cm: pause mode will be activated.

Note: when no user is detected, the device automatically enters sleep mode. When a user is detected again, the device automatically reactivates.

End of cycle

After a cycle is completed, the weekly cycle progress is shown on the display. The letters at the top indicate the days of the week. The yellow bar shows progress: a full bar represents one completed cycle.

To start a new cycle:

1. Press Ⓢ to close the weekly cycle progress screen.
2. Press Ⓢ again to start a new cycle.

Note: If no action is taken within 3 minutes, the device enters sleep mode and the display turns off. Press any button to exit sleep mode and resume.

Factory reset

Press and hold Ⓢ and Ⓢ simultaneously for 10 seconds to perform a factory reset.

Note: Changing the language, time, or date settings requires performing a factory reset. This will erase all stored progress data.

Cleaning and storage

- Clean the product with a dry cloth, never clean it with water!
- Never use any abrasive, corrosive, or scouring cleaning agents, as these may damage the product.
- Store the product in a clean, dry place inside the original packaging.

Troubleshooting

- The device cannot operate in dark environments. If the message “Environment too dark” appears, ensure that the device is placed in a well-lit environment.
- If the device overheats, the near-infrared lights are disabled and the device enters sleep mode. Allow the device to cool down before resuming use.
- If you move within 30 cm of the device, the near-infrared lights are disabled and the message “You’re too close” appears. Ensure the device is not positioned too close to obstacles such as walls.
- The device may not function with certain USB ports. If this occurs, use a USB charger connected to a power outlet. The adapter specifications can be found in the section “Specifications”.
- If a webcam without an infrared filter is used while the device is active, its colours may appear distorted.
- The device may malfunction if it shares its wavelength with other devices (e.g. wireless infrared headsets).
- If an error screen (E + number) appears on the display, please contact your retailer or visit www.sunbooster.health.

Software updates

Go to www.sunbooster.health for further information.

Warranty

See the warranty card for more information on product warranty.

Disposal

This product is subject to European Directive 2012/19/EU. The product (including the USB-C cable) and packaging must not be disposed of with other household waste. The product must be handed in at an authorized collection point or recycling centre.

Specifications

Model:	SunBooster SLS2000
Working distance:	30 – 80 cm (optimal); 81 – 100 cm (reduced efficiency)
Dimensions:	83 × 62 × 56 mm
Operating conditions:	0 – 30 °C at 10 – 95 % RH
Storage conditions:	-40 to 70 °C at 10 – 95 % RH
Weight:	120 g (without packaging and cable)
Max. power consumption:	4,5 W
Cable requirements:	USB cable rated for ≥ 1 A
USB power adapter requirements:	USB-C power delivery 10 W (DC 5 V / 2 A)
Wavelengths:	850 nm (LED) 940 nm (TOF sensor)

EC Declaration of Conformity

SunLED Life Science Research B.V. declares that the above mentioned apparatus fulfils all relevant provisions of

- EMC Directive 2014/30/EU
- RoHS Directive 2011/65/EU & Delegated Directive (EU) 2015/863
- Battery Regulation 2023/1542

In conjunction with the following harmonised standards and where appropriate other technical standards and specifications:

EN ISO 12100:2010; EN 55032:2015/A11:2020; EN 55035:2017/A11:2020;
EN 61000-3-2:2014; EN 61000-3-3:2013; EN 61000-6-1:2007;
EN 61000-6-3:2007/A1:2011; EN IEC 63000:2018

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SunLED Life Science Research B.V. försäkrar att den ovan nämnda apparaten uppfyller alla relevanta bestämmelser i