KLEGU

MASSAGE MAT

MASSAGE MAT

Model: HY-8639

NL · EN

Inhoudsopgave / Table of contents

Nederlands (NL)

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EN – English

This user manual accompanies the **KLEGU** massage mat and contains important information about safety, disposal and use. Read the manual carefully and follow the instructions and warnings to ensure your own safety and the safety of others.

Keep the manual for further reference, even if you are familiar with the product. Make sure to include this manual when passing on this product to third parties.

The images in the manual are for illustrative purposes only and may differ slightly from the product. In that case, KLEGU is not obliged to adjust the manual. Use the actual product as a reference.

1. INTRODUCTION

Thank you for choosing a **KLEGU** product. Read the user manual carefully before use to ensure proper and safe product use. Always follow the instructions to prevent damage or harm.

Any questions or concerns regarding the product, please feel free to contact us by email: Info@klegu.nl.

We hope you enjoy your purchase!

1.1 Certification overview

Our product complies with the highest international safety and quality standards and holds the following certifications:

CE mark: The CE mark (an acronym for the French “conformité européenne”) certifies that a product meets EU health, safety, and environmental requirements, guaranteeing consumer safety.

This certification guarantees that our product is safe, reliable, and compliant with all relevant laws and regulations.

1.2 Declaration of CE conformity

This product complies with the essential requirements and provisions of the following directives:

- 2014/35/EU (Low Voltage Directive - LVD)
- 2014/30/EU (EMC Directive)
- 2011/65/EU (Restriction of Hazardous Substances - RoHS)

The full EU Declaration of Conformity is available upon request from the supplier.

2. EXPLANATION OF SYMBOLS

Throughout this user manual, the following symbols can be used:

INFO: Provides important information or helpful tips to enhance the use and performance of the product.

DANGER: Indicates a hazardous situation that, if not avoided, will result in serious injury, death, or significant property damage.

WARNING: Indicates a potentially dangerous situation that could result in serious injury or

material damage.

CAUTION: Indicates a situation that could result in minor injury or material damage.

3. SAFETY INSTRUCTIONS

WARNING: **READ** the user instructions carefully and especially observe the safety information. If you do not follow the safety instructions and information on proper handling in this manual, we assume no liability for any resulting personal injury or damage to property. Such cases will invalidate the warranty.

Electrical devices can expose people, especially children, to hazards that can cause serious injury or death. These safety instructions are intended to give you the information you need to use the KLEGU massage mat and avoid such hazards. Common sense and respect for the device will help reduce the risk of injury.

All possible hazards cannot be covered here, but we have tried to highlight the important ones. Also, observe the warnings or safety instructions in the individual chapters of this manual.

Contact us by e-mail at Info@klegu.nl if you have questions that remain unanswered by these warnings & safety instructions.

3.1 General device safety

- Do not use the massage mat for any purpose other than its intended use (see the chapter "INTENDED USE").
- The massage mat is designed and intended for private use and is not suitable for commercial or industrial use.
- The massage mat is intended for indoor use only. Do not use the device outdoors.
- Handle the device carefully. Jolts, impacts, or falls can damage the device.
- Place the massage mat on a flat, firm and stable surface. Do not use it on an unstable or slippery surface where the mat could slide away.
- Do not place the massage mat under excessive mechanical stress. Do not stand or jump on the mat or exert point pressure on it.
- Do not fold or kink the massage mat sharply. This can damage the internal massage mechanisms, heating elements and air bags.
- The device must be properly and securely connected. Follow the instructions to prevent any possible damage/injury.
- Never use accessories not recommended by the manufacturer. They can pose a danger to the user and damage the device.
- Never insert any object or fingers into any opening or between the moving massage parts.
- Do not apply excessive force when using the buttons on the remote control.
- Do not rest any body part with your full weight on a single massage point for an extended period.
- Ensure that the moving parts can always move freely and that no objects, clothing or hair are jammed.
- Do not leave the switched-on massage mat unattended. Close supervision is necessary when using near small children.
- Do not use the massage mat in dusty environments.

- Do not use the massage mat when the fabric cover is torn or missing.
- If it is no longer possible to operate the massage mat safely, take it out of operation and protect it from accidental use. A safe operation can no longer be guaranteed if the device:
 - is visibly damaged,
 - is no longer working properly,
 - has been stored for extended periods in poor ambient conditions,
 - or has been subjected to serious transport-related stresses.

You can take the device out of operation by removing the power cord or by cutting off the mains plug of the power cord.

3.2 Persons & use

- The massage mat is not a toy. Keep and store it out of the reach of young children and pets. Don't let them climb or sit on the device.
- Do not use the massage mat on pets.
- Children under the age of 18 should use the device only under the constant supervision of an adult.
- The massage mat can be used by children aged 14 years and above and persons with reduced physical, sensory, or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning the safe use of the device, and they understand the resulting hazards.
- Cleaning and maintenance must not be performed by children.
- Do not fall asleep while using the unit.
- Do not use the massage mat longer than the recommended duration per session and not on the same body area for a long time.
- Do not use the massage mat if you are intoxicated or feeling unwell.

3.3 Electronics

- You may not reconstruct, self-repair, or disassemble the massage mat. This can damage the device.
- Maintenance, modifications, and repairs may only be carried out by the manufacturer, a professional workshop appointed by the manufacturer, or a similarly qualified person.
- Check the device from time to time for damage. Never use the device if it shows signs of damage.
- Do not use any other mains plug or power adapter for this device. This can cause failure and/or damage.
- Ensure your mains socket supplies the correct power voltage to avoid damage to the device.
- The mains socket must be located near the device and be easily accessible.
- To avoid energy loss and potential safety risks, always unplug the mains plug when the device is not in use, during cleaning or maintenance, or if a fault occurs.
- Do not turn the massage mat on and off by plugging and unplugging it. Always use the ON/OFF button on the remote control.
- Never force the mains plug into a mains socket; if it does not fit easily, stop using it. Always ensure the plug is properly connected.
- If the power fails during operation, turn off the device immediately and pull the mains plug from the mains socket.

- If smoke is produced, or there are any unusual sounds or smells, switch off the massage mat immediately and remove the mains plug. The mat should not be used until it has been inspected. Never inhale smoke from a potential device fire; if you have inhaled smoke, consult a doctor.
- Do not extinguish with water in the event of a fire. Smother the flames with a fireproof blanket or a suitable fire extinguisher.
- For safety reasons, always unplug the device in case of a thunderstorm.

3.4 Power cord

- Fully unwind the power cord before use.
- Do not modify the power cord.
- Check the power cord periodically for damage. Never use the power cord if there are any signs of damage.
- Do not touch the power cord if it is damaged. A damaged power cord can cause a deadly electric shock. First power down the respective mains socket, then carefully pull the mains plug from the socket. Never use the device if the power cord is damaged.
- To prevent risks, damaged power cords may only be replaced by the manufacturer, a professional workshop appointed by the manufacturer, or a similarly qualified person.
- Never pull the mains plug from the socket by pulling the power cord. Always pull it using the intended grips.
- Do not carry or pull the massage mat by its power cord.
- Never wrap the power cord tightly around the device. This will cause the cord to wear prematurely and break.
- Ensure that you do not trip over the power cord.
- Do not let the power cord hang over the edge of a table or counter.
- The power cord should not be kinked, trapped, or twisted.
- Ensure that the power cord is not squeezed, bent, damaged by sharp edges, or put under mechanical stress.
- Do not position the power cord where it can be pinched or crushed by doors.
- Never tie a knot in the power cord.
- Avoid excessive thermal stress on the power cord from extreme heat or cold.
- If using an extension cord, ensure that the total power consumption of connected devices does not exceed the rated capacity.
- Children can become entangled while playing with cords, posing a strangulation hazard. Cords are not toys.

3.5 Heat & cooling

- The massage mat may become warm during operation. This is normal; the heating zones (neck and lumbar) are intended to warm up.
- Do not cover the device while it is in operation, as this may cause it to overheat. This poses a risk of fire, electric shock, or injury.
- Ensure that the ventilation around the device is good during use.
- Keep the device away from heat sources and from areas where the ambient temperature is relatively high. Never use the device near naked flames or on hot surfaces.
- The massage mat is designed to work best in ambient temperatures between 0°C and 35°C and to be stored at temperatures between 0°C and 45°C. The device may be damaged if

used or stored outside this range. Do not subject your device to dramatic changes in temperature. Do not store it in a hot car or direct sunlight for extended periods.

WARNING: Do not use the heating function unattended and do not fall asleep with the heating switched on. Prolonged contact with heated surfaces can cause skin irritation or low-temperature burns, even at moderate temperatures, especially in persons with reduced heat sensitivity.

3.6 Water

WARNING: This device is not waterproof.

- Ensure your skin is completely dry before using the massage mat. Using the device with wet skin may result in an electric shock or malfunction.
- Never touch the plug with moist hands while connecting or disconnecting it from the mains socket. Moisture can cause electric shock or short circuits.
- Do not use the massage mat in or near a bath, shower, swimming pool, or any container filled with water.
- Do not use in wet environments, and do not submerge in liquid. Never pour liquids over or near the massage mat. If liquid does get into the housing, unplug it immediately.
- Do not reach for the device if it has fallen into the water. Immediately unplug it and only then remove it from the water.
- Ensure that the device is placed on a dry surface.

3.7 Medical considerations

- Do not use this device together with other therapeutic equipment.
- This device is designed for comfort and relaxation only, not for medical treatment.
- Those with any of the following conditions should use this device with caution and follow their doctor's advice:
 - Individuals with back, neck or spinal conditions (e.g., a herniated/slipped disc or spinal stenosis).
 - Individuals who have recently had surgery on the back, neck or spine.
 - Individuals with whiplash or a recent neck or back injury.
 - Individuals with impaired nerves or reduced skin sensitivity.
 - Individuals with metal, silicone, or plastic implants due to fractures or plastic surgery.
 - Individuals who are not sensitive to heat or have lost their ability to sense temperature.
 - Individuals with back, neck or spinal complications resulting from an accident.
 - Individuals with a history of or current severe conditions (e.g., tendon injuries, arthritic deformities, fractures, or varicose veins).
 - Individuals with edema, ischemic vascular disease, or severe peripheral neuropathy.
 - Individuals with bleeding disorders or coagulation problems.
 - Individuals with malignant tumours or cancer.
 - Individuals with severe heart disease or abnormal blood pressure.
 - Pregnant or postpartum women, especially those with complications in the first trimester or shortly after childbirth.
 - Individuals with bone diseases or fragile or injured bones.
 - Individuals with osteoporosis.
 - Individuals with impaired consciousness.

- Individuals with diabetes-related sensory impairment.
- Individuals with scarred skin.
- Individuals with infected, swollen, burned, or otherwise damaged skin. The device should only be used on healthy skin.
- Individuals who require a quiet environment or who are in poor overall health.
- Individuals with severe fractures, sprains, muscle pain, or acute pain (e.g., headaches).
- Individuals with a body temperature of 38°C (fever) or higher.
- Individuals experiencing acute inflammation, fatigue, chills, unstable blood pressure, or other acute symptoms.
- Patients currently undergoing treatment in medical facilities.
- Individuals with erythema (skin redness).
- Individuals with epilepsy.
- Individuals with active tuberculosis.
- Consult your doctor before use if you are using any of the following electronic medical devices:
 - Cardiac pacemakers or other implanted electronic devices
 - Artificial heart-lung machines or other life-sustaining medical equipment
 - ECG monitors, medical scanners, or similar diagnostic devices
- Do not use the massage mat on children or adolescents who are developing physically.

4. INTENDED USE

The **KLEGU** massage mat is intended to promote relaxation and relieve muscle tension by massaging the neck, back, lumbar area, hips and legs. It uses techniques such as rolling, kneading and shiatsu, combined with heat therapy and air compression, to relieve muscle tension and reduce fatigue. This device is commonly used at home as part of a personal wellness routine, especially after long periods of sitting, standing or physical activity.

The massage mat is designed for non-medical, home use and should not be considered a substitute for medical treatment. People with health conditions should consult a doctor before use and always follow the manufacturer's instructions.

INFO: The massage mat is designed for users approximately 155 – 188 cm tall, so the massage zones align correctly with the neck, back, lumbar area and hips.

There are no other uses for this product.

WARNING: Only use this product as described in this user manual. Any other use is considered improper and may result in personal injury or material damage. KLEGU is not liable for damage or injury resulting from improper or careless use.

5. THE PACKAGE

5.1 Contents

The following components are included in the scope of delivery:

- 1x **KLEGU** massage mat with integrated storage bag (hook-and-loop fastener and carry handle on the back; rolled up, the mat forms a compact bag)

- 1x AC power adapter
- 1x Remote control
- 1x User manual

5.2 Unpacking

Remove everything from the packaging and check for damage. If this is the case, do not use the massage mat. Be very careful when opening the package to avoid damage to the product.

After unpacking, check if the delivery is complete and if all parts are in good condition. If you are missing a part or if a part is damaged, please contact us by email: Info@klegu.nl.

Before use, remove all packaging materials and stickers or labels. Do not remove the rating plate or any safety warnings.

5.3 Safety during unpacking

Do not leave the packaging material lying around carelessly during unpacking. It may become a dangerous playing material for children. Do not dispose of the packaging material until the massage mat has been tried and is working.

WARNING: Plastic bags & foils. To avoid the danger of suffocation, keep plastic bags and foil away from babies and children. Thin foils may cling to the nose and mouth and prevent breathing.

WARNING: Strangulation hazard. Cords longer than 22 cm can be harmful to babies and very young children. Children can become entangled while playing with it. Do not leave the massage mat and power cord lying around during unpacking.

6. TECHNICAL DATA

NOTE: Due to continuing improvements on our product, the technical specifications are subject to change without notice.

Product name	Massage mat
Model no.	HY-8639
Brand	KLEGU
Massage technology	3D stretch and traction massage
Programmes	4 programmes: Auto / Stretch / Twist / Flow
Intensity levels	3 levels: Low / Medium / High
Massage zones	3 zones: Neck / Back / Hip (plus leg and calf compression)
Heat settings	3 settings: 35 °C / 38 °C / 42 °C
Suitable for body height	155 – 188 cm
Suitable surfaces	Bed, sofa or floor
Rated voltage	AC 100 – 240 V, 50/60 Hz (adapter); output DC 12 V / 3 A
Rated power	36 W
Operation method	Remote control

Auto shut-off	After 20 minutes
Weight	4.8 kg
Colour	Anthracite with white
Material	PU leather cover
Safety class	II
Dimensions	Unfolded: 150 × 56 × 5 cm Folded: 37 × 56 × 15 cm
Waterproof	No
Certifications	CE, EMC, LVD, RoHS

7. PRODUCT PARTS

The massage mat consists of the following main parts:

- Neck zone with heating (kneading/twist massage)
- Back zone (stretch and rolling shiatsu massage)
- Lumbar zone with heating and traction
- Hip airbags (air compression)
- Leg and calf airbags (air compression)
- Remote control
- AC power adapter (with cord)

8. THE CONTROL BUTTONS AND FUNCTIONS

NOTE: Remove any protective film from the remote control before using it.

The remote control is connected to the massage mat by a cord. The button functions are described below.

Power (ON/OFF):

Press the power button (the round button at the bottom) to switch the massage mat on or off. The indicator light at the top is lit when the device is on.

Start / pause (Play/Pause):

Use the central button (Play/Pause symbol) to start or pause the massage. During a pause, the air compression stops; the heating keeps working. Press again to resume.

4 Massage modes: Auto / Stretch / Twist / Flow

Select the programme using the buttons around the central start/pause button. On starting, the device begins in Auto by default. Press the desired programme:

- **Auto:** default automatic massage mode.
- **Stretch:** press this button to switch to Stretch mode.
- **Twist:** press this button to switch to Twist mode.
- **Flow:** press this button to switch to Flow mode.

3 Massage zones: Neck / Back / Hip

Use the Neck, Back and Hip buttons to select the zone(s) you want to massage.

Intensity: Low / Medium / High

Use the intensity button to switch between three levels; the three indicator lights show the selected level. Start at the lowest level and increase it gradually for a deeper massage. If you experience discomfort, reduce the intensity.

Heating: Heat1 / Heat2 / Heat3

Use the Heat1, Heat2 and Heat3 buttons to select one of the three heat settings — approximately 35 °C, 38 °C and 42 °C respectively. The heating zones (neck and lumbar) warm up gradually over a few minutes. Press the active heat button again to switch the heating off. Heat therapy can promote blood circulation and help reduce stiffness.

9. HOW TO USE

The **KLEGU** massage mat provides a massage experience for your neck, back, lumbar area, hips and legs.

Instructions for use:

1. Place the massage mat on a flat, firm surface, such as a bed, sofa or the floor.
2. Plug the massage mat into the mains socket.
3. Carefully lie down on the mat and align your neck, back, lumbar area and hips with the corresponding massage zones. The mat is designed for a body height of approx. 155 – 188 cm.
4. Switch on the massage mat using the ON/OFF button on the remote control.
5. Select the desired massage mode, zone, intensity and heating using the remote control.
6. Relax during the massage. After approximately 20 minutes (the auto shut-off time), the device switches off automatically. Press the ON/OFF button to start again.

10. PRODUCT FEATURES

Neck massage with heating:

- The neck zone combines a kneading/twist massage with heating to relieve tension in the neck and shoulders.

Back stretch & rolling shiatsu massage:

- The back zone provides a rolling shiatsu and stretch massage that penetrates deep into the muscle layers and promotes circulation.

Lumbar traction with heating:

- The lumbar zone provides traction and warmth to relieve tension in the lower back.

Hip and leg relaxation (air compression):

- Airbags in the hip and leg zone simulate kneading movements to relieve muscle aches and promote circulation.

Modes & settings:

- 4 automatic massage modes (Auto / Stretch / Twist / Flow), 3 intensity levels and 3 heating levels (35 °C / 38 °C / 42 °C).
- A pause function that lets you pause the air compression while the heating keeps working.

Heating and safety features:

- The heating warms up gradually over a few minutes for safe and comfortable use.
- A built-in timer function with automatic shut-off after 20 minutes for your safety. To restart, press the ON/OFF button again.

Cover:

- A removable cover keeps the massage mat clean and hygienic and can be taken off to make cleaning easier.

11. MASSAGE

11.1 Massage applications

- Relaxation massage to relieve fatigue and muscle tension. Perfect for home use, after a long day.
- Supportive relaxation for various complaints, including:
 - Tension and stiffness in the neck, shoulders and back
 - Lower back pain and lumbar tension
 - General fatigue and poor circulation
 - Stress and tension
 - Sleep problems
 - Tired, heavy legs

11.2 Benefits

Massage therapy is known for promoting relaxation and supporting overall well-being. The massage mat combines modern technology with proven massage techniques.

- The massage mat helps relieve tension in the neck, back and lumbar area.
- The rolling, kneading and warming action penetrates the muscle layers, leaving you feeling rejuvenated and relaxed.
- The air compression of the hip and leg airbags stimulates circulation.
- The massage mat helps relax and prepare your muscles.
- The massage mat relieves tension and reduces stress.
- Self-massage helps improve quality of life and promotes a more balanced body.

11.3 Tips

Usage tips for the best result:

- Wear comfortable, light clothing during use for comfort and hygiene.
- Always start with the lowest intensity and increase it gradually for a deeper massage.
- If desired, place a thin towel between your body and the mat to reduce the intensity or the heat.

Positioning for comfort:

- Lie down on a flat surface and align your body with the massage zones.
- Shift your body slightly up or down so the neck, back and lumbar zones align well.

Recommended use:

- Start with one session per day and build up gradually as your comfort improves.
- Ideal times of use: after work or before going to sleep.

Tips for medical conditions:

- For back or neck complaints, diabetes or similar conditions: start with the lowest intensity and a gentle mode, and consult your doctor.
- Turn off the air pressure or the heating if this is more comfortable.

12. FREQUENTLY ASKED QUESTIONS

What is the recommended duration of one massage?

A single massage usually lasts about 20 minutes. The device switches off automatically after 20 minutes.

What should I do after a massage?

Drink a glass of water within 30 minutes after the massage to promote circulation. Take your time getting up calmly.

Is the massage mat suitable for everyone?

Most people benefit from regular use. However, consult your doctor if you are in poor health, recovering from an illness, or pregnant (see “Medical considerations”).

What body height is the mat suitable for?

The mat is designed for users approximately 155 – 188 cm tall, so the massage zones align well.

13. PROBLEM SOLVING

We have compiled a list of possible causes and solutions. If you still have problems, please feel free to contact us.

WARNING: Do not disassemble the massage mat. Contact us if the cause of the malfunction has not been found.

Problem	Cause	Solution
The massage mat does not work.	The plug is not in the mains socket.	Insert the plug into the mains socket. If the remote control lights up but the mat does not work, please contact us.

The massage mat switches off suddenly.	The plug is not properly seated in the mains socket.	Insert the plug properly into the mains socket.
	The massage mat switches off automatically after 20 minutes of use due to the time limit.	Restart the massage mat.
	The overheating protection has been activated.	Wait until the mat has cooled down completely before using it again.
	Too much force was applied to the massage parts.	Keep your body relaxed during use. Do not stand or jump on the mat. Avoid unnecessary pressure on the internal mechanisms.

14. CLEANING

- Let the massage mat cool down completely before cleaning.
- Switch off the device before cleaning. Always unplug it from the mains socket.
- A slightly damp or dry cloth is recommended for cleaning the exterior. Dry with a soft cloth.
- Be careful that no liquid gets into the remote control or the connections, and do not apply too much pressure to the buttons.
- Do not immerse the massage mat in water or other liquids. Do not rinse the device with water, to prevent water from leaking into the electrical circuit and causing malfunctions.
- Under no circumstances use aggressive cleaning agents, cleaning alcohol, or other chemical solutions.
- Do not use abrasive objects for cleaning.
- The cover is removable and can be taken off to make cleaning easier. Clean the PU leather cover with a slightly damp cloth; do not machine wash, tumble dry, bleach or iron. Refit the cover only once it is completely dry.
- Cleaning must not be performed by children.

15. MAINTENANCE & SERVICE

- Apart from cleaning, the massage mat is otherwise maintenance-free.
- You may not reconstruct, self-repair, or disassemble the massage mat.
- Check regularly for damage and severe wear.

WARNINGS: Do not attempt to open the massage mat for maintenance. Removing the housing or incorrect assembly may cause electric shock, damage, short circuits, or fire.

If you open the housing for any reason, the warranty will be voided. The massage mat contains no user-serviceable parts.

Maintenance, modifications, and repairs may only be carried out by the manufacturer, a professional workshop appointed by the manufacturer, or a similarly qualified person. For maintenance or service, first contact us at Info@klegu.nl.

- For missing accessories such as the user manual, please contact us at the email address above.

- In the event of a defect, return the product in its original packaging to **KLEGU**, together with a description of the defect, proof of purchase, and all accessories.

16. STORAGE & TRANSPORT

- We recommend cleaning the massage mat before storing the device for a long time (see the chapter “CLEANING”).
- Let the device cool down completely before storing.
- Preferably store the massage mat in its original packaging.
- Do not fold or roll up the massage mat too tightly; this can damage the internal parts, heating elements and airbags.
- Do not store any wet parts, as mould or rust may form.
- Store the massage mat in a dust-free place that is not exposed to sunlight and humidity.
- Do not store the massage mat at extremely high or extremely low temperatures. This can damage the device.
- Store the device out of the reach of small children.
- Store the power cord loosely coiled and do not wrap it tightly around the device.
- During transport, protect the device from vibration and shocks; preferably use the original packaging.
- After transport or storage at very low temperatures, the massage mat must acclimatise to room temperature for approx. 30 minutes due to possible condensation. The device can then be used.

17. DISPOSAL & RECYCLING

This product was designed and manufactured with high-quality materials and components, which can be recycled and reused.

17.1 Disposal of the product

WEEE: This product complies with the Waste Electrical and Electronic Equipment Directive 2012/19/EU (WEEE).

The crossed-out wheeled bin symbol indicates that at the end of its life cycle, the product must be disposed of separately from ordinary household waste and taken to a separate waste collection centre for the recycling of electrical and electronic equipment, or returned to the retailer when purchasing a new, similar product.

Correct separate waste collection helps prevent possible negative effects on the environment and health. For more information about recycling this product, contact your local waste collection service or the shop where you bought the product.

When disposing of the massage mat, safely cut through the power cord so that the device can no longer be used.

17.2 Disposal of the packaging material

Follow the instructions and markings on the packaging. Make sure that plastic bags, packaging, etc. are disposed of in an environmentally friendly way and collected for reuse.

18. WARRANTY & LIABILITY

KLEGU warrants that this product will be free from defects in materials and workmanship for one year from the date of purchase. This warranty applies to products used for personal purposes and not for commercial or rental purposes. **KLEGU** will, at its sole discretion, replace or repair this product or any defective part with a new or refurbished product or part. If the product is no longer available, a replacement may be made with a comparable product of equal or greater value. No refund will be granted.

This warranty is valid for the original purchaser from the date of initial purchase and is non-transferable. Proof of purchase is required to obtain warranty support.

This warranty does not cover:

- Negligence or misuse of the product.
- Use with incorrect voltage or current.
- Damage possibly caused by pets.
- Damage resulting from drops or impacts.
- Water damage.
- Use contrary to the operating instructions.
- Disassembly, repair, or modification by anyone other than **KLEGU**.
- Damage caused by accidents or disasters, such as fire, flood, earthquake, war, vandalism, or theft.
- Damage due to improper maintenance and storage.
- Damage caused by incorrect connections.
- Normal wear and tear.

In the case of damage caused by failure to observe this manual, the warranty is void. We accept no liability for consequential damage. We accept no liability for material damage or personal injury caused by improper use or failure to follow the safety instructions.

KLEGU reserves the right to make changes to the manual, colours, specifications, accessories, materials, and models at any time without notice or obligation.

All rights reserved. No part of this publication may be reproduced, distributed, or transmitted in any form or by any means without the prior written permission of the publisher, except in the case of brief quotations embodied in critical reviews and certain other non-commercial uses permitted by copyright law.

KLEGU

V.O.F Klegu

Burseman 13

7156JM Beltrum

The Netherlands

Info@klegu.nl

Made in China

Versie / Version: **1.0 – 04-06-2026**