FITAGE

USER MANUAL

ENGLISH	5
NEDERLANDS	18
FRANÇAIS	32
DEUTSCH	47
ESPAÑOL	62
ITALIANO	75

SCAN THE QR CODE

to view the digital manual
in all languages

EN WE WANT TO LISTEN TO OUR CUSTOMERS

Your opinion is everything for FITAGE. Your experience drives us to help contributing to improve people's health and fuels our community. If you like your FITAGE product, please share your experience with us!

1. Go to Amazon / Bol.com / fitage.com
2. Click "Your Orders" and search for the FITAGE product.
3. Click on "Write A Review" and tell us about your experience!

ANY PROBLEMS?

If you encounter any problems, please contact us via Email. We are available:
on Monday - Friday
09.00-17.00 (+1 GTM)
in English, Dutch or German.
Support Email: support@fitage.nl

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EN: Scan the QR code and discover more about our FITAGE community!

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FITAGE

GET IN SHAPE

Introducing the FITAGE Massage Gun Premium, a groundbreaking device crafted to enhance your health and well-being. This massage gun is designed to deliver deep tissue heated or cooled massage and promote relaxation, aiding in quicker recovery and improved performance. Whether you're an athlete striving to enhance your abilities or someone seeking relief from everyday stress and tension, the FITAGE Massage Gun Premium is the ideal solution to support your goals.

With its advanced features and easy-to-use design, this massage gun is essential for anyone aiming to reclaim their health. Don't hesitate, invest in your wellness today with the FITAGE Massage Gun Premium!

SPECIFICATIONS AND CERTIFICATIONS

Product Name	FITAGE Massage Gun Premium
Product Model	FITAGE-22
Message heads	8
Motor type	Brush Motor
Power supply	USB-C with EU Adapter
Input voltage	5V
Net weight	620 g
Rotation time	1800 – 2800 rpm
Noise level	40- 55db
Certifications	CE, RoHS

The product and its packaging are made from materials that can be recycled, which helps reduce waste and protect the environment. Please dispose the packaging properly by utilizing local facilities designated for collecting paper, cardboard, and light-weight packaging.

Proper disposal of your massage gun is essential to protect the environment and human health. Reach out to local facilities to dispose this electronic device. The Waste Electrical and Electronic Equipment (WEEE) Directive provides guidelines for the responsible handling and disposal of electronic devices.

WHAT'S INCLUDED

- 1x FITAGE Massage Gun Premium
- 8x Massage Heads in Mesh Bag
- 1x Massage Ball
- 1x Instructions Manual (EN/NL/FR/DE)
- 1x USB-C Cable with EU Power Plug
- 1x Suitcase

INSTRUCTIONS FOR USAGE

Percussion therapy delivers precise pulses of pressure to the muscles and tissues, easing muscle pain, stiffness, and soreness while enhancing flexibility. By targeting both superficial and deep muscle fibers, percussive stimulation offers thorough relief.

The FITAGE Massage Gun Premium is designed for deep muscle relaxation and is ideal for post-workout physical massages. It helps alleviate muscle stiffness and pain, reduces muscle fatigue, and enhances blood circulation to improve tissue flexibility.

Purpose and Benefits:

Suitable for physical massage, skincare, or post-exercise recovery.

1. Accelerates warm-up and recovery by delivering targeted pulses of pressure to the muscle tissue.
2. Reduces muscle soreness and stiffness with a powerful high-torque motor featuring Quiet Glide™ technology. Can be used to activate, soothe, or loosen muscles and fascia throughout the entire body.
3. Improves range of motion and flexibility. Percussive stimulation offers comprehensive relief by targeting both superficial and deep muscle fibers.

TOUCH DISPLAY

1. Heat/Cold Level Indicator
2. Decreasing speed button (-)
3. Increasing speed button (+)
4. Speed indicator (1 - 30 gears)
5. Cold switch
6. Heat switch
7. ON/OFF button
8. Battery indication

GENERAL OPERATING INSTRUCTIONS

When the device is turned off, flip the power switch located at the bottom of the product to turn it on. This will keep the display screen illuminated. To power off the device, slide the bottom switch to the "OFF" position, and the display screen will turn off.

TO ACTIVATE THE HEATING MODE:

- Install the hot/cold head.
- Turn on the device's power switch to illuminate the display screen.
- Press the heating button once:
blue color bar indicates a temperature of 38°C.
- Press the heating button again:
orange color bar indicates a temperature of 40°C.
- Press the heating button once more:
red color bar indicates a temperature of 42°C.
- There are three temperature settings in total.

TO ACTIVATE THE COOLING MODE:

- Install the hot/cold head.
- Turn on the device's power switch to keep the display screen illuminated.
- Press the cooling button once:
blue color bar indicates a temperature of 20°C.
- Press the cooling button again:
orange color bar indicates a temperature of 18°C.
- Press the cooling button once more:
red color bar indicates a temperature of 16°C.
- There are three temperature settings in total.

TO ACTIVATE THE REGULAR MASSAGE MODE:

- Install any of the regular massage heads.
- Turn on the device's power switch to keep the display screen illuminated.
- Use the "+" button to increase the intensity level of the product.
- Use the "-" button to decrease the intensity level.

HOW TO USE THE MASSAGE GUN

1. Start massaging the muscles by setting the massage gun to the lowest setting.
2. Place the massage gun on the muscle and hover over it with the massage gun. Gradually increase the speed if you want the massage to be stronger. You don't have to apply a lot of pressure; the massage gun does most of the work.
3. Massage back and forth across the muscle and avoid keeping the massage gun in the same spot for longer than five seconds. If you are massaging a more tense area, you can hold a little longer. Be sure not to get too close to harder bone areas.

DIFFERENT MASSAGE HEADS

ROUND HEAD

Palm, Back,
Buttocks
and Legs.

LONG TAPER

Eases deep
tissue pain
in areas like the
knees, joints,
palms, and soles.

FLAT HEAD

Back, Buttocks,
Legs, Calves.

MUSHROOM HEAD

Joints, Sensitive
Muscles

GYRO HEAD

Joints, Sensitive
Muscles

DISPERSED HEAD

Shoulders,
Buttocks

SEMI-CIRCLE HEAD

Legs, Buttocks,
Arms

HEAT & COLD HEAD

Heating and
cooling function

BENEFITS

EFFECTIVE AND RELIEVING


IMPROVES BLOOD
CIRCULATION

REJUVENATES

RELIEVES PAIN

NORMAL MESSAGE HEAD ASSEMBLY


1. Press the message head gently into the message gun

2. Gently pull out and remove the message head

HEAT AND COLD MESSAGE HEAD ASSEMBLY

This massage head can be assembled using a plug-in method. Secure the red vertical bar of the cold and hot head onto the piston notch. Ensure that it is inserted fully to the bottom.

Fix the red vertical bar on the piston notch

CHARGING

The charging port is situated at the bottom part of the product handle.

1. Charge the battery fully for two to four hours before using it for the first time.
2. Use the provided special USB cable to connect to the TYPE-C socket. Use a 5-volt charger (5V and 2A is recommended).

3. The battery power indicators show the charging and battery status. When all the battery power indicators on the massage gun are lit up, it indicates that the battery is fully charged.
4. The battery can be charged at any time and at any battery level, under specified environmental conditions.
5. Place the massage gun on charge if the battery level is less than 25%. It is not recommended to let the battery power drop below 25% for increased longevity.
6. The average running time is three to six hours after the device is charged, depending on the vibration frequency, function use and applied pressure during use.
7. Charging the device via a computer is not permitted.

MAINTENANCE

1. Clean the device with a slightly damp towel, then wipe and dry it with a soft cloth.
2. Store the device in its packing box if it won't be used for an extended period.
3. Charge the device at least once a month if it will not be used for a long time.
4. Before packing and storing the massage device for long periods, clean it thoroughly. It's advisable to take it out for charging at least once a month to prevent moisture damage.
5. The product is suitable for transportation by vehicle, train, ship, and air.
6. Secure the product during transportation to prevent violent vibrations, collisions, rolling, or falling. Handle it with care, keep it upright, and protect it from rain.

FREQUENTLY ASKED QUESTIONS

The massage gun won't turn on?

1. Make sure you have properly charged the device using the included power charger. The charger light will turn full on once the battery is fully charged.
2. Check that you press longer than two seconds on the On/Off button.
3. Excessive Temperature: If the massage gun's operating environment exceeds 40°C, it may activate overheating protection, preventing normal startup.

Solution: Allow it to cool to below 40°C for 2 hours before attempting to start it.

4. Battery Capacity is "0": If the battery is completely drained and left uncharged for an extended period, it may activate low power protection, causing the massage gun to fail to start.

Solution: Keep it turned off, connect a charger to initiate charging mode. Unplug the charger when fully charged, then turn on the massage gun for normal operation.

5. If the massage gun still won't turn on, please contact our customer support team at support@fitage.nl.

I can't charge my massage gun?

1. Ensure that the charging port is properly plugged in.
2. Excessive Temperature of the Massage Gun: If the main body of the massage gun becomes too hot, charging may be disabled even when the charger is connected.

Solution: Allow it to cool to below 40°C for 2 hours, then attempt charging again.

3. Excessive Temperature of the Charger: If the operating environment exceeds 40°C, the charger may activate overheating protection, preventing normal charging of the massage gun. **Solution:** Allow it to cool to below 40°C for 2 hours, then attempt charging again.

The massage head slips off and falls?

Check if the massage head is properly installed (refer to this manual for detailed instructions).

It suddenly stops after continuous operation for 10-15 minutes?

This may be due to the long-term running protection function activating, switching from running mode to standby mode. Restart by briefly pressing the power button; this is a normal occurrence (refer to the manual for detailed instructions).

Is it safe to use a massage gun for neck pain?

With caution, yes, you can use a massage gun on your neck. However, it's essential to assess your tolerance level first. Applying too much pressure could potentially exacerbate neck pain.

Is it safe to use a massage gun during pregnancy?

If you're pregnant, always consult your doctor before using a massage gun or any other therapeutic device.

Can a massage gun alleviate lower back pain?

Yes, a massage gun can be effective for relieving lower back pain. However, it's important to consult your doctor before using one.

How frequently should you use a massage gun?

It varies for each person. However, using it too often, directly on bones, or for too long in one area could exacerbate the issue. If you're new to massage guns, gauge your tolerance by starting with no more than one session per day.

SAFETY WARNING

To prevent damage, such as malfunctions, electric shock, injury, fire, death, and device damage, carefully read all instructions, warnings, and reminders provided in the product's operation manual before using it.

1. To prevent damage to your property or injury to yourself and to ensure that the product continues to function properly over time, please read and follow these instructions carefully while using the product.
2. This unit is not suitable for individuals who are 12 years old or younger.
3. In order to minimize the risk of injury, do not directly massage the skin with the device. Wear clean and dry clothing over the massage area and use gentle pressure for no more than 60 seconds to prevent abrasions.
4. Use high speed and high pressure settings only on soft tissue areas of the body to prevent pain and discomfort. Do not use high speed and pressure on the head or any bones to avoid injury.
5. Keep your fingers and hair away from the telescopic rod and motor heat dissipation hole to avoid pinching injury.
6. Use the device in a dry environment with normal room temperature to maintain its long-term capability.
7. Handle the device gently and avoid dropping or hitting it, as this can damage the device.
8. This device is for household use only and should not be used commercially.
9. Do not run or charge the device in an unmanned state.
10. Do not use the device to diagnose or treat any medical issues, and consult a doctor before using it if you have any health problems.

11. Do not use the device if you are pregnant, have diabetes, neuropathy, retinal damage, wear a pacemaker, or have undergone surgery.
12. If the device malfunctions, check the battery level first and charge it if necessary. Do not attempt to open or modify the device. Use the original charger and remove it after the battery is fully charged.
13. To avoid damaging the device, do not put it in water or let water enter the openings. If the device becomes wet, use a dry cloth to clean it and do not use soap or any other cleaning agents. Also, remember to clean the motor heat sink frequently to keep it free of dust and debris.
14. Do not use the unit continuously for more than one hour.
15. To prevent any damage to the device, do not use or charge it in high-temperature environments. Also, make sure to remove the charger from the unit before using it and after it has fully charged.
16. Charging the product using a computer is prohibited.
17. Use caution if you are sensitive to heat or cold, as the device may produce heat and cold with its hot & cold massage head.

WARRANTY

We provide a two-year repair warranty on our product, with exceptions including:

- Damage resulting from incorrect operation or transportation by the user.
- Unauthorized tampering or repair of the equipment.
- Failure to consult the instruction manual for assistance.
- Overuse of the product beyond its recommended capacity.
- Damage due to inadequate storage or maintenance by the customer.
- Lack of proof of purchase provided by the customer.

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