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
User Guide

Electric Teeth Cleaner

Product Overview

1. Working Tip

2. Cleaning Head

3. LED Light

4. Mode Indicator

5. Mode Area

6. Mode Button

7. Power Button

8. DC Charging Port

Flat Cleaning Head

Suitable for removing stains and polishing

Hook-Like Cleaning Head

Suitable for removing large tartar

Needle-Like Cleaning Head

Suitable for removing small tartar

Precautions

1. The device is for adults use only. Keep out of reach of children to avoid injury.

2. Do not soak the body of the device in water.

3. The device has a built-in rechargeable battery. Do not throw into fire or charge, use, or leave in a high temperature environment.

4. Do not touch the working tip while the device is on.

5. Do not use the device if it is visibly damaged.

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
Precautions

6. Rinse the body of device or wipe with wet cloth if there is any dirt on its surface.

7. Place the device in cool shadow. Avoid direct sunlight or high temperature.

8. Keep the device out of reach of children to avoid any hazard.

9. The device has a built-in rechargeable battery. To maintain its service life, please recharge the device every 3 months if it is not in use for a long time.

10. Never attempt to open or repair the device yourself. If any malfunction happens or it doesn't work, please contact the customer service center.

Application Methods

1. Choose and assemble a suitable cleaning head.

2. Clean the working tip with running water or alcohol.

3. Press the Power Button to turn the device on. The LED light will turn on automatically. The default mode is the Soft Mode.

4. Press the Mode Button to choose a mode. There are 3 modes: Soft Mode, Normal Mode and Strong Mode.

5. Look in a mirror, and apply the working tip to the needed part on your teeth to remove any tartar or stain.

6. After use, press the Power Button to turn the device off. The LED light will go out automatically. Gargle with clean water.

7. Rinse the cleaning head under water. Wipe dry or wait for air dry, then put the device in a clean place.

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