Parvale®

LED FACE MASK

User Manual

Beste Klant,

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Valued Customer,


Thank you for choosing our **PARVALE™ LED Face Mask**.

We sincerely appreciate your trust in our product and are dedicated to providing you with the best possible experience. Before using the mask, carefully read this manual for proper usage and operation. Keep this manual safe for future reference to ensure optimal use and care of your **LED Face Mask**.


Product Overview:


PARVALE™ LED Face Mask


Remote Control


Package Includes:


PARVALE™ LED
Face Mask x 1


Remote Control
x 1


Eye Protectors
x 2

Self-buckle
Straps x 2

Charging Cable x 1

Adapter x 1

User Manual x 1

Technical Specifications:

LED Lamp Bead:	240
Operating Time:	10/ 15/ 20/ 25 mins
Service Temperature:	0 - 40°C
Adapter Input:	100 - 240V
Adapter Output:	5V DC 3000MA
Frequency:	0.3A 50/ 60Hz
Wavelengths Of Light:	Red Light 630nm Blue Light 460nm Yellow Light 590nm Near-Infrared Light 850nm
Net Weight:	0.32 kg


Introduction:

The **PARVALE™ LED Face Mask** provides an advanced light therapy experience to improve skin health and rejuvenation. With cutting-edge technology, this mask delivers targeted photon energy that helps address various skin concerns, including fine lines, wrinkles, and uneven texture. Suitable for home use, the mask offers a convenient and effective solution for those seeking a non-invasive treatment that enhances overall skin radiance.

Features:

- Multi-Mode Light Therapy for Personalized Treatment
- 3 Intensity Settings for Customizable Energy Levels
- Infrared Light Control for Deep Skin Penetration
- Adjustable Treatment Duration with 5-Time Options
- 60 Luminous Matrices with 240-beam Light Sources
- Ergonomic Design for Comfortable Fit
- Safe and Easy to Use with User-Friendly Controls
- Suitable For Targeting Various Face Areas

Photon Interpretation:

1. Red Light: 630nm

Anti-ageing, promotes collagen, relieves the relaxation state, and increases elasticity.

2. Blue Light: 460nm

Germicidal acne bacteria, reduce inflammation and scars, increase skin oxygen, and brighten skin luster.

3. Yellow Light: 590nm

Balance skin texture, stimulate red blood ball to reduce red and fine lines, effectively dilute acne marks, and maintain skin stability.

4. Infrared Light: 850nm


Repair the skin, accelerate the skin metabolism ability, improve blood circulation, and increase skin moisture.

Getting Started:

Carefully unbox the product and confirm that all the parts listed in the **“Package Includes”** are available in the box.

Preparation:


1. Remove the dustproof film from the face mask.

⚠ Warning! When removed, the product can not be returned.


2. Attach the eye protectors to the face mask.

Note: The eye protectors are labelled with "**L**" for the left side and "**R**" for the right side. The letter is visible from the back of the mask for correct placement.


3. Secure the self-buckle straps on the mask by snapping them into the strap connectors.

Note: The straps are adjustable to ensure a secure and comfortable fit for any head shape. Lengthen or shorten them as needed to achieve the optimal fit.


Operation Instructions:

1. **Power Button:** Long press for 1-2 seconds to turn the LED mask on or off.
2. **E Button (intensity increase):** Short press to adjust the intensity. The mask features 3 intensity levels to customize brightness based on your preference.
3. **T+ Button (Time Increase):** Press to extend the duration of the treatment session.
4. **T- Button (Time Decrease):** Press to reduce the treatment duration. The mask supports 5-time settings: 10, 15, 20, 25, and 30 minutes.
5. **NIR Button (NIR-Infrared Activation):** Use this button to activate or deactivate the NIR-Infrared lamps as needed.
6. **M Button (Mode Selection):** Press to switch between 4 light source modes to customize your therapy session for different treatments.



How to Use?

1. Rinse your face to remove makeup, sunscreen or skin care products and dry it to ensure the LED light can penetrate your skin effectively.


2. Charge the remote control by following the instructions provided in the **"Charging Instructions"** section.

3. Connect the mask to the remote control by inserting its plugs into the designated slots and turning on the LED light therapy.



4. After preparing the LED face mask according to the **"Preparation"** section, position it securely on your face, then adjust the strap buckles for a comfortable and snug fit.

- Power on the mask using the remote control.

- Follow the **“Operation Instructions”** section to select your desired therapy mode, energy level, and duration. Start with low settings (gear 1-2) to build skin tolerance.

- Turn off the mask and disconnect it from the remote control after use.

Tips:


- Avoid increasing the light intensity, time, or frequency without assessing your skin's reaction.
- Use the masks for 15-25 minutes per session, 3-4 times per week, based on your skin's needs.

How to Charge the Device?


The LED Face mask is powered by the remote control, so the remote control must be charged before the LED mask is used.

Charging Instructions:

1. Unplug the LED face mask from the remote control.


2. Insert the charging cable into the designated charging port on the remote control.

3. Connect the adapter to the charging cable and plug it into a standard electrical outlet.

Note: The screen indicator will flash while charging and turn off when fully charged.

4. Once fully charged, disconnect the adapter from the outlet and unplug the charging cable from the remote control.

Note: For first use, fully charge the remote control for 6 hours. Ensure the device is fully charged before use for optimal performance.

Charging Indicators:

1. **Flashing Light:** The device is currently charging.
2. **Light off:** Charging is completed.

Contraindications:

- Do not use the product if you are pregnant or breastfeeding.
- This product is not suitable for children.
- If you have a physical, sensory, or mental illness, or are undergoing treatment for dermatological diseases (such as lupus erythematosus, photosensitive eczema, or albinism), consult with a healthcare professional before use.
- **Photosensitive disorders:** If you suffer from any condition that causes sensitivity to light, do not use this product.
- Avoid using this product on the same day as products containing photosensitive ingredients.
- If you have any genetic eye conditions or are susceptible to light sensitivity, do not use the mask.
- Do not use the product if you experience any severe skin reactions to light or experience discomfort while using it.

Note: Always consult a specialist before using the product if you have any underlying medical conditions.

Safety Instructions:

- Use the product strictly as instructed in the manual. Improper use may lead to equipment failure or accidents.
- Do not use the product for purposes other than those specified in the manual.
- Protect your eyes during use. Wear the provided eye mask and keep your eyes closed during the session.
- If you experience any discomfort while using the product, stop use immediately.
- This product is not waterproof. Avoid exposing it to water or moisture. If it comes into contact with water, do not use it.
- Keep the product out of reach of infants and children.
- Do not use the product if any parts are damaged.
- Do not dismantle, modify, or repair the product without authorization.
- Use the product 3-4 times per week, with each session lasting 15-25 minutes.
- Always perform the LED light sensitivity test before use to ensure compatibility with your skin.
- The LED mask should lie flat on the skin during use, without pinching or applying pressure to specific areas of the face. The mask should not be fastened too tightly, even by using the head straps. Too tight a fit may cause local overheating, friction or pressure points, which in some cases may result in skin reactions such as hyperpigmentation, especially in sensitive or pigmented skin types.

Frequently Asked Questions:

a. How do I use the LED Face Mask?

Simply adjust the straps for a comfortable fit, connect the mask to the remote control, and follow the operation instructions to customize the settings based on your skin's needs.

b. How often should I use the LED mask?

We recommend using the mask 3-4 times a week for optimal results, with each session lasting between 15-25 minutes. You can also use our usage schedule in the first phase: it is as follows; 1st week every 3 days 10 min. 2nd week every 2 days 10 min. 3rd week every 2 days 10-15 min. 4th week every day 10-15 min. When you have the desired result you can maintain it by using the mask 1-2 times a week for 10-15 min each time. The intensity we often set medium and can be turned up depending on how the reaction is from the skin, sometimes it is not necessary to use a higher intensity as the light can penetrate the skin well.

c. Can I use the mask while pregnant or breastfeeding?

No, the LED mask is not suitable for use during pregnancy or breastfeeding.

d. What are the benefits of using the LED Light Therapy Mask?

The mask promotes skin rejuvenation, help reduce fine lines and wrinkles, improve skin texture, and enhance overall radiance through non-invasive photon therapy.

e. Can I use the LED mask with makeup or skincare products on my skin?

It is recommended to remove all makeup, creams, and lotions before using the LED masks for the best results.

f. How do I know when the mask is fully charged?

The charging indicator will flash during charging and turn off once the mask is fully charged.

g. Is the product waterproof?

No, the LED mask are not waterproof. Avoid exposing them to water or moisture. Clean the mask without cleaning agents and use a slightly damp cloth. After this, immediately use a dry cloth; do not use sharp cleaning brushes or sponges. Use only a slightly damp cloth.

h. Can I use the mask if I have sensitive skin?

If you have sensitive skin, start with the lowest intensity setting and perform a test run for 1 minute. If discomfort occurs, stop use immediately and consult with a dermatologist.

i. Can you use multiple colors in one session?

Each color has its own effect, the effects of the colors can be seen in the manual. The red light is most often used together with the Infrared mode (by the way, this color cannot be seen with the naked eye). For example, red light is often used when one wants a reduction of fine lines around the eyes and a firming of the skin. It is possible to use multiple colors for other problems, however, we recommend using a maximum of two or three colors in succession. For each color, use the build-up phase first!

j. How do I clean the mask?

The mask is made of medical silicone and we recommend cleaning the mask with a slightly damp cloth, do not use detergent as it may irritate the skin. To do this, first disconnect it from the power source. Use a damp cloth and then dry it with a dry cloth. Do not use harsh detergents and brushes with metal or nylon bristles. Also, do not use sharp objects when cleaning. Before using the mask, cleanse your face thoroughly and make sure you do not have any cream or makeup on while using the mask.

k. When can I expect results?

With most of our clients, changes are seen after the first 4 weeks, a reduction in symptoms is then certainly desirable. You can keep repeating the usage schedule until the desired results are achieved. Once the result is visible you can maintain it by using the mask 1-2 per week.

Warranty:

We offer a 12-month warranty from the date of purchase for the LED Face Mask. This warranty covers any defects in materials or workmanship under normal use. Should you experience any issues with your product during this period, please contact our customer service team for assistance.

What's Covered:

- Defects in material or workmanship
- Functional issues arising from manufacturing defects

What's Not Covered:

- Damage caused by misuse, improper care or accidents
- Normal wear and tear
- Unauthorized repairs or modifications

To claim your warranty, please retain your proof of purchase and contact us with the issue details. Our team will guide you through the next steps for repair or replacement.

Note: The warranty does not cover damages caused by water exposure, overheating, or any use outside of the intended purpose as outlined in the manual.

Parvale®