EN-DE-NL-FR-ES-IT-DK-SE-PL-TR-CS-SK

BASIC FIRST AID MANUAL

Compatible with DIN 13157 = ELRO EHBO2500 | with DIN 13169 = EHBO5000

EHBO1500 + EHBO5000

Home is not a place, it's a feeling

Basic First Aid Manual

- Compatible with DIN 13157 = EHBO2500
- Compatible with DIN 13169 = EHBO5000

Initial Steps

1. Ensure your own safety.
2. Assess the scene: What happened? Is there danger?
3. Call emergency services (112) if needed.
4. Use disposable gloves before touching any wound.

Key First Aid Kit Contents & Usage

- Sterile compresses: Cover open wounds.
- Adhesive bandages: For small cuts or abrasions.
- Wound dressings (small/large): Control bleeding.
- Triangular bandages: Use as a sling or to stabilize limbs.
- Adhesive tape: Fix dressings securely.
- Gauze bandages: Hold dressings in place or support injuries.
- Emergency blanket: Prevent hypothermia or shock.
- Scissors: Cut clothing or bandages safely.
- Tweezers: Remove splinters or debris.
- Face masks: For hygiene during close contact.
- Moist cleaning wipes: Clean intact skin.
- First aid manual: For basic guidance.

Note: DIN 13169 kits contain double quantities for larger workplaces.

COMMON EMERGENCY SITUATIONS

Bleeding

- Apply direct pressure using a sterile dressing.
- Elevate the affected limb if possible.
- Secure with gauze or a bandage.

Burns

- Cool with lukewarm water for at least 10 minutes.
- Do not use ointments or creams.
- Cover with a sterile, non-stick dressing.

Sprain or Bruise

- Rest the injured area.
- Apply a cold pack (wrapped in cloth).
- Use a compression bandage if needed.

Unconscious but breathing

- Place the person in the recovery position (on their side).
- Monitor breathing continuously.
- Call 112.

No breathing (Cardiac Arrest)

- Call 112 immediately.
- Start CPR: 30 chest compressions followed by 2 rescue breaths (if trained).
- Use an AED (Automated External Defibrillator) if available.

Additional Tips

- Check the contents of your first aid kit every 6 months.
- Replace any used or expired items immediately.
- Store the kit in a clearly marked, easily accessible location.
- Encourage staff to take certified First Aid and CPR training.

Symbolen

Recycling and disposal: The WEEE symbol means that this product and its batteries must be disposed of separately from household waste. When this product reaches the end of its service life, take it to a designated waste collection point in your area to ensure safe disposal or recycling. Protect the environment and public health, use natural resources responsibly!

Read the user manual and keep it in a safe place for future use and maintenance..

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