
8-SENSOR BASIC

Composition Scale

Manual / Handleiding / Handbuch / Manuel

Index

Forword	4
Display structure	6
Place & use	7
Read before use	8
Safety warnings	8
Body composition	10
App download	15
Instructions for use	16
100% satisfaction guarantee	23

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Foreword

This scale is now your personal health assistant. In order to guarantee an accurate measurement and a long service life of the scale, please read this manual before using and keep it for further reference.

This Sacoma 8-Sensor Basic Composition Scale could help you keep track of your personal health data. It records and analyzes the information from your scale, then transfer the data to your mobile device via wireless transmission technology. You can view your health data at any time of the day by downloading the free app.

Thank you for your order!

Thank you for choosing our 8-Sensor Basic Composition Scale!
Your opinion matters. If you love our Composition Scale as much as we do, we would love to hear your feedback in a review.

Problems with your order?


The 8-Sensor Basic Composition Scale is finished with care and is delivered to your home in perfect condition. If this is not the case, you can contact our customer service via: info@sacoma.nl
(Response within 12 hours on business days)

Discover more products online

Sacoma |


Display structure

Place & use

- Place the scale on an even and dry surface. As soon as the LED screen shows 0.0 kg, you can step on the scale.
- The lower part of the scale provides conductance. Therefore, make sure that your feet are always symmetrical on the scale.

Read before use

The instruction manual must be considered an integral part of the product and must be kept for its entire life. If the product is transferred to another owner, its documentation must also be transferred in its entirety.

To ensure the safe and correct use of the product, the user must carefully read the instructions and warnings contained in the manual insofar as they provide important information regarding safety, user and maintenance instructions. Should the instruction manual get mislaid or you require additional information or clarifications, please contact us by email at info@sacoma.nl.

Body composition scales that measures the weight and calculates the body composition (percentage of body fat, water, muscle, and BMI) and the basal metabolic rate based on height, age and gender. Compatible with bluetooth 4.0 and higher technology, it allows you to automatically send the data to your smartphone so you can monitor your progress at any time.

Safety warnings

- Before using, make sure the device is intact with no visible damage. If you have any doubts, do not use the device and contact the retailer.
- Keep the plastic bag for packaging out of reach of children: suffocation hazard.
- This device must only be used for its intended purpose and in the manner indicated in the instructions. All other uses are considered improper and therefore hazardous.
- The manufacturer cannot be considered responsible for any damage due to improper or incorrect use.
- This device can be used by children on their own from the age of 8 upwards.

- Individuals with reduced physical, sensory or mental capacities or unskilled people must only use the product if previously instructed in its use, in safety and only if informed of the dangers connected to the product itself.
- Children must not play with the device.
- Cleaning and maintenance operations must not be performed by children without adult supervision.
- In the event of a failure and/or malfunction, do not tamper with the device. Always contact your retailer for repairs.
- Handle the scale with care, protect it against knocks, extreme temperature fluctuations, humidity, dust, direct sunlight and sources of heat. ☐
- Do NOT allow to penetrate into the body of the scale. This appliance is not waterproof. Do NOT leave and do NOT use the scale in water-saturated environments; humidity rates greater than 85% or contact with water or other liquids may cause corrosion or compromise proper scale operation.
- Do NOT step onto the scale with wet feet, or if the surface of the scale is wet: slipping hazard.

Attention! Before using this appliance

- The body fat, water, and skeletal muscle percentages and the basal metabolism calculated by the scale must only be considered to be indicative data. For further information, consult a doctor or dietician.
- Through the electrodes located on the platform, the scale sends a weak, low intensity electric current that passes through the body.
- The scale is not dangerous in any way.
- The scale is not recommended for the following individuals as the data detected by the scale may not be accurate:
 - people with a pacemaker or other electronic medical device
 - people with fever
 - children under 12 years of age
 - pregnant women
 - individuals undergoing dialysis, with oedema in the legs or osteoporosis
 - people who take cardiovascular medication

Safe use of batteries

- Remove the batteries if the scale is not used for long periods of time and keep them in a well-ventilated and dry place, at room temperature.
- Do NOT recharge the batteries if they are not rechargeable.
- Do NOT recharge rechargeable batteries using methods other than those indicated in the manual or using unsuitable equipment.
- NEVER expose the batteries to sources of heat or direct sunlight.
- Failure to observe these instructions may result in damage to and/or the explosion of the batteries.
- Do NOT throw the batteries into an open flame.
- The removal or replacement of batteries must always be performed by an adult.
- Keep the batteries far from the reach of children: swallowing batteries carries the risk of death.
- If batteries are swallowed, immediately contact a doctor.
- The acid contained in the batteries is corrosive. Avoid contact with the skin, eyes or garments.

Body composition

The body composition is made up of:

- lean mass: all components made up without fat (water, muscle, connective tissue, and other non-fatty organic components), bones, bowels, organs, and blood;
- fatty mass: all fats that can be extracted from adipose tissue and other tissues/organs;
- bone mass.

Through the electrodes located on the platform, the scale sends a weak, low intensity electric current that passes through the body. The different masses respond differently to the current: lean mass provides less resistance than the opposite. This body resistance is called bioelectrical impedance. The relative percentage of the various body masses varies based on gender, age, and physical stature.

Body fat %

Body fat percentage is the amount of body fat as a proportion of your body weight. The body fat percentage indicated in the table below represents the total fat that is calculated by adding essential fat and storage fat. Essential fat is necessary for various bodily functions. The percentage is 3% in men and 12% in women.

Storage fat gives energy, protects against physical trauma and insulates the body to preserve heat. These percentages vary with the age and gender of the individual as shown in the table below. In Fitdays™, there will be an automatic comparison of your measured body fat percentage against a range chart for your reference.

Bodyfat Table					
	Thin	Athletes	Fitness	Acceptable	Too high
Women	<14%	14–21%	21–25%	25–32%	>=32%
Men	<6%	6–13%	13–17%	17–25%	>=25%

*The data in the table is for references only.

Percentage of Body Water

The % of water in the body differs based on gender: in healthy adult women it ranges from 45% to 60% while in men it ranges from 50% to 65%. The percentage of water is greater in childhood and decreases with age. Body weight variations can occur in a short period of time (hours), a medium period (days) or over a long period (weeks).

Small variations in weight over a few hours is normal and due to water retention since extracellular water (which includes the interstitial fluid, plasma, lymph and transcellular fluid) is the only part capable of varying so rapidly.

When weight losses occur in just a few days, it is very likely that it is part of the total water. In Fitdays™, there will be an automatic comparison of your measured body water % against a range chart for your reference.

% Body water

	Men	Women
Normal	50-65%	45-60%

Body mass index

The body mass index (BMI) is the relation of body weight to height. The classification given below is correlated to the risk that body weight exerts on health.

When BMI is between 18.5 and 24.9 it means that the weight

BMI Table

BMI Value

Underweight	< 18.5
Normal	18.5 – 24.9
Overweight	25.0 – 29.9
Obese	>30

does not statistically impact health and is to be considered normal. Any deviation, above or below this BMI interval should be considered a health risk. In this case, we recommend consulting a physician.

Getting in the habit of following a balanced diet helps improve quality of life, promote health and prevent illness, increase vitality, and maintain ideal weight.

It is important to know that BMI does not allow accounting for weight variations due to muscle mass, bone mass, water content or fat accumulation. Therefore, a simple BMI calculation alone should not be misinterpreted. The evaluation of BMI should be made by a doctor or other specialists such as a nutritionist or dietician. In Fitdays™, there will be an automatic comparison of your measured BMI against a range chart for your reference.

Muscle Mass Percentage

Muscle mass is the weight of your muscles in your body. There are three different types of muscle: smooth, striated skeletal, and striated cardiac.

Smooth muscle is made of smooth fibres. These are found in the walls of the digestive tract, respiratory tract, and genital tract and in the wall of the venous and arterial blood vessels. Human striated skeletal muscle is in the skeletal and facial muscles.

Cardiac striated muscle is found in the heart. Muscles play an important part in the overall health, including maintaining regular body functions. Muscle mass starts to decrease with age. Keeping track of muscle metrics overtime, along with other body composition metrics, can provide indications of your health fitness.

This scale captures and displays in the Fitdays™ app 3 muscle metrics: muscle mass (lbs/kg), muscle rate (all muscle types) and skeletal muscle (%). On average, 40% of body mass is made up of skeletal muscles. The chart below shows a range of your muscle rate (all muscle types) based on age and gender.

% Muscle mass normal range		
Age	Men	Women
20–39	75–89%	63–75.5%
40–59	73–86%	62–73.5%
60–79	70–84%	60–72.5%

*The data in the table is for reference only.

Bone Mass Information

Bone mass is the weight of bone mineral in your body. It is important that you develop and maintain healthy bones by having a balanced diet and plenty of exercise. Keep track of your bone mass over time for any changes in the long term. Your bone mass tends to decrease after the age of 30. Below is a chart to show average estimated bone mass based on weight and gender.

In Fitdays™, there will be an automatic comparison of your measured bone mass against a range chart for your reference also.

Bone Mass – Average bone mass			
Weight	Male	Weight	Female
<143 lbs	5.9 lb	<110 lbs	4.3 lb
143–209 lbs	7.3 lb	110–165 lbs	5.3 lb
>209 lbs	8.1 lb	>165 lbs	6.5 lb

*The data in the table is for references only.

Recommendations for use

- Always weigh yourself on the same scale at the same time, preferably when you wake up while naked.
- It is important to weigh yourself with dry, bare feet positioned correctly on the electrodes (see figure).
- Also avoid contact between the legs so the scale operates correctly.
- Wait about 15 minutes after rising so water can distribute itself through your body.
- Scale precision could vary in the following cases:
 - a) high level of alcohol
 - b) high level of caffeine or drugs
 - c) after heavy meals
 - d) after rigorous sports activities
 - e) during illness
 - f) when pregnant
 - g) when dehydrated
 - h) with wet or unclean feet
 - i) swollen legs

Technical specifications

Capacity: 180 kg / 400 lb

Division: 50 g / 0.1 lb

4 disposable 1.5V AAA carbon-zinc batteries

LED display

Operating frequency: 2402–2480 MHz, effective radiated power: <0.9 dBm

Environmental conditions for use: +10°C +40°C; RH ≤85%

Environmental storage conditions: +10°C +50°C; RH ≤85%

App download

After you downloaded the App,
do not forget to turn on your bluetooth.

Instructions for use

This scale can be used in 2 ways:

1. Body Weight Measurement only
2. Body Weight and Body Composition Measurement with the Fitdays™ app.

1. Body weight measurement only

(does not require pairing to smartphone application)

1. Place the scale on a flat, sturdy surface.

Note: An uneven or disjointed floor (e.g., rug, carpet, linoleum) may compromise the accuracy of the reading.

2. With the scale off, step onto it with your feet parallel and your weight evenly distributed. The scale will begin measuring. Do not move until the scale has measured your weight. The weight will flash 3 times.

Note: Do not lean on anything while the scale measures your weight to avoid altering the reading.

3. Step off the scale. The display will continue to show the reading for a few seconds before automatically turning off.

2. Body weight and body composition measurement with the Fitdays™ App

(requires pairing to smartphone application)

When it is paired successfully, all body composition metrics will be available to be viewed on the Fitdays™ app.

First time connecting the scale to your smartphone

The scale sends body composition data to the smartphone and make it available on the Fitdays™ app. During the following operation, we recommend keeping the smartphone near the scale to keep the connection active.

List of compatible devices:

- iOS 8.0 and above
- Android™ version 5.0 and above

1. Turn on bluetooth on your smartphone settings.
2. Download the free Fitdays™ app from App Store or Play Store™. You can also scan the QR code with your smartphone's camera to direct you to the app in the Apps store.
3. Launch the Fitdays™ app.
4. Follow the instructions in the app to pair the Sacoma scale and to create a new user.

If the scale is paired successfully with the app, you will see the weight measuring on the app.

- If the weight is not measuring on the app, please make sure the device is paired by trying the following:
- In the app, go to **Device**. You should see an entry of "**Sacoma Scale**". If not, tap on **ADD Device** to add a device.
- Step on the scale and tap on **Scale > Searching bluetooth for pairing**. If the scale is found, it will appear with an option to connect. Tap on "**connect**" and it will confirm with "**DEVICE IS PAIRED**".

If a device is already linked to the app and you wish to add another one, go to [Account] > [My Device] in the lower right corner. Then, tap the [+] icon in the upper right corner to access the device management screen. Once the app detects the new device, tap the scale icon to complete the pairing process.

- In order to ensure the accuracy of measurement data, please fill in your basic information correctly (height, birthday, gender).
- The icon of Fitdays™ will automatically change form to Sacoma after the device once connected to your smart phone.

How to measure your body composition

1. Launch the Fitdays™ app.
2. Login your own account.
3. Select Device page, and choose your current using device.
4. Press Measure key.
5. Step onto the scale, taking care to position your bare feet on the electrodes, and hold the measuring hand bar correctly (please check the picture on the right or on page 1).
6. After few seconds, you can find your body composition data in the app of your smartphone.

Maintenance

Battery replacement

When the "LO" or battery symbol appears on the display, replace the batteries with the same type.

1. Open the battery compartment by pulling the tab on the cover.
2. Remove the worn batteries.
3. Insert the new batteries according to the polarity markings.
4. Close the battery compartment.

Maintenance

- Store the device in a cool, dry area.
- Store the device in a horizontal position
- Only clean the scale with a soft, dry cloth.
- Never use chemicals or abrasive products

Adding Users to the Fitdays App

The Fitdays App allows multiple users to track their measurements individually. Follow the steps below to add a new user.



Adding a New User

1. Open the Fitdays App and go to the Account section.
2. Tap on your profile name in the top-left corner. A menu will appear showing the current users.
3. Select "Add User".
4. Enter the required information for the new user:
 - Nickname
 - Date of Birth
 - Height
 - Target Weight
 - Select Mode – choose the appropriate measurement mode.
5. Tap Confirm to save the new user.

Switching Between Users

1. Open the Fitdays App and go to the Account section.
2. Tap on the profile name in the top-left corner.
3. Select the user you want to switch to from the list.
4. The selected user's data will now be displayed in the app.

Each user's data is stored separately, ensuring accurate tracking for multiple individuals.

Troubleshooting

Problem	Possible cause	Solution
The scale doesn't turn on.	You have not stepped onto the scales properly.	Step onto the scale with your feet parallel. Note: the scale will turn on only after stepping onto it.
	The batteries have run out.	Replace the batteries.
"LO" or the battery symbol is shown on the display.	The batteries are low.	Replace them.
The display shows "Err".	Weight greater than the scale capacity.	Step off the scale immediately so as not to damage it and check the maximum capacity.
The "Fitdays+™" app doesn't display any data.	The connection has not been enabled in the smartphone settings.	Connect your smartphone.
	The smartphone went into standby during weighing.	Keep the app active on your smartphone during weighing to check that data is sent instantly.
The "Fitdays+™" app only shows your body weight.	You did not step on the scale with bare feet or did not contact the electrodes correctly.	Step onto the scale with bare feet and make sure your feet and hand contact the electrodes correctly.
	The skin on your feet is too dry.	Slightly dampen your feet.

All weights are approximate having been rounded up or down to the nearest whole number. We guarantee accuracy within $\pm 1\%$ + 1 division.

DISPOSAL PROCEDURE

The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment (Dir. 2012/19/EU-WEEE).

At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m² without having to purchase a new, similar device.

This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same.

Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

To correctly dispose of batteries (Dir. 2013/56/EU, the Waste Batteries and Accumulators Regulations 2009) do not throw them into household waste but rather dispose of them as special waste at specialised recycling centres. For more information about disposal of used batteries, contact the store where you purchased the equipment containing the batteries, the municipality or the local waste disposal service.

This product is meant for household use. The conformity, certified by the CE-marking on the device, refers to the Directive 2014/30/EU on electromagnetic compatibility and on Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment. This product is an electronic device that was tested to assure, according to the current technical knowledge, that it does not interfere with other nearby located equipment (electromagnetic compatibility) and that it is safe if used according to the advices given in the instructions for use. In case of operation anomalies, do not keep using the device and, if necessary, refer to the manufacturer directly.

The manufacturer, Sacoma S.p.A., declares that the device is compliant with Directive 2014/53/EU.

Warranty

This appliance is guaranteed for 2 years from the time of delivery of the goods, or another longer term envisaged by the national legislation of the consumer's residence. This provision complies with Italian and European legislation. The Sacoma products are designed for home use and must not be used in public venues. The warranty only covers manufacturing defects and does not apply if the damage is caused by an accidental event, incorrect use, negligence or misuse of the product. Use only the accessories supplied; the use of different accessories may result in invalidity of the warranty. Do not open the unit for any reason; in the case of opening or tampering, the warranty is definitively voided. This warranty does not apply to parts subject to wear or to the batteries when supplied. After 2 years from delivery, or another longer term envisaged by the national legislation of the consumer's residence, the warranty expires; in this case, the technical assistance interventions will be carried out against a fee. Information on technical assistance, whether under warranty or for a fee, can be requested by contacting info@sacoma.nl.

No payment will be due for repairs or replacements of products that fall within the terms of the warranty. In the event of faults, contact the retailer. Do NOT send the appliance directly to Sacoma. All the operations under warranty (including those of replacement of the product or part thereof) will not prolong the duration of the original period of warranty of the product replaced. The manufacturer declines any liability for any damage that may, directly or indirectly, be caused to persons, property or animals as a result of the non-observance of all the requirements established in the relevant instructions manual and concerning, especially, warnings relating to installation, use and maintenance of the appliance. Sacoma, in its constant commitment to improving its products, is entitled to changing without any notice, in whole or in part, its products in relation to production requirements, without this entailing any liability for Sacoma towards its dealers.

100% satisfaction guarantee

We guarantee all our products are designed to the highest standard, we base our ethos on customer satisfaction. This means the customer experience and journey you have had with us should be exceptional. If you have experienced our customer service or product, please let us know by leaving us an honest review.

INFO@SACOMA.NL

COPYRIGHT © 2021 SACOMA. ALL RIGHTS RESERVED.

The material in this publication is protected under International and Federal Copyright Laws and Treaties, and as such, any unauthorized reprint or use of this material is strictly prohibited. No part of this book may be reproduced or transmitted in any form without the written permission of the author, except for the inclusion of brief quotations in a review. Reproduction or translation of any part of this work without the permission of the copyright holder is against the law.

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COPYRIGHT © 2025 SACOMA. ALLE RECHTEN VOORBEHOUDEN.
COPYRIGHT © 2025 SACOMA. ALLE RECHTE VORBEHALTEN.
COPYRIGHT © 2025 SACOMA. TOUS DROITS RÉSERVÉS.
Le fabricant, Sacoma S.p.A., déclare que l'appareil est conforme à la Directive 2014/53/EU.