
 **NORTHWALL**

User Manual / Handleiding
Northgun G1

 **GOLD SERIES**

Maintenance

The battery in this device is built in and is not intended to be removed or replaced by the user. If you experience any issues related to the power supply, please contact our support team directly.

Do not attempt to open, repair, or disassemble the device yourself. Any unauthorized modification or tampering with the device's internal components will automatically void the product warranty.

Warnings

Please read all safety instructions carefully before using this device.

Do not immerse the device in water or any other liquid.

Do not use the device continuously for more than 30 minutes. Allow the device to rest for at least 5 minutes before resuming use.

Do not use the device on or near the head or genital area.

This device is not intended to diagnose, treat, cure, or prevent any medical condition, and should not be used as a substitute for professional medical treatment.

Consult your physician before use if you are experiencing unexplained pain, undergoing medical treatment, or using any medical devices or equipment.

Stop using the device immediately if you experience discomfort, irritation, or unusual side effects.

Do not attempt to open, repair, or modify the device. Unauthorized tampering may result in damage, injury, and voiding of product warranty.

CE & RoHS Compliant

This device is fully compliant with CE and RoHS standards, ensuring it meets applicable European safety, health, environmental, and consumer protection requirements.

- CE certification confirms that the product complies with relevant EU directives and regulations related to safety, electromagnetic compatibility, and user protection.
- RoHS compliance ensures that the device is manufactured without the use of hazardous substances restricted under the Restriction of Hazardous Substances Directive.

The product has been tested and manufactured in accordance with internationally recognized quality and safety standards.

Disclaimer

Please consult your physician before using the Northgun G1 massage gun if you have any health concerns or if any of the following conditions apply to you: spondylolisthesis, spondylolysis, spondylosis, recent joint replacements, pregnancy, diabetes, epilepsy, migraines, hernias, the use of an IUD, having a pacemaker, recent surgeries, or the presence of metal pins, implants, or metal plates in the body.

For a safe and comfortable experience, always use the device as instructed and listen to your body. If you experience any discomfort or uncertainty during use, stop immediately and consult a medical professional.

Main Structure

IMPORTANT:

The massage gun shuts down after 10 min. of use as a safety measure. You can turn it back on by pressing the ON/OFF button.

Massage Heads

Use the heat/cold massage head to warm up muscles, reduce soreness, and support recovery before or after exercise.

Use the ball head for large muscle groups to provide a comfortable, full-body massage and reduce overall muscle tension.

Use the bullet head to target deep tissue, trigger points, and small muscle areas with precise massage pressure.

Use the fork head to massage around the spine, neck, and Achilles tendon without applying pressure directly on bones.



Use the air cushion head for sensitive areas and a softer massage experience with reduced impact on muscles and joints.

Use the stainless steel massage head for a firm, targeted massage and a smooth glide over the skin. Ideal for hot and cold therapy.

Operation

● When the Northgun G1 is turned off, screw in the massage head in a clockwise motion. You can remove the massage head by rotating it counterclockwise.

● Press the On/off button  for 2 seconds to turn on the massage gun. It will start at the first intensity level. Press the On/off button again to select the next gear. There are 5 intensity levels. Press the On/off button  for 2 seconds to turn off the massage gun. The massage gun shuts down after 10 min. of use as a safety measure. You can turn it back on by pressing the On/off  button.

● The 5 LED lights above the On/off  button indicate the 5 intensity levels.

● When the massage gun is in standby mode (level 0), it will automatically shut down after 20 sec.

Charging

Before using the device for the first time, fully charge the battery. A complete charge takes approximately 3 hours.

Make sure the device is powered off before charging.

Connect the charging cable to the USB-C charging port located at the bottom of the device.

Plug the other end of the cable into a 5V/2A USB power source, such as:

- a phone charger
- a USB adapter
- a USB port on your laptop or PC

Note: Using a charger with a lower output than 5V/2A will increase charging time. Most standard phone chargers are suitable for charging the device.

When charging is complete:

- The battery icon will stop flashing.
- All four LEDs will light up.

The battery can be recharged at any time, regardless of the remaining charge level.

Technical Specifications

Model: NG-G1: Northgun G1 Massage Gun

Speed Range: 2000–3000 RPM

Amplitude: 11 mm

Rated Power: 25 W

Dimensions: 244 × 147 × 73 mm

Battery Voltage: 14.6 V

Battery Capacity: 2600 mAh

Charging Time: Approximately 180 minutes

Warranty

Northwall provides a 2-year warranty from the date of purchase against manufacturing defects. If a defect occurs within the warranty period, we will repair or replace the massage gun free of charge. For warranty claims or support, please contact us at: info@northwallproducts.com

The warranty will be void under the following conditions:

- Damage caused by improper use or incorrect storage.
- Damage resulting from attempts to repair, modify, or disassemble the device.

Use the Northwall massage gun before, during, or after your workout to support muscle performance and recovery. The targeted percussion therapy helps stimulate muscle responsiveness, improve mobility and range of motion, and relieve tension in tight or fatigued muscles.

Sessions	Goal	Intensity Level	Pressure	Treatment Time	Movement
Muscle activation & warm-up	Prepare and activate muscles prior to exercise	Level 2 or 3	Light	25-30 sec.	Move steady along the muscle
Mid-workout mobility	Support recovery between exercises	Level 4 or 5	Moderate to high	10-15 sec.	Move slowly over the trained muscle
Pain relief	Helps relieve pain and improve mobility	Level 1	Very light	25-30 sec.	Move slowly over the affected area
Recovery	Reduce muscle tension and support recovery	Level 1 or 2	Moderate	90-120 sec.	Move slowly along the muscle

Every Northwall product is backed by a 2-year warranty and supported by our dedicated customer support team. If you experience any issues with your product, we are here to help. Depending on the situation, we will repair, replace, or refund your product under warranty. Simply contact us by email and our team will take care of the rest.

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