
NORTHWALL

Foot Massager

Handling
User Manual

DISCLAIMER

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DISCLAIMER

READ THIS BEFORE USE

- Do not use this massager near water.
- Keep the ventilation grille free of dust and dirt.
- Do not remove any screws or attempt to open the device yourself.
- For adult use only. Do not use if you have any injuries. Consult your doctor if in doubt.

DANGER: To limit the risks of electric shock, injury and fire, use this device only in accordance with the instructions below:

- Suitable for adults only.
- Never leave the Northwall foot massager switched on unattended.
- Never use this device in any way other than as directed.
- Carefully check the device for defects before each use.
- Never use water near this massager.
- Check treated body parts regularly and stop immediately in case of pain or discomfort.

GUARANTEE

If a manufacturing or a material defect of the Northwall massage device occurs within two years of purchase, Northwall will repair or replace the product free of charge. In the event of a defect, please contact us via email: info@northwallproducts.com

CONTACT

For questions about the Northwall foot massager or any other Northwall product, you can always contact the Northwall customer service. We are happy to assist you and will do our best to find a solution.

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
Specifications

1. Place the foot massager on a flat surface.
2. Plug the device into a power outlet.
3. Slide your feet into the device and press the power button.
4. The massage starts at the lowest intensity for 15 minutes.
5. Press the timer button to adjust the massage duration between 5 and 30 min.
6. Press the intensity button to switch between 5 different intensities.
7. To switch the heating function on or off, press the heating button. The heat indicator light will turn on or off.
8. Press the program button to choose between 4 different programs. Program 4 is the relax mode.
9. The relax mode can also be activated by pressing the relax button. The relax indicator light turns on or off. In relax mode, the kneading massage is switched off.

Programs

Mode 1: The massage heads rotate inwards and the airbags are active.


Mode 2: The massage heads rotate outwards and the airbags are active.

Mode 3: The massage heads rotate inwards without airbags.

Mode 4: Relax mode: only the airbags and heat are active. The relax indicator light illuminates.

Attention

Always disconnect the foot massager from the power outlet after use. You can neatly store the cable on the underside of the device.

Maintenance

Always ensure that the device is disconnected from the power outlet before cleaning it with a slightly damp cloth.

The inner sock can be easily unzipped and can be hand-washed at a maximum of 30°C.