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CE REP

Components

Device	1 pc	Charging cable	1 pc
User manual	1 pc	Scaler tips	3 pcs
Protective cap	1 pc	Dental kit	1 pc
Disassembly wrench	1 pc	Alcohol wipes	2 pcs

Specification:

Charging Current: $\leq 0.5A$ Operating frequency: $38kHz \pm 2kHz$
Charging Time: 3 hours Power: $< 1.67W$
Size: $25.5 \times 25.5 \times 210mm$


You will not feel any vibrations when holding the device; the vibrations are only noticeable when the tip makes contact with the teeth.

FAQ

- Q Why does the device sometimes turn off during use?
- A If the battery indicator flashes or the device shuts off immediately after powering on, the battery is nearly depleted. Charge the device before use.
- Q Is it normal for the tip and the device to warm up slightly during use?
- A Yes, the internal operation may cause a slight increase in temperature. This is normal.
- Q Is it normal for the device to make noise?
- A A slight noise may occur if the tip is loose or due to resonance. Tighten the tip using the included wrench.

User Guide

- For safety reasons, the device is not fully charged upon delivery. Charge it fully before first use.
- Remove the components from the packaging.

- Attach the scaler tip and secure it using the disassembly wrench.
- Press and hold the on/off button for 1.5 seconds to turn the device on. Repeat to turn it off.

- The device switches off automatically after 10 minutes
- If the battery indicator flashes or the device shuts off immediately, the battery is empty. Charge the device.

Using the Scaler Tip

- Place the tip gently against the teeth and move it vertically and horizontally along the gum line.
- Adjust the duration of use based on the condition of your teeth.

This device may be used by children aged 8 and above, and by persons with reduced physical, sensory, or mental capabilities, provided they are supervised or have received instructions. Children must not play with the device. Cleaning and maintenance must not be carried out by children without supervision.

The device must only be connected to a safe extra-low voltage supply as indicated on the device.

Not suitable for persons with a pacemaker. The tip must not come into contact with the gums or soft tissue.

In the event of a serious incident involving this medical device, it must be reported to the manufacturer immediately.

Clinical benefit: Removal of tartar from the tooth surface and maintenance of oral hygiene.

Note:

- Consult a dentist before using this device.
- Report serious incidents to the manufacturer and the relevant authorities.

Cleaning: Remove the tip after each use, rinse with clean water, and allow to dry. Replace the tip regularly. Do not store the device in a humid environment.

Operating environment: Use in a dry indoor space at room temperature.

Disposal: Dispose of the device, accessories, and packaging in accordance with local regulations.

Intended use: Thorough removal of plaque and tartar, and maintenance of oral hygiene.

REP Europe GmbH
Prinzenallee 740549, Düsseldorf,
Germany info@rep-europe.de

Komponenten

Gerät	1 Stk	Ladekabel	1 Stk
Bedienungsanleitung	1 Stk	Scaler-Spitzen	3 Stk
Schutzkappe	1 Stk	Dental-Set	1 Stk
Demontageschlüssel	1 Stk	Alkoholtücher	2 Stk

Specification:

Charging Current: $\leq 0.5A$ Operate Frequency: $38kHz \pm 2kHz$
 Charging Time: 3H Power: $< 1.67W$
 Size: $25.5 \times 25.5 \times 210 \text{ mm}$

You won't feel the vibrations when holding it in your hand; vibrations can only be felt when the device makes direct contact with the teeth

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