

Sinocare Safe-Accu

BLOOD GLUCOSE TEST STRIP

IVD For self-testing



- 0.6 μ L Blood Sample
- Fast result within 10 seconds
- Automatic applying blood sample

in vitro diagnostic use only.
capillary and venous whole blood



European Authorized Representation
Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands
Tel: (31) (0) 70 345-8570


Manufacturer:



SteriLance Medical (Suzhou) Inc.
No.168 PuTuoShan Road
New District, 215153
Suzhou, Jiangsu, P.R. China

Instructions for Use:

1. Load your lancing device and twist off the protective cap.
2. Consult the Instructions for use of your lancing device.
3. After use, please dispose of the used lancet in an approved sharps container according to the local sharp disposal regulations.



Made in China



Version : 01