



WORLD CLASS MONITOR

A&D is a Japanese manufacturer of WORLD CLASS blood pressure monitors with industry-leading accuracy and ease-of-use.



Wrap cuff comfortably around wrist and secure. Position the monitor at heart level.

Press START button to automatically inflate the cuff and begin measurement.

Remain still until your measurement appears on the digital display.

WHO Classification Indicator

User & Clock set button

Clock & Calendar



- IT**
- Guida alla corretta posizione del braccio per un'accurata misurazione
 - %IB/AFib per indicare la frequenza di aritmia Cardiaca/Rivelazione AFib
 - Segnalazione corretta applicazione del bracciale
 - Rilevazione movimento del paziente
 - Valutazione della pressione secondo le linee guida O.M.S.
 - Orologio e calendario
 - Memoria per 2 utilizzatori (60 x 2) + calcolo del valore medio

- ES**
- Guía de posición correcta para mediciones exactas
 - %IB/AFib para mostrar frecuencia en la arritmia cardiaca / detección AFib
 - Error del ajuste del manguito
 - Error de movimiento
 - Indicador de Presión Sanguinea
 - Reloj y Calendario
 - Memoria de 60 registros para 2 usuarios + visualización de promedio de mediciones

- ZH**
- 正確量度血壓位置指引
 - %IB/AFib指標顯示不規則心律或心房纖維性顫動發生之頻率
 - 手袖佩戴錯誤提示
 - 運動誤差提示
 - WHO分類顯示測量結果
 - 時鐘與日曆功能
 - 雙用家記憶 (60x2) + 平均讀數



Technology You Can Trust
Blood Pressure Monitor
UB-533

- Correct Position Guidance
- 2 Users



- Correct Position Guidance
- 2 Users

- EN**
- Correct position guidance for accurate measurement
 - %IB/AFib to show frequency of Irregular Heartbeat / AFib Detection
 - Cuff fit error indicator
 - Movement error indicator
 - WHO blood pressure classification indicator
 - Clock & calendar
 - 2 user memory (60 x 2) + average reading

- FR**
- Indicateur de bonne position pour une mesure précise
 - %IB/AFib pour indiquer la fréquence de l'arythmie cardiaque / Détection AFib
 - Indicateur d'erreur de position du brassard
 - Indicateur de mouvement du bras
 - Diagnostic de la tension selon l'OMS
 - Fonction Horloge et Calendrier
 - 2 mémoires utilisateurs (2x60) + la moyenne



Wrist cuff
13.5-21.5cm
5.3"- 8.5"

A&D Company, Limited
1-243 Asahi, Kitamoto-shi, Saitama 364-8585 JAPAN
<http://www.aandd.jp>

EC REP Emergo Europe B.V.
Prinsessegracht 20 2614 AP The Hague The Netherlands
Telephone: [31] (0) 70 345-8570 Fax: [31] (0) 70 346-7299

CE 0123 Designed by A&D Medical, Tokyo Japan, Made in China



UB-533

