

UK SAFE SUPER LUBE CONDOMS

Standard condoms with an anatomical shape and reservoir. Provided with extra lubricant for a pleasant, natural feeling during sex.

Length 190 mm **Nominal width** 56 mm

NL SAFE SUPER LUBE CONDOOMS

Standaard condooms met een anatomische vorm en reservoir. Voorzien van extra glijmiddel voor een prettig, natuurlijk gevoel tijdens de seks.

Lengte 190 mm **Nominale breedte** 56 mm

FR SAFE SUPER LUBE PRÉSERVATIFS

Préservatifs standard avec une forme anatomique et un réservoir. Fourmi avec un lubrifiant supplémentaire pour une sensation agréable et naturelle pendant les rapports sexuels.

Longueur 190 mm **Largeur nominale** 56 mm

DE SAFE SUPER LUBE KONDOME

Standardkondome mit anatomischer Form und Reservoir. Mit extra Gleitmittel für ein angenehmes, natürliches Gefühl beim Sex.

Länge 190 mm **Nominale Breite** 56 mm

ES SAFE SUPER LUBE CONDOMES

Condomes ultra delgados para sentir aún más. Tan Seguros como los condones usuales.

Largo 190 mm **Ancho nominal** 56 mm

Distributed by **CLS Healthcare BV**
Maziestraat 22A | 2514 GT | The Hague
The Netherlands
www.enjoysafe.com

Manufactured by: Thai Nippon Rubber Industry PLC
789/139 Moo 1 | Pinthong Industrial Estate
Nongkhram | Sriracha | Chonburi 20110 | Thailand

EC REP Obelis S.A. | Boulevard Général Wahnis 53 | 1030 Brussels | Belgium

EXP: 06/2023

LOT: Z10180703

www.enjoysafe.com

Follow us   #enjoysafe

USE SAFE CONDOMS

This package contains condoms that are important for your health. Use them to reduce the risk of pregnancy and most STDs. However, no method can offer a 100% guarantee. In the EU, condoms can only be sold if they carry a CE mark. This means these condoms comply with strict requirements. The use of condoms is at one's own risk. Safe BV cannot be held liable for any damage.

CONDOM SAFETY

Subject to correct use, condoms offer an added protection from pregnancy and most STDs. However, no method can offer a 100% guarantee. In the EU, condoms can only be sold if they carry a CE mark. This means these condoms comply with strict requirements. The use of condoms is at one's own risk. Safe BV cannot be held liable for any damage.

PAY ATTENTION

Always keep condoms in the original packaging, at a cool and dry place and out of reach of sunlight. No contraceptive offers 100% protection from STDs, HIV or pregnancy. Oil-based lubricants or massage oil can damage condoms. Therefore, only use lubricants that are appropriate.

INSTRUCTIONS FOR USE

Make sure the use-by date has not expired. Open the packaging and take out the condom. Be careful with sharp nails and/or objects! Squeeze the tip of the condom to remove any air. Check to which side the condom rolls down and then roll condom over the stiff penis after the foreskin has been pulled back. After ejaculation, pull the penis back before it gets soft. While doing so, hold the rim of the condom to prevent it from slipping off. Dispose of the condom. Never flush a condom down the toilet.

NEVER REUSE CONDOMS!

EXP: 06/2023

LOT: Z10180703

CE 0120

LATEX

2

NO FIRE

NO SHARP OBJECTS

NO OIL

NO MASSAGE OIL

NO REUSE

EN ISO 4074:2002



SAFE CONDOMS | SUPER LUBE | EXTRA LUBRICANT | 36