

Kernel M^{Multi}Check
Cholesterol

Reagent Composition for each
MultiCheck TestStrip:

- Cholesterol Oxidase 200
- Phosphotungstic Acid 2.5 mg
- Phosphomolybdic Acid 2.5 mg

Kernel M^{Multi}Check Cholesterol 10



Test strip
MultiCheck
 Cholesterol 10
 For use with
 MultiCheck Serial Reader

CE Marked
 Cholesterol (total) in
 human whole blood with
 MultiCheck Serial Reader

Kernel M^{Multi}Check Cholesterol 10

Sample For MultiCheck
 Only use with Kernel MultiCheck
 Serial Reader
 45-100 mg/dL



Kernel
 Instrumentation
 Europe