

precautions after hair removal

- It is recommended to avoid contact with water within 6 hours after the hair removal, then moisturizing water can be used to moisturize the irradiated area.
- Do not use perfumes, masks, or cosmetics containing frosted beads after skin depilation. Avoid UV exposure for 24 hours after depilation, and apply sunscreen if necessary.
- Sunbathing is likely to cause skin darkening or sunburn. Irradiation at this time will cause redness, swelling, and inflammation. Therefore, during the use of this product, be careful not to be exposed to the sun, and the time to repair the sunburn is not less than 1 week.
- Artificial tanning is not allowed during the hair removal operation.
- If you need plastic surgery after hair removal (only for the hair removal area), consult your doctor.
- If you feel unwell, consult your surgeon.

Cleaning and maintenance

In order to ensure the best results and the longest service life of the product, please clean the product before and after each epilation to avoid long-term accumulation of hair and dust. If you do not clean it properly, the product will not get its proper effect.

before starting cleaning, the device must be switched off and unplugged to allow the machine to

use a soft, clean and dry cloth to wipe and clean the mirror surface of the light outlet.

Keep the device in a dry, dust-free environment at -10~60°C.



Do not wash the product and accessories in faucets and dishwashers.

Do not use steel balls, abrasive cleaners (toothpaste, etc.) or corrosive liquids (gasoline, motor oil,

, etc.) to clean the product.

If the water does not clean well, you can use a few drops of high concentration

(not more than 96%) to clean the above parts.



9

User manual of hair removal instrument

Thank you for purchasing the product. For convenience, please read and keep this manual carefully before use.