



Wheelchair and Scooter Bags



VA137ST/VA137SST

Fitting and Maintenance Instructions

INTRODUCTION

Thank you for deciding to purchase the Wheelchair or Scooter Bag from Aidapt. When used correctly they are designed to give many years of reliable, trouble-free service.

Please examine your product for any visible damage before using/installing. If you see any damage or suspect a fault, please do not use, but contact your supplier.

INTENDED USE

Wheelchair/ Scooter Bags are suitable for use with most standard wheelchairs/scooters. They provide the user and/or carer with extra storage.

FITTING INSTRUCTIONS

Overlap the straps so that the hook and loop fabric fastens the two sections together to secure the bag to the armrest of the wheelchair/scooter.

CLEANING

Clean your Wheelchair or Scooter Bag using a mild detergent with a slightly damp, soft cloth only. Only use a mild detergent or disinfectant; never use an abrasive cleaner.

After cleaning, leave the bag somewhere warm and dry to air. Do not leave in direct contact with a heat source.

CARE & MAINTENANCE

Please make a safety check of the product at regular intervals or if you have any concerns.

IMPORTANT INFORMATION

Please exercise common sense and do not take any unnecessary risks when using this product; as the user you must accept liability for safety when using the product.

Please do not hesitate to contact the person who issued this product to you or the manufacturer (detailed below) if you have any queries concerning the assembly/use of your product.

**Aidapt Bathrooms Ltd, Lancots Lane, Sutton Oak, St Helens, WA9
3EX**

Telephone: +44 (0) 1744 745 020 • Fax: +44 (0) 1744 745 001

**Web: www.aidapt.com
Email: sales@aidapt.co.uk**