

ACCU-CHEK® Safe-T-Pro Plus





200 Sterile, single-use lancing devices with 3 depth settings

Penetration depth / Stechtiefe / Profondeur de piqûre / Profondità di puntura / Profundidad de punción / Profundidade de penetração / Prikdiepte / Instiksdybde / Stikkedybde / Insticksdjup / Pistosyvvyys / Βάθος τρυπήματος / Głębokość nakłucia / Hĺoubka vpichu / Hĺbka vpichu / Szúrásmélység / Globina vboda / Dubina uboda / Adâncime de înterpare / Дълбочина на убождане / Torkesügavus / Duršanas dziļums / Įdūrimo gylis:


1.3 mm; 1.8 mm; 2.3 mm

Gauge / Durchmesser der Nadel / Diamètre de l'aiguille / Diametro dell'ago / Diámetro de la aguja / Diâmetro da agulha / Diameter van de naald / Lancettens diameter / Nålens diameter / Neulan läpimitta / Διάμετρος αιχμής / Średnica igły / Průměr jehly / Priemer ihly / Άτμέρο / Premer igle / Promjer igle / Diametrul acului / Диаметър на иглата / Lantseti läbimõõt / Adatas diametrs / Adatos skersmuo:

Ø 0.63 mm (23 G)



Consult instructions for use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consúltense las instrucciones de uso / Consultar as instruções de utilização / Raadpleeg de gebruiksaanwijzing / Se brugsanvisningen / Se bruksanvisningen / Katso käyttöohjetta / Συμβουλευθείτε τις οδηγίες χρήσης / Patrz instrukcja obsługi / Čtěte návod k použití / Pozri návod na použitie / Nézze meg a használati útmutatót / Upoštečajte navodila za uporabo / Pogledajte upute za uporabu / A se consulta instrucțiunile de utilizare / Вижте инструкцията за работа / Vaadake kasutusjuhendit / Skatiet lietošanas instrukciju / Žiūrėti naudojimo taisyklės



STERILE R

CE 0123



Roche