

**SPECIFICATIONS**  
**ROMED MEDICAL UNDERPADS**  
**REF: UP60X60, UP60X90**



*Medical device conform MDR (EU) 2017/745*

**INTENDED USE**

The Romed disposable medical underpads are used to collect moisture and wetness from lying patients in hospitals or at nursing homes, for example postoperatively, with incontinency or after giving birth.

**MATERIAL**

Spun bond nonwoven fabric top layer, moderate absorbent fabric mid layer with superabsorbent polymer (SAP) and waterproof polyethylene (PE) back layer. It can absorb 1500ml liquid in saturation situation.

**DIMENSIONS**

Reference no.	Weight	Length	Width	SAP content
UP60X60	50 gram	60 cm	60 cm	4 gram
UP60X90	60 gram	60 cm	90 cm	5 gram

**SHELF LIFE**

5 years

**PACKING**

10 pieces in a polybag and 20 polybags (=200 pieces) in an export carton.